## InTouch

CLINICAL UPDATE

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- Weiser TG, Haynes AB, Molina G, et al. Size and distribution of the global volume of surgery in 2012. *Bull World Health Organ*. 2016;94(3):201-209F.
- World Health Organization. Global Guidelines for the Prevention of Surgical Site Infection. 2016; <http://www.who.int/gpsc/ssi-guidelines/en/> Accessed 10th November 2016.
- Gaines S, Luo JN, Gilbert J, Zaborina O, Alverdy JC. Optimum Operating Room Environment for the Prevention of Surgical Site Infections. *Surgical infections*. 2017;18(4):503-507.
- Magid B, Murphy C, Lankiewicz J, Lawandi N, Poulton A. Pricing for safety and quality in healthcare: A discussion paper. *Infection, Disease & Health*. 2018;23(1):49-53.
- World Health Organization. Global action plan on antimicrobial resistance. 2015; [http://apps.who.int/iris/bitstream/handle/10665/193736/9789241509763\\_eng.pdf?sequence=1](http://apps.who.int/iris/bitstream/handle/10665/193736/9789241509763_eng.pdf?sequence=1). Accessed 4th April 2018, 2018.
- Chauveaux D. Preventing surgical-site infections: Measures other than antibiotics. *Orthopaedics & Traumatology: Surgery & Research*. 2015;101(1, Supplement):S77-S83.
- Australian College of Operating Room Nurses. Standards for Perioperative Nursing in Australia. In: 14th ed. Adelaide, South Australia: Australian College of Operating Room Nurses; 2017.
- Association of Operating Room Nurses. Perioperative Standards and Recommended Practices For Inpatient and Ambulatory Settings. In: *Sterilization*. Denver, CO: Association of Operating Room Nurses; 2013.
- Anderson DJ, Podgorny K, Berríos-Torres SI, et al. Strategies to Prevent Surgical Site Infections in Acute Care Hospitals: 2014 Update. *Infection Control & Hospital Epidemiology*. 2014;35(S2):S66-S88.
- Loveday HP, Wilson JA, Pratt RJ, et al. epic3: national evidence-based guidelines for preventing healthcare-associated infections in NHS hospitals in England. *J Hosp Infect*. 2014;86 Suppl 1:S1-70.
- Berríos-Torres SI, Umscheid CA, Bratzler DW, et al. Centers for disease control and prevention guideline for the prevention of surgical site infection, 2017. *JAMA Surgery*. 2017;152(8):784-791.
- Hebden JN, Murphy C. Minimizing ambiguity to promote the translation of evidence-based practice guidelines to reduce health care-associated infections. *Am J Infect Control*. 2013;41(1):75-76.
- Dunn D. Examining the Significance of Infection Prevention Practices. *AORN Journal*. 2017;105(5):442-445.
- Spruce L. Back to Basics: Sterile Technique. *AORN Journal*. 2017;105(5):478-487.
- Parvizi J, Barnes S, Shohat N, Edmiston CE, Jr. Environment of care: Is it time to reassess microbial contamination of the operating room air as a risk factor for surgical site infection in total joint arthroplasty? *Am J Infect Control*. 2017.
- Septimus EJ, Schweizer ML. Decolonization in prevention of health care-associated infections. *Clinical Microbiology Reviews*. 2016;29(2):201-221.
- Young HL, Reese S, Knepper B, Miller A, Mauffrey C, Price CS. The effect of preoperative skin preparation products on surgical site infection. *Infect Control Hosp Epidemiol*. 2014;35(12):1535-1538.
- Beldame J, Lagrave B, Lievain L, Lefebvre B, Frebourg N, Dujardin F. Surgical glove bacterial contamination and perforation during total hip arthroplasty implantation: When gloves should be changed. *Orthopaedics & Traumatology: Surgery & Research*. 2012;98(4):432-440.
- Daeschlein G, Kramer A, Arnold A, Ladwig A, Seabrook GR, Edmiston CE, Jr. Evaluation of an innovative antimicrobial surgical glove technology to reduce the risk of microbial passage following intraoperative perforation. *Am J Infect Control*. 2011;39(2):98-103.
- National Health and Medical Research Council. Australian Guidelines For The Prevention and Control of Infection In Healthcare. 2010. <https://www.nhmrc.gov.au/guidelines/publications/cd33>. Accessed 6th July 2014.
- U.S. Food and Drug Administration. Guidance for Industry and FDA Staff: Medical Glove Guidance Manual. 2008; <https://www.fda.gov/downloads/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/UCM428191.pdf>. Accessed 17 January, 2018.
- de Oliveira AC, Gama CS. Evaluation of surgical glove integrity during surgery in a Brazilian teaching hospital. *American Journal of Infection Control*. 2014;42(10):1093-1096.
- Tlili MA, Belgacem A, Sridi H, et al. Evaluation of surgical glove integrity and factors associated with glove defect. *American Journal of Infection Control*. 2018;46(1):30-33.
- Goldman AH, Haug E, Owen JR, Wayne JS, Golladay GJ. High Risk of Surgical Glove Perforation From Surgical Rotatory Instruments. *Clin Orthop Relat Res*. 2016;474(11):2513-2517.
- Makama JG, Okeme IM, Makama EJ, Ameh EA. Glove Perforation Rate in Surgery: A Randomized, Controlled Study To Evaluate the Efficacy of Double Gloving. *Surgical infections*. 2016;17(4):436-442.
- Misteli H, Weber WP, Reck S, et al. Surgical glove perforation and the risk of surgical site infection. *Archives of surgery*. 2009;144(6):553-558; discussion 558.
- Li X, Li M, Li J, et al. Glove perforation and contamination in fracture fixation surgeries. *Am J Infect Control*. 2017;45(4):458-460.
- Jid LQ, Ping MW, Chung WY, Leung WY. Visible glove perforation in total knee arthroplasty. *J Orthop Surg (Hong Kong)*. 2017;25(1):2309499017695610.
- Harnoss JC, Partecke LI, Heidecke CD, Hubner NO, Kramer A, Assadian O. Concentration of bacteria passing through puncture holes in surgical gloves. *Am J Infect Control*. 2010;38(2):154-158.
- Hubner NO, Goerdts AM, Stanislawski N, et al. Bacterial migration through punctured surgical gloves under real surgical conditions. *BMC infectious diseases*. 2010;10:192.
- Leitgeb J, Schuster R, Yee BN, et al. Antibacterial activity of a sterile antimicrobial polyisoprene surgical glove against transient flora following a 2-hours simulated use. *BMC Surg*. 2015;15:81.

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31. Leitgeb J, Schuster R, Yee BN, et al. Antibacterial activity of a sterile antimicrobial polyisoprene surgical glove against transient flora following a 2-hours simulated use. *BMC Surg*. 2015;15:81.
32. Assadian O, Kramer A, Ouriel K, et al. Suppression of surgeons' bacterial hand flora during surgical procedures with a new antimicrobial surgical glove. *Surgical infections*. 2014;15(1):43-49.
33. Mischke C, Verbeek JH, Saarto A, Lavoie MC, Pahwa M, Ijaz S. Gloves, extra gloves or special types of gloves for preventing percutaneous exposure injuries in healthcare personnel. *The Cochrane database of systematic reviews*. 2014;3:CD009573.
34. American Academy of Orthopaedic Surgeons. Preventing the Transmission of Bloodborne Pathogens. 2012; [https://www.aaos.org/uploadedFiles/PreProduction/About/Opinion\\_Statements/advistmt/1018%20Preventing%20the%20Transmission%20of%20Bloodborne%20Pathogens.pdf](https://www.aaos.org/uploadedFiles/PreProduction/About/Opinion_Statements/advistmt/1018%20Preventing%20the%20Transmission%20of%20Bloodborne%20Pathogens.pdf). Accessed 4th April, 2018.
35. Johnson J, Osborne S. Surgical hand antisepsis, gowning and gloving - Not necessarily the way it has always been done. *ACORN*. 2016;29(2):52-54.
36. Landelle C, Kubilay Z, Damani N, et al. Do different intraoperative glove practices reduce surgical site infections? A systematic review. *Antimicrobial Resistance and Infection Control*. 2015;4(Suppl 1):O30-O30.
37. Lipsett PA. Surgical site infection prevention—what we know and what we do not know. *JAMA Surgery*. 2017;152(8):791-792.
38. Ali S, Wilson APR. Effect of poly-hexamethylene biguanide hydrochloride (PHMB) treated non-sterile medical gloves upon the transmission of *Streptococcus pyogenes*, carbapenem-resistant *E. coli*, MRSA and *Klebsiella pneumoniae* from contact surfaces. *BMC infectious diseases*. 2017;17(1):574.
39. Kahar Bador M, Rai V, Yusof MY, Kwong WK, Assadian O. Evaluation of the efficacy of antibacterial medical gloves in the ICU setting. *J Hosp Infect*. 2015;90(3):248-252.
40. Leitgeb J, Schuster R, Eng AH, et al. In-vitro experimental evaluation of skin-to-surface recovery of four bacterial species by antibacterial and non-antibacterial medical examination gloves. *Antimicrob Resist Infect Control*. 2013;2(1):27.31. Leitgeb J, Schuster R, Yee BN, et al. Antibacterial activity of a sterile antimicrobial polyisoprene surgical glove against transient flora following a 2-hours simulated use. *BMC Surg*. 2015;15:81.
32. Assadian O, Kramer A, Ouriel K, et al. Suppression of surgeons' bacterial hand flora during surgical procedures with a new antimicrobial surgical glove. *Surgical infections*. 2014;15(1):43-49.
33. Mischke C, Verbeek JH, Saarto A, Lavoie MC, Pahwa M, Ijaz S. Gloves, extra gloves or special types of gloves for preventing percutaneous exposure injuries in healthcare personnel. *The Cochrane database of systematic reviews*. 2014;3:CD009573.
34. American Academy of Orthopaedic Surgeons. Preventing the Transmission of Bloodborne Pathogens. 2012; [https://www.aaos.org/uploadedFiles/PreProduction/About/Opinion\\_Statements/advistmt/1018%20Preventing%20the%20Transmission%20of%20Bloodborne%20Pathogens.pdf](https://www.aaos.org/uploadedFiles/PreProduction/About/Opinion_Statements/advistmt/1018%20Preventing%20the%20Transmission%20of%20Bloodborne%20Pathogens.pdf). Accessed 4th April, 2018.
35. Johnson J, Osborne S. Surgical hand antisepsis, gowning and gloving - Not necessarily the way it has always been done. *ACORN*. 2016;29(2):52-54.
36. Landelle C, Kubilay Z, Damani N, et al. Do different intraoperative glove practices reduce surgical site infections? A systematic review. *Antimicrobial Resistance and Infection Control*. 2015;4(Suppl 1):O30-O30.
37. Lipsett PA. Surgical site infection prevention—what we know and what we do not know. *JAMA Surgery*. 2017;152(8):791-792.
38. Ali S, Wilson APR. Effect of poly-hexamethylene biguanide hydrochloride (PHMB) treated non-sterile medical gloves upon the transmission of *Streptococcus pyogenes*, carbapenem-resistant *E. coli*, MRSA and *Klebsiella pneumoniae* from contact surfaces. *BMC infectious diseases*. 2017;17(1):574.
39. Kahar Bador M, Rai V, Yusof MY, Kwong WK, Assadian O. Evaluation of the efficacy of antibacterial medical gloves in the ICU setting. *J Hosp Infect*. 2015;90(3):248-252.
40. Leitgeb J, Schuster R, Eng AH, et al. In-vitro experimental evaluation of skin-to-surface recovery of four bacterial species by antibacterial and non-antibacterial medical examination gloves. *Antimicrob Resist Infect Control*. 2013;2(1):27.

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Ansell Healthcare Products LLC  
111 Wood Avenue South  
Suite 210  
Iselin, NJ 08830 USA

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Ansell Healthcare Europe NV  
Riverside Business Park  
Blvd International 55  
1070 Brussels, Belgium

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Ansell Services (Asia) Sdn. Prima  
6, Prima Avenue,  
Block 3512, Jalan Teknokrat 6  
63000 Cyberjaya, Malaysia

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Ansell Limited  
Level 3, 678 Victoria Street  
Richmond, Vic, 3121  
Australia

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© 2021 Ansell Limited. All Rights Reserved.