

# BIOCLEAN

## ANTI-FOG SPRAY

For Autoclavable Goggles

# BAFS

### About

To address any fogging issues with our autoclavable goggles an application of BioClean Anti-fog Goggle Spray is recommended. The procedure for the successful application of the spray is as follows:

1. Autoclave the BioClean autoclavable goggles as usual
2. Just before donning, apply one or two sprays to the INSIDE of the lens
3. Using a sterile cleanroom wipe, dry the excess liquid from the inside surface of the lens with a spreading motion, covering the whole of the inside of the lens
4. Ensure the treated side of the lens is dry after wiping
5. Repeat the process after every autoclave cycle

### Specifications

**MATERIAL:**

Non-hazardous ingredients

**LENGTH:**

60ml

**COLOUR:**

Colourless

**STERILIZATION:**

Gamma Irradiated

### Features

- Reduces fogging
- Airless pump action spray
- Alcohol free
- Quick and easy application
- Suitable for all BioClean autoclavable goggles

**ALCOHOL FREE**

**QUICK AND EASY APPLICATION**



Billy Bioclean©



## Quality Standards

---

- Manufactured in a facility holding ISO 9001:2015 and PPER 2016/425 Module D

## Sterilization

---

Method: Gamma irradiation  
Minimum Dose: 25kGy

## Shelf Life & Storage

---

Three (3) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

## TO ORDER

### Re-order code

BAFS-601

### Size

60ml

---

### Packing

One bottle per sealed inner PE bag; one inner bag per second PE bag; 10 double bags per sealed outer PE bag; 10 outer bags per lined carton (100 bottles)

For additional information visit us at [www.bioclean.com](http://www.bioclean.com), or call us at +44 1638 663338

**Nitritex Ltd**  
Ground Floor, 15 Kings Court  
Willie Snaith Road  
Newmarket, Suffolk  
CB8 7SG  
United Kingdom

BAFS/1809/PDS12

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2019 Ansell Limited. All Rights Reserved.