# ISOLATOR & RABS GLOVES

## INSTRUCTIONS FOR USE

Available in 32 languages

EN	DE	FR	ES	JA	IT	NL
PL	DA	EL	FI	SV	CS	HU
LV	HR	BG	SL	ET	LT	MT
RO	SK	GA	PT	TR	MS	ZH
TH	AR	NO	KO			

### Europe, Middle East and Africa Region

Ansell Healthcare Europe NV  
Riverside Business Park  
Blvd International, 55  
1070 Brussels, Belgium

### UK Importer / Authorized Rep:

Ansell (U.K.) Limited,  
Building C, Willerby Hill Business Park,  
Willerby, Hull, United Kingdom  
HU10 0FE

Nitrilex Ltd,  
Ground Floor, 15 Kings Court, Willie Snaith Road,  
Newmarket, Suffolk, United Kingdom  
CB8 7SG

[ansell.com](http://ansell.com)


Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2023 Ansell Limited. All Rights Reserved.

V.1 2022-12

## CONTENTS

<b>SYMBOLS &amp; PICTOGRAMS.....</b>	<b>1</b>
<b>EN: INSTRUCTIONS FOR USE.....</b>	<b>2</b>
<b>DE: GEBRAUCHSANLEITUNG.....</b>	<b>3</b>
<b>FR: MODE D'EMPLOI.....</b>	<b>4</b>
<b>ES: INSTRUCCIONES DE USO.....</b>	<b>5</b>
<b>JA: 使用説明書.....</b>	<b>6</b>
<b>IT: ISTRUZIONI PER L'USO.....</b>	<b>7</b>
<b>NL: GEBRUIKSAANWIJZING.....</b>	<b>8</b>
<b>PL: INSTRUKCJE UŻYTKOWANIA.....</b>	<b>9</b>
<b>DA: BRUGSANVISNING.....</b>	<b>10</b>
<b>EL: ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....</b>	<b>11</b>
<b>FI: KÄYTTÖOHJEET.....</b>	<b>12</b>
<b>SV: BRUKSANVISNING.....</b>	<b>13</b>
<b>CS: NÁVOD NA POUŽITÍ.....</b>	<b>14</b>
<b>HU: HASZNÁLATI ÚTMUTATÓK.....</b>	<b>15</b>
<b>LV: LIETOŠANAS INSTRUKCIJA.....</b>	<b>16</b>
<b>HR: UPUTE UPORABU.....</b>	<b>17</b>
<b>BG: ИНСТРУКЦИИ ЗА УПОТРЕБА.....</b>	<b>18</b>
<b>SL: NAVODILA ZA UPORABO.....</b>	<b>19</b>
<b>ET: KASUTUSJUHISED.....</b>	<b>20</b>
<b>LT: NAUDOJIMO INSTRUKCIJOS.....</b>	<b>21</b>
<b>MT: STRUZZJONIJET MT GHALL-UŽU.....</b>	<b>22</b>
<b>RO: INSTRUCTIUNI DE UTILIZARE.....</b>	<b>23</b>
<b>SK: NÁVOD NA POUŽÍVANIE.....</b>	<b>24</b>
<b>GA: TREORACHÚSIDE.....</b>	<b>25</b>
<b>PT: INSTRUÇÕES DE UTILIZAÇÃO.....</b>	<b>26</b>
<b>TR: KULLANIM TALIMATLARI.....</b>	<b>27</b>
<b>MS: ARAHAN PENGGUNAAN.....</b>	<b>28</b>
<b>ZH: 使用说明.....</b>	<b>29</b>
<b>TH: คำแนะนำในการใช้งาน.....</b>	<b>30</b>
<b>AR: إرشادات الاستخدام.....</b>	<b>31</b>
<b>NO: BRUKSANVISNING.....</b>	<b>32</b>
<b>KO: 사용 설명서.....</b>	<b>33</b>

## SYMBOLS & PICTOGRAMS

(01)

ABCDE

(02)

VIRUS

(03)

ABCDEFGHIJKLMNPST

(04)

(05)

(06)

(07)

**UK  
CA** 0321 **UK  
CA** 8515

(08)

(09)

(10)

(11)

(12)

(13)

(14)


(15)

(16)


(17)

(18)

(19)

(21)

(22)

(23)

(24)

(25)

**CA XX.XXX****STERILE|R**

(28)

**LOT****BATCH****INT****REF**

(29)

(30)

(31)

(32)

**RETURN TO CONTENTS**

# EN INSTRUCTIONS FOR USE

**USE:** This Instruction for Use is to be used in combination with the specific information that appears on the first packaging. Gloves are designed as a single use two way barrier protection against cross contamination and also protect the hands mainly against chemical splash risks and comply with the standards shown by the relevant pictograms. **EXPLANATION OF SYMBOLS & PICTOGRAMS THAT MAY APPEAR ON GLOVES/PACKAGING (01) EN ISO 21420:2020 / ISO 15523-1:2016** – Please read the Instructions for Use, prior to using the products, or contact Ansell for more information. If a level X is mentioned under any of the pictograms, this means this test is not applicable and glove is not designed and therefore not to be used for this specific hazard. **(02) EN 388:2016+A1:2018 – A B C D E** = Protection from mechanical risks: A: Abrasion resistance (performance levels 0 to 4) B: Blade cut resistance (performance levels 0 to 5) C: Tear resistance (performance levels 0 to 4) D: Puncture resistance (performance levels 0 to 4) E: TDM ISO EN 13997 cut resistance (performance levels A to F) **(03) EN ISO 374-5:2016** – Protection against bacteria, fungi and virus. **(04) EN ISO 374-1:2016+A1:2018 TYPE A, B OR C** – Protection against chemical hazards: Type A = chemical breakthrough time > 30 minutes against at least 6 chemicals as per list defined in EN ISO 374-1:2016 Type B = chemical breakthrough time > 30 minutes against at least three chemicals as per list defined in EN ISO 374-1:2016 Type C = chemical breakthrough time > 10 minutes against at least one test chemical as per list defined in EN ISO 374-1:2016 (no code underneath the pictogram) A = methanol, B = acetone, C = acetonitrile, D = dichloromethane, E = carbon disulfide, F = toluene, G = diethylamine, H = tetrahydrofuran, I = ethyl acetate, J = n-heptane, K = sodium hydroxide, 40%, L = sulphuric acid, 96 %, M = nitric acid, 65%, N = acetic acid, 99%, O = ammonia, 25%, P = hydrogen peroxide, 30%, S = hydrofluoric acid, 40%, T = formaldehyde, 37%. **(05) ISO 18889:2019** – Protection against pesticides – Glove suitable when the potential risk is relatively low. These gloves are not suitable for use with concentrated pesticide formulations and/or for scenarios where mechanical risks exist. The pesticide shall not have the possibility to penetrate between the garment sleeve and the glove. If the overlap is less than approximately 50 mm between the glove and the sleeve, a glove with a longer length should be used. For pesticide resistance data, the duration of the test is not based on actual use time since the permeation test is an accelerated test in which the surface of the specimen is in constant contact with the test chemical. Although the duration of the exposure may be for a longer period during field application with a dilute formulation, the entire surface is not in constant contact with the test chemical. **WARNING!** Chemical resistance data has been assessed under laboratory conditions and relates only to the chemical tested. It can be different if used in a mixture. For gloves equal or longer than 400 mm, the chemical resistance data is based from samples taken 80 mm from the end of the cuff. The data may not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. Check that the gloves are suitable for the intended use. Conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Chemical permeation data, tested per EN 16523-1:2015, and degradation data, tested per EN 374-4:2013, are available upon request and/or via [ansell.com](http://ansell.com), through the Ansell product page/downloads criteria/chemical recommendation guides. **(06) EN 421:2010** – Protection against radioactive contamination. **(07) CE** = Product is compliant and certified to the requirements of the European Regulation on Personal Protective Equipment 2016/425. When MD is mentioned with the CE mark, this is also a Class I product under the Medical Device Regulation 2017/745. **(08) UKCA** = Product is compliant and certified to the requirements to Regulation 2016/425 on personal protective equipment, as amended to apply in GB. For Great Britain: Type-examination certificate (Module B) and conformity to type certificates based on quality assurance of the production process (Module D) for CE marking are used as the basis for applying a UKCA. **(09) UKCA 0321** = Product is compliant and certified to the requirements to Regulation 2016/425 on personal protective equipment, as amended to apply in GB. For Great Britain: Type examination certificate (Module B) and, where applicable, Supervised product checks (Module C2) or Conformity to type based on quality assurance of the production process (Module D) by Satra Technology Centre, Wyndham Way, Telford Way, Kettering, Northamptonshire, NN16 8SD, UK. To obtain the EU-or UK Conformity Declaration, please go to: [www.ansell.com/regulatory](http://www.ansell.com/regulatory) **(10) UKCA 8515** = Product is compliant and certified to the requirements to Regulation 2016/425 on personal protective equipment, as amended to apply in GB. For Great Britain: Type examination certificate (Module B) and, where applicable, Supervised product checks (Module C2) or Conformity to type based on quality assurance of the production process (Module D) by Centexbel International Ltd, 8 Northumberland Avenue, London, WC2N 5BY, UK **(11) FOOD PICTOGRAM** = Suitable for contact with foodstuffs. Products carrying this pictogram are in conformity with the European Regulations 1935/2004 and 2023/2006 as well as with all applicable National Regulations for Food-contact materials. **(12) KOHSA** = Product is compliant and certified to the requirements of the Korean Occupational Health & Safety Act legislation for PPE. **(13) UKRSEPRO** = Certified in Ukraine **(14) NON-STERILE** **(15) NOT MADE WITH NATURAL RUBBER LATEX** **(16) MADE WITH NATURAL RUBBER LATEX** **(17) SINGLE USE ONLY** **(18) KEEP AWAY FROM SUNLIGHT** **(19) KEEP DRY** **(20) AVOID OZONE** **(21) MAXIMUM STORAGE TEMPERATURE, 40°C** **(22) KEEP AWAY FROM HEAT AND RADIOACTIVE SOURCES** **(23) DO NOT USE IF PACKAGE IS DAMAGED** **(24) EXPIRY DATE** **(25) MANUFACTURED DATE** **(26) MANUFACTURER** **(27) CA XX.XXX** = Certificate of Approval, as certified to the requirements of the Brazilian Regulation (whereas xx.xxxx refers to the certificate number). **(28) STERILIZED USING IRRADIATION** **(29) LOT NUMBER** **(30) BATCH NUMBER** **(31) ANSELL INTERNAL NUMBER** **(32) CATALOG NUMBER** For more detailed information on the product's performance, please consult Ansell.

**PRECAUTIONS FOR USE:** 1. Before usage, inspect the gloves for any defects or imperfections such as holes, pinholes and tears. If the gloves are ripped or punctured during use, dispose of them immediately. If in doubt, do not use the gloves, get a new pair. 2. It is essential to keep all chemicals from contact with the skin, even if they are thought to be harmless. Ensure the chemicals cannot enter via the cuff. Remove the glove immediately if contaminated by a concentrated spill of pesticides. 3. Contaminated gloves should be cleaned or washed or wiped dry before removal. Avoid touching contaminated surfaces with bare hands. 4. Gloves should not come in contact with a naked flame nor used for protection against heat. 5. Gloves shall not be used for protection against ionising radiation. 6. Gloves suitable for contact with foodstuffs may show some migration against specific foodstuffs. Please obtain advice from Ansell or consult the Ansell Food Conformity declaration to know if specific restrictions apply and for which specific foodstuffs the gloves can be used. 7. If gloves are being used in explosive environments, please ensure they meet the EN 16350 requirements. Persons wearing these gloves should be properly earthed, e.g. by wearing adequate footwear & clothing. Warning: the gloves shall not be unpacked, opened, adjusted or removed whilst in flammable or explosive atmospheres. The electrostatic properties of the gloves might be adversely affected by ageing, wear, contamination and damage and might not be sufficient for oxygen enriched flammable atmospheres where additional assessments are necessary 8. For Sterile Products – Contents of the pouch are sterile unless the pouch is opened or damaged. Do not use if the pouch is opened or damaged. 9. For Medical Use - Lifetime of the Glove – For normal use, following routine inspection, Ansell recommends Examination gloves to be changed every 15 minutes or during patient care if moving from a contaminated body site to another body site and after caring for a patient.

**INGREDIENTS / HAZARDOUS INGREDIENTS:** Some gloves might contain ingredients which are known to be a possible cause of allergies in sensitised persons, who may develop irritant and/or allergic contact reactions. If allergic reactions should occur, obtain medical advice immediately. Warning: If gloves contain natural latex, this would be mentioned on the packaging. In that case, **THIS PRODUCT MAY CAUSE ALLERGIC REACTIONS, INCLUDING ANAPHYLACTIC RESPONSES** to sensitised people. For more information, please contact Ansell.

**CARE INSTRUCTIONS: STORAGE:** Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials. **CLEANING:** The gloves are single-use only and not designed to be laundered nor to be reused. Re-use of the glove could cause post contamination and postoperative complication. Cleaning and re-sterilization has not been validated for this product by the manufacturer. Product integrity may be compromised by any cleaning or re-sterilization process used. **DISPOSAL:** Used products which have been in contact with chemicals or contaminated with infectious or other hazardous materials such as residual pesticides should be disposed after each working shift and not reused. They should also be disposed once they show any signs of degradation during usage, such as tearing, holes, discoloration and weakening of the gloves. Dispose of according to Local Authority Regulations. Landfill or incinerate under controlled conditions. If product has been used in a clinical setting after use, the product should be incinerated or disposed as per the clinical waste disposal protocol.

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