

	EN	ES	IT	NL	PT
AR	FR	DE	CS	PL	TR

## INSTRUCTIONS FOR USE

EN



For more information,  
visit: [www.ansell.com](http://www.ansell.com)

### Europe, Middle East and Africa (EMEA) Region

Ansell Healthcare Europe NV  
Riverside Business Park  
Blvd International, 55  
1070 Brussels, Belgium  
☎ +32 2 528 74 00  
📠 +32 2 528 74 01

### North America Region

US Legal Manufacturer/Fabricant:  
Ansell Healthcare Products LLC  
2301 Robb Drive, Reno, NV 89523, USA  
☎ +1 732 345 5400  
📠 +1 800 800 0445

### Australia

Ansell Limited  
Level 3, 678 Victoria Street  
Richmond, Vic, 3121  
☎ +61 1800 337 041  
📠 +61 1800 803 578

Ansell, ® and ™ are trademarks owned  
by Ansell Limited or one of its affiliates.  
© 2016–2019 Ansell Limited.  
All Rights Reserved.

#### A. Use

This Instruction for Use note is to be used in combination with the specific information that is mentioned on the packaging enclosure.

The GAMMEX® Cut Resistant Liner is a single use glove liner intended to provide ANSMSEA 105 Cut Level Protection 2 against cuts when used with an inner and outer surgical glove. The glove liner should be worn between two surgical gloves.

These glove liners are designed to protect the hands against mainly cut risks. They are in conformity with, and are marked per the requirements, of the European Regulation 2016/425. They also comply with the applicable harmonised European Standards.

Please ensure the glove liners are used only for the designated purposes.

Gloves are available in sizes: 5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9 and 9.5

Explanation of symbols & pictograms:



A B C D E  
EN 388:2016  
Protection from  
mechanical risks

A: Abrasion resistance  
(performance levels 0 to 4)  
B: Blade cut resistance  
(performance levels 0 to 5)  
C: Tear resistance  
(performance levels 0 to 4)  
D: Puncture resistance  
(performance levels 0 to 4)  
E: TDM ISO EN 13997 cut resistance  
(performance levels A to F)



EN 420:2003 + A1:2009

Please read the Instructions for Use,  
prior to using the glove liners or  
contact Ansell for more information.  
To obtain the EU-Conformity  
declaration, please use  
[www.ansell.com/regulatory](http://www.ansell.com/regulatory)



Product is compliant and certified to the requirements of the European PPE Regulation 2016/425

**Warning:** the performance levels claimed for the glove liners are based on tests performed on the palm area of the glove liners.

PPE EU -Type examination certificate (Module B) by Centexbel Belgium (I.D. 0493), Technologiepark 7, B-9052 Zwijnaarde. Gloves marked with CE 2797 = EU Conformity to type based on quality assurance of the production process (Module D) by BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands.

#### B. Precautions for use

1. NEVER USE THESE GLOVE LINERS WITH CHEMICALS.
2. Glove liners should not be used when there is a risk of entanglement with moving machine parts.
3. Before usage, inspect the glove liners for any defects or imperfections. Avoid wearing damaged glove liners.
4. The glove liners should not come in contact with a naked flame.
5. No glove liner will provide complete protection against cuts. Use with care when handling sharp objects.
6. The glove shouldn't be used against serrated blades.

#### C. Ingredients / Hazardous ingredients

Some glove liners might contain ingredients which are known to be a possible cause of allergies in sensitised persons, who may develop irritant and/or allergic contact reactions. If allergic reactions should occur, obtain medical advice immediately.

#### D. Care instructions

**Storage:** Keep away from direct sunlight; store in a cool dry place. Keep away from ozone sources or naked flame.

**Cleaning:** This glove is for single-use only. These glove liners are not designed to be laundered.

#### E. Disposal

Used glove liners may be contaminated with infectious or other hazardous materials. Dispose of according to Local Authority Regulations. Landfill or incinerate under controlled conditions.



XS	S	M	L	XL
5½-6	6½-7	7½-8	8½-9	9½-10