ME&DI

INSTRUCTIONS FOR USE

ANSELL SURGICAL Cat IIa GLOVES

Available in 26 languages

v.1 2022-09

EN	ES	IT	NL	DA	NO	FR
DE	PT	EL	FI	UA	SV	HU
LV	PL	BG	SL	CZ	ET	LT
MT	RO	SK	TR	MS		

North America

US Legal Manufacturer/Fabricant:
Ansell Healthcare Products LLC
2301 Robb Drive,
Reno, NV 89523, USA
+1 732 345 5400
+1 800 800 0445

EU & UK Legal Manufacturer

Ansell Healthcare Europe NV
Boulevard International 55
Brussels
B-1070
Belgium

Middle East and Africa (EMEA)

Ansell Healthcare Europe NV
Boulevard International 55
Brussels
B-1070
Belgium

Latin America and Caribbean

Ansell Commercial México S.A. de C.V.
Blvd. Bernardo Quintana No. 7001-C,
Q7001 Torre II, Suites 1304,
1305 y 1306. Col. Centro Sur, C.P. 76079
Querétaro, Qro. México
+52 442 296 2050

Brazil

Ansell Brazil Ltda.
Rua das Figueiras 474 – 4º Andar,
Santo André, SP - Brazil
CEP 09080-300
CNPJ: 03.496.778/0001-21

Asia Pacific

Ansell Global Trading Center (Malaysia) Sdn Bhd
Prima 6, Prima Avenue,
Block 3512, Jalan Teknokrat 6
63000 Cyberjaya, Selangor, Malaysia
+60 3 8310 6688
+60 3 8318 6699

UK Responsible Person for Medical Products

Ansell (U.K.) Limited
Building C, Willerby Hill Business Park
Willerby, Hull, United Kingdom
HU10 6FE

Australia

Authorised Sponsor ANZ: Ansell Limited
Level 3, 678 Victoria Street,
Richmond, Vic, 3121
+61 1800 337 041
+61 1800 803 578

Switzerland

Authorised Representative for Switzerland:
Cosanum AG
Brandstrasse 28
8952 Schlieren
Switzerland

Applicable to Great Britain

UK CA 0321 For compliance as a PPE product

UK CA 0086 For compliance as a Medical device product

UK IMPORTERS

Nitritex Limited, Ground Floor, 15 Kings
Court, Willie Snaithe Road, Newmarket,
Suffolk, CB8 7SG, United Kingdom

Ansell (U.K.) Limited, Block C, Willerby Hill
Business Park, Willerby, Hull, HU10 6FE,
United Kingdom

ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2022 Ansell Limited. All Rights Reserved.

CONTENTS

SYMBOLS & PICTOGRAMS	1
EN: ANSELL SURGICAL CAT IIa GLOVES.....	2
ES: GUANTES QUIRÚRGICOS ANSELL CAT IIa.....	4
IT: GUANTI CHIRURGICI ANSELL DI CAT IIa.....	6
NL: ANSELL CHIRURGISCHE HANDSCHOENEN CAT IIa.....	8
DA: ANSELL OPERATIONSHANDSKER KAT IIa.....	10
NO: ANSELL KIRURGISKE HANSKER KATEGORI IIa.....	12
FR: GANTS CHIRURGICAUX ANSELL – CATÉGORIE IIa.....	14
DE: OP-HANDSCHUHE DER KAT. IIa VON ANSELL.....	16
PT: LUVAS CIRÚRGICAS DE CATEGORIA IIa DA ANSELL.....	18
EL: ΧΕΙΡΟΥΡΓΙΚΑ ΓΑΝΤΙΑ ANSELL ΚΑΤΗΓΟΡΙΑΣ IIa.....	24
FI: ANSELLIN KAT. IIa LEIKKAUSKÄSINEET.....	26
UA: ХИРУРГІЧНІ РУКАВИЧКИ ANSELL КАТЕГОРІЇ ІІА.....	28
SV: ANSELL OPERATIONSHANDSKAR KAT. IIa.....	30
HU: ANSELL IIa. KAT. SEBÉSZETI KESZTYŰK.....	32
LV: ANSELL ĶIRURĢISKO CIMDU KATEGORIJA IIa.....	34
PL: RĘKAWICE CHIRURGICZNE ANSELL KATEGORII IIa.....	36
BG: ХИРУРГИЧЕСКИ РЪКАВИЦИ ANSELL, КАТЕГОРИЯ ІІА.....	38
SL: KIRURŠKE ROKAVICE ANSELL RAZREDA IIa.....	40
CZ: ANSELL CHIRURGICKÉ RUKAVICE CAT IIa.....	42
ET: ANSELLI KIRURGILISED KINDAD KATEGOORIA IIa.....	44
LT: „ANSELL“ IIa KATEGORIJS CHIRURGINĖS PIRŠTINĖS.....	46
MT: INGWANTI KIRURĠIĊI ANSELL KAT IIa.....	48
RO: MĂNUȘI CHIRURGICALE CAT IIa ANSELL.....	50
SK: CHIRURGICKÉ RUKAVICE ANSELL KAT. IIa.....	52
TR: ANSELL CERRAHİ KAT IIa ELDİVENLER.....	54
MS: SARUNG TANGAN PEMBEDAHAN CAT IIa ANSELL.....	56

SYMBOLS & PICTOGRAMS

(01)

(02)

EN 388:2016
+A1:2018

ABCDE

(03)

EN ISO 374-5: 2016

VIRUS

(04)

EN ISO 374-1

ABCDEFGHIJKLMNPST

(05)

EN 421:2010

(06)

(07)

(08)

(09)

(10)

(11)

(12)

(13)

(14)

(15)

(16)

(17)

(18)

(19)

(20)

A. USE: This Instruction for Use is to be used in combination with the specific information that appears on the first packaging. Gloves are designed as a two way barrier protection against cross contamination at the surgical site and also to protect the hands mainly against chemical splash risks and comply with the standards shown by the relevant pictograms. **EXPLANATION OF SYMBOLS & PICTOGRAMS (01) EN ISO 21420:2020 / ISO 15223-1:2016** – Please read the Instructions for Use, prior to using the products, or contact Ansell for more information. If a level X is mentioned under any of the pictograms, this means this test is not applicable and glove is not designed and therefore not to be used for this specific hazard. **(02) EN 388:2016+A1:2018 – A B C D E** = Protection from mechanical risks: A: Abrasion resistance (performance levels 0 to 4) B: Blade cut resistance (performance levels 0 to 5) C: Tear resistance (performance levels 0 to 4) D: Puncture resistance (performance levels 0 to 4) E: TDM ISO EN 13997 cut resistance (performance levels A to F) **(03) EN ISO 374-5:2016** – Protection against bacteria, fungi and virus. **(04) EN ISO 374-1:2016 TYPE A, B OR C** – Protection against chemical hazards: Type A = chemical breakthrough time > 30 minutes against at least 6 chemicals as per list defined in EN ISO 374-1:2016 Type B = chemical breakthrough time > 30 minutes against at least three chemicals as per list defined in EN ISO 374-1:2016 Type C = chemical breakthrough time > 10 minutes against at least one test chemical as per list defined in EN ISO 374-1:2016 (no code underneath the pictogram) A = methanol – B = acetone – C = acetonitrile – D = dichloromethane – E = carbon disulfide – F = toluene – G = diethylamine – H = tetrahydrofuran – I = ethyl acetate – J = n-heptane – K = sodium hydroxide, 40% – L = sulphuric acid, 96 % – M = nitric acid, 65% – N = acetic acid, 99% – O = ammonia, 25% – P = hydrogen peroxide, 30% – S = hydrofluoric acid, 40% – T = formaldehyde, 37%.

B. WARNING! Chemical resistance data has been assessed under laboratory conditions and relates only to the chemical tested. It can be different if used in a mixture. For gloves equal or longer than 400 mm, the chemical resistance data is based from samples taken, 80 mm from the end of the cuff. The data may not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. Check that the gloves are suitable for the intended use. Conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Chemical permeation data, tested per EN 16523-1:2015, and degradation data, tested per EN 374-4:2013, are available upon request and/or via ansell.com, through the Ansell product page/downloads criteria/chemical recommendation guides. **(05) EN 421:2010** – Protection against radioactive contamination. **(06) CE** = Product is compliant and certified to the requirements of the European Regulation on Personal Protective Equipment 2016/425 and Medical Device Directive 93/42/EEC or European Regulation 2017/745. **(07) MEDICAL DEVICE (08) UKCA** = Product is compliant and certified to the requirements of the Personal Protective Equipment Regulation 2016/425, as brought into UK law and amended and UK Medical Devices Regulation 2002 (SI 2002 No. 618 as amended) PPE Type examination certificate (Module B) and Supervised product checks (Module C2) or Conformity to type based on quality assurance of the production process (Module D) by Centexbel Belgium (I.D. 0493), Technologiepark 70, B-9052 Zwijnaarde, for EU, and SATRA Technology Centre (I.D. 0321), Wyndham Way, Telford Way, Kettering, Northamptonshire, NN16 8SD, United Kingdom (Module B, Module C2 or Module D), for UK. For compliance as a Medical device product, the CE and UKCA mark is followed by a four digit code which refers to the identification number of the Notified / Approved Body used For EU: I.D. 2797 – BSI Netherlands. For UK: I.D. 0086 – BSI UK. To obtain the EU-or UK Conformity Declaration, please go to: www.ansell.com/regulatory **(09) SINGLE USE ONLY (10) A SINGLE STERILE BARRIER SYSTEM, STERILISED USING IRRADIATION (11) UA.TR.116** = Medical device Class IIa under Technical Regulation on medical devices, approved by Decree of Cabinet of Ministers of Ukraine No. 753 on October 2nd., 2013 under conformity assessment procedure by the Notified Body “UNI-CERT” LLC, with reference (UA.TR.116) **(12) MADE WITH NATURAL RUBBER LATEX (13) NOT MADE WITH NATURAL RUBBER LATEX (14) DO NOT USE IF PACKAGE IS DAMAGED (15) KEEP DRY (16) KEEP AWAY FROM SUNLIGHT (17) AVOID OZONE (18) LOT NUMBER (19) MANUFACTURING DATE (20) EXPIRY DATE** For more detailed information on the product’s performance, please consult Ansell.

C. PRECAUTIONS FOR USE 1. Before usage, inspect the gloves for any defects or imperfections such as holes, pinholes and tears. If the gloves are ripped or punctured during use, dispose of them immediately. If in doubt, do not use the gloves, get a new pair. 2. Do not reverse the gloves. 3. It is essential to keep all chemicals from contact with the skin, even if they are thought to be harmless. Ensure the chemicals cannot enter via the cuff. 4. Contaminated gloves should be cleaned or washed or wiped dry before removal. Avoid touching contaminated surfaces with bare hands. 5. Gloves should not come in contact with a naked flame nor for protection against heat. 6. Gloves shall not be used for protection against ionising radiation nor for use in containment enclosures. 7. Gloves suitable for contact with foodstuffs may show some migration against specific foodstuffs. Please obtain advice from Ansell or consult the Ansell Food Conformity

declaration to know if specific restrictions apply and for which specific foodstuffs the gloves can be used. 8. Contents of the pouch are sterile unless the pouch is opened or damaged. Do not use if the pouch is opened or damaged. 9. For Medical Use - Lifetime of the Glove – For normal use, following routine inspection, Ansell recommends Surgical gloves be changed every 60-90 minutes. For gloves used with chemotherapy, we recommend they are changed every 30 minutes, or immediately if damaged.

D. INGREDIENTS / HAZARDOUS INGREDIENTS: Some gloves might contain ingredients which are known to be a possible cause of allergies in sensitised persons, who may develop irritant and/or allergic contact reactions. If allergic reactions should occur, obtain medical advice immediately. Warning: If gloves contain natural latex, this would be mentioned on the packaging. In that case, **THIS PRODUCT MAY CAUSE ALLERGIC REACTIONS, INCLUDING ANAPHYLACTIC RESPONSES** to sensitised people. For more information, please contact Ansell.

E. CARE INSTRUCTIONS: STORAGE: Keep away from direct sunlight; store in a dry place at 25°C to a maximum of 35°C and keep in the original packaging. Keep away from ozone sources. If gloves are properly stored, as indicated above, they won't lose their performances and won't change the glove characteristics significantly. If gloves could be affected by ageing or storage, the expiry date is mentioned on the packaging materials. **CLEANING:** The gloves are single-use only and not designed to be laundered nor to be reused. Re-use of the glove could cause post contamination and postoperative complication. Cleaning and re-sterilization has not been validated for this product by the manufacturer. Product integrity may be compromised by any cleaning or re-sterilization process used.

F. DISPOSAL: After use, the product should be incinerated or disposed as per the clinical waste disposal protocol. Used products which have been in contact with chemicals or contaminated with infectious or other hazardous materials should be disposed after each working shift and not reused. They should also be disposed once they show any signs of degradation during usage, such as tearing, holes, discoloration and weakening of the gloves.

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