



Ansell

SAFETY SOLUTIONS CATALOG

Veterinary

www.ansell.com

CONTENTS

03

WHY ANSELL? ANSELL ADVANTAGE

04

SUSTAINABILITY AT ANSELL

05

PROPRIETARY TECHNOLOGIES

06

VET SAFETY SOLUTIONS

07

HAND PROTECTION SOLUTIONS

HAND FATIGUE

- MICROFLEX® Ultraform™ UF-524
- MICROFLEX® XCEED® XC-310
- MICROFLEX® NEOGARD™ C52

ALLERGY PREVENTION

- MICRO-TOUCH® NITRAFREE™ Blue

CHEMICAL PROTECTION

- AlphaTec® Solvex® 37-175
- MICRO-TOUCH® Nitrile

EXTENDED PROTECTION

- MICRO-TOUCH® Nitrile E.P.
- MICRO-TOUCH® Nitrile Sterile
- MICRO-TOUCH® Central Supply – SPD Plus ULNCS

VALUE PRODUCTS

- MICROFLEX® 92-134 Versatility
- MICRO-TOUCH® MICRO-THIN™ Nitrile
- MICROFLEX® COBALT™ N19

PUNCTURE RESISTANCE

- HyFlex® 11-755
- RINGERS® R068

13

SURGICAL PROTECTION SOLUTIONS

- GAMMEX® PI Hybrid (Surgical)
- GAMMEX® Non-Latex Sensitive (Surgical)
- ENCORE® Latex Textured (Surgical)
- ENCORE® Latex Micro (Surgical)
- ENCORE® Latex Acclaim (Surgical)
- ENCORE® Perry Style 42™ PF (Surgical)
- ENCORE® SENSI-TOUCH® PF (Surgical)

17

SANDEL® MEDICAL SOLUTIONS

SHARPS SAFETY SOLUTIONS

- Z-Friction Drape™ Neutral Zone
- SharpsRink™ Neutral Zone
- Z-Tray™ & Stretch-A-Tray Hands-Free Transfer Trays

STAFF ERGONOMICS

- Z-Slider™ Animal Transfer Sheet
- ErgoPlus™ Anti-Fatigue Mats
- Ergo-Step™ Stool
- Sit-Stand™ Stool

TABLE STRAPS

- STAT-STRAP™ Table Straps

DISPOSABLE SHEETS

- Draw Sheets

FLUID AND TRIP MANAGEMENT

- Trip-No-More™ Cord Cover
- DriFloor™ Absorbent Floor Pad

INFECTION CONTROL

- STAT-SORB™ Solidifier

SPECIMEN MANAGEMENT

- Specimen Handling System

SURGICAL SITE MARKING

- PREP-RESIST™ Skin Markers
- TIME OUT® Skin Markers & Permanent Markers

MOPS

- Rayon Mop
- STAT-MOP™
- Microfiber Mops

WASTE BAGS

- STAT-BAG™ Heavy-duty Waste Bags
- Z-Suture Saddle Bag Waste Bags

WHY ANSELL?

Our mission is to provide innovative and effective solutions for safety, well-being, and peace of mind, no matter who or where you are.

Every day, more than 14,000 people in 58 countries design, manufacture, and market products on which millions of workers and healthcare professionals rely on.

MARKET LEADERSHIP



No. 1 or 2 market share
in key segments
globally



Safety
assessments
for every industry



Sell
~12 billion gloves
per year



Protect more than
10 million workers
each day

ANSELL ADVANTAGE

Protection

- Protecting workers across continuum of care
- Broad portfolio base to meet all application needs, head to toe
- Advanced chemical protection

Innovative

- Worker Experience Innovation: Technology innovations focused on customer Comfort, Performance, and Protection
- Deep understanding of the ever-changing regulatory landscape, experts to agencies and ability to educate customers

Reliable

- Owned Manufacturing excellence
- Manufacturing & Warehousing Diversification allows for supply chain flexibility and mitigates supply disruptions
- Competitively priced product portfolio



High Quality

- Drives worker satisfaction & safety
- Clear differentiation (AQL, coating technology, in-house & external testing, etc.)
- Best-in-class controls at each stage of manufacturing process

Trusted Partner

- Superior customer service
- Clinical resources and education
- Guardian safety assessments
- 125-year legacy and long-term focus

Presence

- Global leader protecting more than 10 million workers every day
- Strong focus on corporate social responsibility

ANSELL SUSTAINABILITY FOCUS AREAS AND ACHIEVEMENTS IN FY22

Sustainability is the foundation of a new purpose fully integrated into our business strategy at Ansell: Thinking of People and Planet First.

People

We are dedicated to investing in people as a part of our sustainability efforts

- In FY22, we focused on employee health & safety, labor rights, recruitment & employee engagement, diversity & inclusion, and community engagement & investment
- 87% close out rate for FY22 SMETA non-conformances at our plants¹
- 74% of workers at our plants are represented by worker councils, worker committees, unions, or covered by collective bargaining agreements



In FY22, we recorded our **lowest MTI date in 10 years of 0.161**, and **improved LTI rate of 0.051**. Employees submitted 10,511 safety improvement ideas, an **increase of 50%** in engagement rates



Trained 49 suppliers on modern slavery and labour rights, as part of our Supplier Management Framework



Roll-out of Supplier Management Framework to drive meaningful change in our supply chain

Planet

We are committed to achieve Net Zero within our operations by 2040 (Scope 1 & 2)



Zero Waste to Landfill

Waste to landfill diversion rate exceeding 99% against Intertek's Zero Waste to Landfill criteria

- Plants in Brazil, Thailand and China were certified as Zero Waste to Landfill, bringing a total of 5 certified plants
- 96% of waste generated from our plants were diverted from landfills



Water Withdrawals

Renewed water stewardship commitments with a goal of reducing 35% of water withdrawals by FY25²

- We implemented Reverse Osmosis in two of our plants, Bangkok and ATL
- Ansell partnered with the World Resources Institute (WRI) Aqueduct Alliance to conduct a water risk assessment at our manufacturing plants



Energy Consumption

- **45% renewables** in our energy mix across our plants
- **100% renewable grid electricity** in Xiamen, Portugal, Kedah A & B, and Lithuania

Product

In FY22 we focused on reducing the carbon footprint of our products and completed Life Cycle Assessments to understand our products' baseline impacts



Polybag made from bio-blend material³

- Biodegradation by 90% within 180 days



3.5% reduction in shipping volume

- Completed Phase 1 of standardization and optimization of shipper case dimensions⁴



Digitization of IFU in NA/LAC regions

- Reduction of 55 g of CO₂ emissions per shipper case



Plastic-free packaging for selected industrial products⁴

- Reduction of inner packaging material by 70%
- 82% reduction in GHG emissions



EMEA/APAC IFU right sizing

- 43% paper and GHG emissions reduction
- 122 g of paper reduction per shipper case, or 281 g reduction in CO₂ emissions⁵ (calculated for most common shipper case of 144 gloves)

Collaborations



Joined the Renewable Thermal Collaborative, a leading coalition for executing solutions, including knowledge sharing and advocating policy, to decarbonise thermal energy generation.



Ansell joined leading global players in the rubber glove industry as a founding member of the Responsible Glove Alliance. The collaboration aims to drive transformational change on labor rights, and both prevent and remediate forced labour in the Malaysian glove industry.



Ansell was awarded a silver medal by EcoVadis, placing it in top 25% of all companies assessed and the top 6% of manufacturers within our applicable industry.⁶



Ansell is a proud new member of Climate Group's RE100 and EP100. We are committed to these missions and initiatives for a low-carbon future.



Achieved full alignment with the recommendations of TCFD.



First Lead Member from Oceania to join the CDP Supply Chain Program on Climate Change and Water Security, a critical step in supplier engagement by improving availability of data and establishing suppliers' commitments to move towards emissions reduction.

¹ As at August 31, 2022 ² Using FY20 as a baseline ³ Achieving +90% biodegradation within 180 days in industrial composting facilities and is certified compostable by TUV Austria, meeting EN13432 standards ⁴ For selected mechanical styles ⁵ Ansell's Life Cycle Assessment (LCA), scenario analysis ⁶ EcoVadis classified Ansell in "other manufacturing n.e.c. (nowhere else classified)" industry



In our pursuit of developing solutions for healthcare safety, we've always looked beyond what is currently feasible and set our sights on what should be feasible. The results are superior gloves and safety devices with ground breaking, proprietary technologies, some of which have been patented, that offer you and your patients optimal protection in any healthcare environment and for any user need.

ANSELL PROPRIETARY TECHNOLOGIES



A.R.T™ Pinhole Barrier Technology

This technology is proprietary process, integrated into the manufacturing line that helps add an extra level of robustness by minimizing possible defects that may later result in pinholes.



DERMASHIELD™ Donning Aid Technology

This technology is a proprietary inner coating for fast and easy damp or dry hand donning whether it is the first don or during intraoperative changes.



ERGOFORM™ Ergonomic Design Technology

This technology enables Ansell to design safety solutions that support musculoskeletal health during repetitive tasks to improve worker performance.



HYDRASOFT™ Skin Moisturizing Technology

This technology is an odorless, water-based coating applied on the inside of the glove that helps retain moisture and rehydrate dry skin to improve the wearers' skin condition.



PI-KARE™ Skin-friendly, non-sensitizing PI Technology

PI-KARE™ Technology enables the elimination of standard chemical accelerators known to cause allergic contact dermatitis for low dermatitis potential; and, has resulted in the world's first full range of skin-friendly, non-sensitizing PI gloves that are US FDA-510(k) approved.



SENSOPRENE™ Ultra Thin Neoprene Technology

This technology delivers a glove up to 30% thinner than Ansell's standard neoprene surgical gloves so wearers' experience a high degree of sensitivity without compromising on durability and the optimum skin protection needed from allergies and sensitivities caused by latex or chemical accelerators.



Whether in the clinic or in the field, veterinary professionals face all kinds of hazards. Beyond the typical safety risks like exposure to pathogens, veterinary professionals also face other hazards that can endanger both themselves and the animals they treat. Using the right PPE and safety equipment is crucial to keep yourself protected against those various hazards.

However, constant using the wrong kind of PPE and safety equipment can increase the risk of injuries. At Ansell, we go beyond regulatory compliance, offering diverse safety solutions tailored to protect you against the challenges from skin allergies to muscle overexertion. We aim to deliver quality products with good value while also improving overall safety and productivity.

Primary Hazards in Vet Industry



SKIN ALLERGIES




CHEMICAL EXPOSURE



MUSCLE OVEREXERTION



SHARPS INJURIES





HAND PROTECTION SOLUTIONS

HAND FATIGUE

Our gloves equipped with ERGOFORM™ Ergonomic Design Technology provides comfort and support to your hands. It is proven to reduce the amount of hand muscle effort while performing repetitive tasks.



NBR

Nitrile

MICROFLEX® Ultraform™ UF-524

Comfortable, thin nitrile gloves that offer affordable protection with an ergonomic design

- Designed with ERGOFORM™ Ergonomic Design Technology to support musculoskeletal health
- Thin design for superior tactile sensitivity
- Soft nitrile formulation forms to hands for comfort during extended wear
- Packed 300/dispenser to reduce waste

External Glove Surface : Textured Fingers
Finger Thickness : 0.10 mm / 3.9 mil
Palm Thickness : 0.06 mm / 2.4 mil

Glove Length : 230 mm / 9.1 inches
Freedom From Holes : 1.5 AQL
 (Inspection Level I)

Ordering Information:

XS	XS/S	S	S/M
UF-524-XS	UF-524-XSS	UF-524-S	UF-524-SM
M	M/L	L	XL*
UF-524-M	UF-524-ML	UF-524-L	UF-524-XL



[SCAN OR CLICK HERE TO LEARN MORE](#)

300 gloves per box ; 10 boxes per case | * 250 gloves per box ; 10 boxes per case / 2500 gloves per case



NBR

Nitrile

MICROFLEX® XCEED® XC-310

Ergonomically designed to reduce hand fatigue

- Designed with ERGOFORM™ Technology to support musculoskeletal health and increase worker productivity
- Thin mil design for exceptional tactile sensitivity
- Strong formulation prevents against rips and tears
- Smartly packaged 250 gloves per box with easy, single glove dispensing to reduce waste
- Low 0.65 AQL for allowable pinholes means exceptional barrier integrity

External Glove Surface : Textured Fingers
Finger Thickness : 0.11 mm / 4.3 mil
Palm Thickness : 0.07 mm / 2.8 mil

Glove Length : 240 mm / 9.5 inches
Freedom From Holes : 0.65 AQL
 (Inspection Level I)

Ordering Information:

XS	S	M	L	XL*
XC-310-XS	XC-310-S	XC-310-M	XC-310-L	XC-310-XL



[SCAN OR CLICK HERE TO LEARN MORE](#)

250 gloves per box ; 10 boxes per case / 2500 gloves per case
 * 230 gloves per box ; 10 boxes per case / 2300 gloves per case



NE

Neoprene

MICROFLEX® NEOGARD™ C52

Ergonomically designed with exceptional grip

- Designed with ERGOFORM™ Technology to support musculoskeletal health and increase worker productivity
- Certified for exceptional wet grip

External Glove Surface : Textured Fingers
Finger Thickness : 0.13 mm / 5.1 mil
Palm Thickness : 0.10 mm / 3.9 mil

Glove Length : 245 mm / 9.6 inches
Freedom From Holes : 1.5 AQL
 (Inspection Level I)

Ordering Information:

XS	S	M	L	XL
C520	C521	C522	C523	C524



[SCAN OR CLICK HERE TO LEARN MORE](#)

100 gloves per box ; 10 boxes per case / 1000 gloves per case

ALLERGY PREVENTION SOLUTIONS

Severe reaction to chemical and latex allergies can be fatal. We have a portfolio of gloves that are designed for comfort and reliable protection from Type I Latex Allergies and Type IV Chemical Allergies.



NBR
Nitrile

MICRO-TOUCH® NITRAFREE™ Blue

Provides comfortable and reliable protection with low dermatitis potential

- Manufactured without sulfur-based chemical accelerators, helping to protect users from Type I latex allergies and Type IV chemical allergies
- Thin and light weight formulation provides enhanced tactile sensitivity and comfort for longer wear time
- Unique berry blue color provides contrast, while masking blood and other stains
- Tested for use with chemotherapy drugs

External Glove Surface : Textured Fingers
Finger Thickness : 0.10 mm / 3.9 mil
Palm Thickness : 0.07 mm / 2.8 mil

Glove Length : 240 mm / 9.5 inches
Freedom From Holes : 1.5 AQL
(Inspection Level I)

Ordering Information:

XS	S	M	L	XL
313028060	313028070	313028080	313028090	313028100

200 gloves per box ; 10 boxes per case / 2000 gloves per case
* 180 gloves per box ; 10 boxes per case / 1800 gloves per case



SCAN OR [CLICK HERE](#)
TO LEARN MORE

CHEMICAL PROTECTION

When handling cleaning chemicals and other solvents, there's a high risk of injury if your PPE breaches. Our chemical protection solutions give you a robust barrier so you can keep on working without worry.



NBR
Nitrile

AlphaTec® Solvex® 37-175

Durable, comfortable chemical gloves, offering versatile hand protection

- Cotton flocking: Moisture-absorbent, for greater comfort
- 17 mil/0.46mm thickness: For comfortable, dexterous safety gloves
- Nitrile compound composition: EN 374 Type A chemical protection

Cuff Style : Gauntlet
Construction : Unsupported
Palm Thickness : 0.46 mm / 18.1 mil

Glove Length : 330 mm / 13 inches
Grip Design : Sandpatch

Ordering Information:

6	7	8	9	10	11
37-175-6	37-175-7	37-175-8	37-175-9	37-175-10	37-175-11

12/polybag; 12 polybags/case



SCAN OR [CLICK HERE](#)
TO LEARN MORE



NBR
Nitrile

MICRO-TOUCH® Nitrile

Provides exceptional fit and reliable protection

- Exceptional comfort and fit
- Tested for use with over 55 chemicals and approved for use with chemotherapy drugs†
- A better than industry standards AQL delivers consistent, high quality production

External Glove Surface : Textured Fingers
Finger Thickness : 0.09 mm / 3.5 mil
Palm Thickness : 0.07 mm / 2.8 mil

Glove Length : 240 mm / 9.5 inches
Freedom From Holes : 1.5 AQL
(Inspection Level I)

Ordering Information:

XS	S	M	L	XL	XXL*
6034300	6034301	6034302	6034303	6034304	6034305

200 gloves per box ; 10 boxes per case
* 180 gloves per box ; 10 boxes per case / 1800 gloves per case



SCAN OR [CLICK HERE](#)
TO LEARN MORE

† See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

EXTENDED PROTECTION SOLUTIONS

Veterinary professionals may treat aggressive animals, or may be exposed to bodily fluids or bloodborne pathogens. You need more robust protection that gives you sufficient protection. Our range of extended protection solutions offer protection over the wrist and forearm while providing a strong barrier.



NBR
Nitrile

MICRO-TOUCH® Nitrile E.P.

Durable design with an extended cuff for high-risk protection

- Extended cuff (295mm/11.6in) provides added protection over the wrist and forearm
- Excellent chemical and puncture resistance
- Tested for use with chemotherapy drugs†

External Glove Surface : Textured Fingers
Finger Thickness : 0.16 mm / 6.3 mil
Palm Thickness : 0.11 mm / 4.3 mil

Glove Length : 295 mm / 11.6 inches
Freedom From Holes : 1.5 AQL
(Inspection Level I)

Ordering Information:

XS	S	M	L	XL
6034050	6034051	6034052	6034053	6034054

100 gloves per box ; 10 boxes per case / 1000 gloves per case



[SCAN OR CLICK HERE TO LEARN MORE](#)



NBR
Nitrile

MICRO-TOUCH® Nitrile Sterile

Nitrile, sterile pair packed, powder-free examination glove featuring extended cuff

- Packed as a pair of gloves per pack
- 11.6" extended cuff for additional protection
- Gamma sterilized minimum 25kGy to provide sterility
- Tested for use with chemotherapy drugs†

External Glove Surface : Textured Fingers
Finger Thickness : 0.16 mm / 6.3 mil
Palm Thickness : 0.11 mm / 4.3 mil

Glove Length : 290 mm / 11.6 inches
Freedom From Holes : 1.5 AQL
(Inspection Level I)

Ordering Information:

S	M	L	XL
6034151	6034152	6034153	6034154

1 pair per pouch ; 4 boxes per case / 200 pairs per case



[SCAN OR CLICK HERE TO LEARN MORE](#)



NBR
Nitrile

MICRO-TOUCH® Central Supply – SPD Plus ULNCS

These long-sleeve disposable gloves provide 33% longer, up-to-the-elbow protection in high-risk situations

- 16in / 400mm long to extend protection to the elbow
- Tested against 10 commonly used chemicals in the Central Supply/SPD Department
- Robust 9mil thickness at the finger for enhanced protection against tears and puncture
- Disposable to minimize cross-contamination occurrences
- Offer cost savings and improved protection versus fluid resistant cuff sleeves

External Glove Surface : Textured Fingers
Finger Thickness : 0.25 mm / 9.8 mil
Palm Thickness : 0.17 mm / 6.7 mil

Glove Length : 400 mm / 15.7 inches
Freedom From Holes : 1.5 AQL
(Inspection Level I)

Ordering Information:

S	M	L	XL	XXL
ULNCS-1	ULNCS-2	ULNCS-3	ULNCS-4	ULNCS-5

50 gloves per box; 10 boxes per case / 500 gloves per case



[SCAN OR CLICK HERE TO LEARN MORE](#)

† See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

VALUE PRODUCTS

Veterinary professionals should not have to choose between value and quality. These glove offers tactility and comfortable protection across numerous applications. A product that provides value and quality without compromising your safety.



NBR
Nitrile

MICROFLEX® 92-134 Versatility

Offers tactility and comfortable protection across numerous applications

- Made of a soft nitrile formulation for comfortable fit and feel
- Lightweight construction for tactile sensitivity and fine motor control
- Designed for a broad range of applications across all industries
- Textured on the fingertips for secure handling of tools and materials

External Glove Surface : Textured Fingers
Finger Thickness : 0.11 mm / 4.3 mil
Palm Thickness : 0.07 mm / 2.8 mil

Glove Length : 240 mm / 9.5 inches
Freedom From Holes : 1.5 AQL
(Inspection Level I)

Ordering Information:

XS	S	M	L	XL
92134060	92134070	92134080	92134090	92134100

100 gloves per box ; 10 boxes per case / 1000 gloves per case



[SCAN OR CLICK HERE TO LEARN MORE](#)



NBR
Nitrile

MICRO-TOUCH® MICRO-THIN™ Nitrile

Nitrile, powder-free examination glove with a thin film formulation for unrestricted tactile sensitivity

- Excellent tactile sensitivity is ideal for palpating veins and starting IVs
- The textured surface of the glove promotes sure handling of materials and instruments
- 300 gloves per box reduces restocking and reorders

External Glove Surface : Textured Fingers
Finger Thickness : 0.10 mm / 3.9 mil
Palm Thickness : 0.06 mm / 2.4 mil

Glove Length : 230 mm / 9.1 inches
Freedom From Holes : 1.5 AQL
(Inspection Level I)

Ordering Information:

XS	S	M	L	XL
6034310	6034311	6034312	6034313	6034314

300 gloves per dispenser ; 10 dispensers per case / 3000 gloves per case



[SCAN OR CLICK HERE TO LEARN MORE](#)



NBR
Nitrile

MICROFLEX® COBALT™ N19

Multipurpose nitrile exam glove offering consistent comfort and durability

- Fully textured for enhanced grip
- Consistent, dependable protection
- Enhanced tactile sensitivity
- Tested for use with chemotherapy/cytotoxic drugs

External Glove Surface : Fully Textured
Finger Thickness : 0.11 mm / 4.3 mil
Palm Thickness : 0.10 mm / 3.9 mil

Glove Length : 240 mm / 9.5 inches
Freedom From Holes : 1.5 AQL
(Inspection Level I)

Ordering Information:

XS	S	M	L	XL
N190	N191	N192	N193	N194

100 gloves per box ; 10 boxes per case / 1000 gloves per case



[SCAN OR CLICK HERE TO LEARN MORE](#)

PUNCTURE RESISTANCE

When handling larger animals where risk of puncture or abrasions are much greater, tougher protection is needed. Our portfolio of rugged gloves are designed to protect you against large forces.



PU

Polyurethane

HyFlex® 11-755

Ultralight ANSI A5/EN ISO E-level cut-resistant gloves with touchscreen compatibility

- ANSI A5/EN ISO E cut protection: 4x greater cut resistance*
 - Ultralight formulation: All-day comfort and dexterity ensured
 - Touchscreen-compatible design: Improved productivity and safety
- *versus standard ANSI A2/EN ISO B-rated gloves*

Construction	: Knitted	Glove Length	: 225 - 285 mm
Cuff Style	: Knitwrist	Finishing	: Palm Dipped

Ordering Information:

6	7	8	9	10	11	12
11-755-6	11-755-7	11-755-8	11-755-9	11-755-10	11-755-11	11-755-12

12 pairs/bag, 12 bags/carton



[SCAN OR CLICK HERE TO LEARN MORE](#)



NBR

Nitrile

RINGERS® R068

Fully coated safety gloves designed for impact, cut and liquid protection

- Patented TPR design: For PPE gloves with strong impact protection
- Fully dipped nitrile coating: Liquid resistance and hand protection
- Touchscreen compatibility: For consistent protection and usability

Construction	: Knitted	Grip Design	: Sandy Nitrile
Cuff Style	: Knitwrist	Finishing	: Fully Dipped

Ordering Information:

8	9	10	11	12	13
R068-8	R068-9	R068-10	R068-11	R068-12	R068-13

10 pairs/bag, 6 bags/carton



[SCAN OR CLICK HERE TO LEARN MORE](#)



SURGICAL PROTECTION SOLUTIONS

SURGICAL PROTECTION

Reducing the transmission of microorganisms between vets and animals in surgical environments is paramount. Our portfolio of sterile surgical gloves minimize the risk of cross-contamination that are instrumental during surgical procedures.



Polyisoprene



Neoprene



GAMMEX® PI Hybrid

Polyisoprene and neoprene blend surgical glove delivering the best of both materials: exceptional comfort and superior durability

- Hybrid technology: Polyisoprene's comfort and neoprene's durability
- Low-allergenic formulation: Reduced risks to a range of allergens*
- Semi-transparent design: Fast and easy breach detection when double gloving

**free from natural rubber latex proteins, as well as Zinc 2-mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG) and cetylpyridinium chloride (CPC)*

External Glove Surface : Smooth with Micro-Textured Finish
Powder Content : Powder-Free
Cuff Style : Beaded with SUREFIT™ Technology

Freedom From Holes : 0.65 AQL
(Inspection Level I)

Ordering Information:

5.5	6	6.5	7
340063055	340063060	340063065	340063070
7.5	8	8.5	9
340063075	340063080	340063085	340063090



[SCAN OR CLICK HERE TO LEARN MORE](#)

50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case



Neoprene



GAMMEX® Non-Latex Sensitive

Neoprene, chemical accelerator-free surgical glove offering unsurpassed combination of sensitivity and durability

- Non-latex, chemical accelerator-free formulation: Optimal allergy protection
- Ultra-thin technology: 30% thinner* for sensitive yet durable gloves
- Smooth surface with micro-textured finish: Secure instrument grip and handling

**when compared with standard Ansell neoprene surgical gloves*

External Glove Surface : Smooth with Micro-Textured Finish
Powder Content : Powder-Free
Cuff Style : Beaded with SUREFIT™ Technology

Freedom From Holes : 0.65 AQL
(Inspection Level I)

Ordering Information:

5.5	6	6.5	7
20277255	20277260	20277265	20277270
7.5	8	8.5	9
20277275	20277280	20277285	20277290



[SCAN OR CLICK HERE TO LEARN MORE](#)

50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case



Natural Rubber Latex



ENCORE® Latex Textured

Latex, powder-free surgical glove with a textured finish for exceptional instrument handling and grip

- Textured finish: Secure grip of instruments under dry or damp conditions
- Optimized thickness: Ideal balance of tactile sensitivity and protection
- Donning aid technology: Easy donning and doffing for dry or damp hands

External Glove Surface : Textured Fingers and Palm
Powder Content : Powder-Free
Cuff Style : Beaded with SUREFIT™ Technology

Freedom From Holes : 0.65 AQL
(Inspection Level I)

Ordering Information:

5.5	6	6.5	7
5785000	5785001	5785002	5785003
7.5	8	8.5	9
5785004	5785005	5785006	5785007



[SCAN OR CLICK HERE TO LEARN MORE](#)

50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case

NRL

Natural Rubber Latex



ENCORE® Latex Micro

Latex, powder-free surgical glove with exceptional dexterity and tactile sensitivity for surgical finesse

- Micro thin: Made 20% thinner* for enhanced precision performance
 - Brown color: Reduces glare of bright surgical lights to prevent eye fatigue
 - Smooth with micro-textured finish: Unique surface offers secure grip
- *when compared with Ansell's standard latex gloves*

External Glove Surface : Smooth with Micro-Textured Finish **Freedom From Holes** : 0.65 AQL
Powder Content : Powder-Free **(Inspection Level GI)**
Cuff Style : Beaded with SUREFIT™ Technology

Ordering Information:

5.5	6	6.5	7
5787000	5787001	5787002	5787003
7.5	8	8.5	9
5787004	5787005	5787006	5787007



[SCAN OR CLICK HERE TO LEARN MORE](#)

50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case

NRL

Natural Rubber Latex



ENCORE® Latex Acclaim

Latex, powder-free surgical glove with a comfortable, natural fit and feel

- Soft yet robust formulation: Optimized comfort and tactile sensitivity
- Ergonomic design: Reduced hand muscle effort during surgery
- Micro-textured external finish: Safe and secure grip of medical instruments

External Glove Surface : Smooth **Freedom From Holes** : 0.65 AQL
Powder Content : Powder-Free **(Inspection Level I)**
Cuff Style : Beaded with SUREFIT™ Technology

Ordering Information:

5.5	6	6.5	7
5795000	5795001	5795002	5795003
7.5	8	8.5	9
5795004	5795005	5795006	5795007



[SCAN OR CLICK HERE TO LEARN MORE](#)

50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case

NRL

Natural Rubber Latex



ENCORE® Perry Style 42™ PF

Latex powder-free surgical gloves, perfectly combining protection and tactile sensitivity

- Optimal glove thickness: A balance of protection and tactile sensitivity
- Anatomic design: Enhanced comfort and reduced hand muscle effort
- Polymer-coated internal surface: Quick and easy donning

External Glove Surface : Smooth with Micro-Textured Finish **Freedom From Holes** : 0.65 AQL
Powder Content : Powder-Free **(Inspection Level I)**
Cuff Style : Beaded

Ordering Information:

5.5	6	6.5	7
5711100PF	5711101PF	5711102PF	5711103PF
7.5	8	8.5	9
5711104PF	5711105PF	5711106PF	5711107PF



[SCAN OR CLICK HERE TO LEARN MORE](#)

50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case



Natural Rubber Latex



ENCORE® SENSITOUCH® PF

Latex, powder-free surgical glove with a smooth finish and strong formulation for everyday protection

- Ideal comfort and fit
- Polymer coated for fast and easy donning
- Features a smooth finish with the addition of a beaded cuff to help prevent cuff roll-down

External Glove Surface : Smooth with Micro-Textured Finish **Freedom From Holes** : 0.65 AQL (Inspection Level I)
Powder Content : Powder-Free
Cuff Style : Beaded

Ordering Information:

5.5	6	6.5	7
7821PF	7822PF	7823PF	7824PF
7.5	8	8.5	9
7825PF	7826PF	7827PF	7828PF



[SCAN OR CLICK HERE TO LEARN MORE](#)

50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case



A photograph of a veterinary team in a clinical setting. Three people wearing blue scrubs, green hairnets, and face masks are performing a procedure on a light-colored dog lying on a stainless steel table. One person is using a large white container to pour liquid over the dog's head. Medical equipment, including a monitor and IV stands, is visible in the background. A red mat is on the floor.

SANDEL® MEDICAL SOLUTIONS

SHARPS SAFETY SOLUTIONS

O.R.'s are areas where risk of sharps injuries are present. Our line-up of safety solutions reduce the risk of sharps injuries while providing the performance vets need.



Z-Friction Drape™ Neutral Zone

Instrument drape for use as a neutral zone to reduce the risk of sharps injuries

- Helps facilitate hands-free transfer of sharps
- Bright orange for high visibility
- Non-magnetic
- Non-slip

Ordering Information:

Z-Friction Drape™ Neutral Zone (Sterile) 10"x16"	1151
--	------



[SCAN OR CLICK HERE TO LEARN MORE](#)



SharpsRink™ Neutral Zone

A non-magnetic, reusable neutral zone designed to support the hands-free transfer of sharps instruments

- Non-magnetic, non-slip design will help hold instruments where you place them
- Ideal for facilities looking to reduce their environmental impact through the implementation of reusable products

Ordering Information:

SharpsRink™ Neutral Zone (Non-Sterile)	1148-NS
--	---------



[SCAN OR CLICK HERE TO LEARN MORE](#)



Z-Tray™ & Stretch-A-Tray Hands-Free Transfer Trays

Disposable instrument trays designed to minimize the risk of sharps injuries through the hands-free transfer of sharps

- Provides a designated area to safely place and retrieve sharp instruments
- Bright orange and highly recognizable
- Adhesive on bottom for the optional use as a non-movable neutral zone
- Stretch-A-Tray expands to accommodate a variety of instrument lengths

Ordering Information:

Z-Tray™ Hands-Free Transfer Tray (Sterile)	1145
Stretch-A-Tray Hands-Free Transfer (Sterile)	1147



[SCAN OR CLICK HERE TO LEARN MORE](#)

STAFF ERGONOMICS

Ergonomics doesn't just stop at the hands. The back and legs of vet staff need to be supported to prevent overexertion and even musculoskeletal disorders



2102



2103 & 2103-H



2101



2104 & 2104-H



2101-H

Z-Slider™ Animal Transfer Sheet

Disposable animal transfer sheet designed to prevent disabling back injuries to staff

- Reduces the risk of physical strain to back, shoulders, neck and arms caused by moving animals
- No weight limit
- Available in a variety of point-of-use dispensers for improved staff compliance
- Disposable and recyclable

Ordering Information:

Z-Slider™ - Single Pack (Non-Sterile)	2101
Z-Slider™ XL (Non-Sterile) 45" x 53"	2151
Z-Slider™ - Roll (Roll Dispenser Box) (Non-Sterile)	2102
Z-Slider™ Wall Dispenser (Large) (Non-Sterile)	2103
Z-Slider™ Wall Dispenser (Small) (Non-Sterile)	2104
Point-of-Use Clear Wall Dispenser Single pack	2101-H
Point-of-Use Clear Wall Holder 12"x20.5"	2103-H
Point-of-Use Clear Wall Holder 11.5"x11"	2104-H



[SCAN OR CLICK HERE TO LEARN MORE](#)



ErgoPlus™ Anti-Fatigue Mats

Designed to reduce the stress and strain on muscles and joints caused from standing for extended procedures

- Bounce-back memory restores energy and reduces strain on the feet, legs and back
- Proven to reduce spinal compression and increase flexibility*
- Premium padding for additional comfort
- Beveled edges help prevent tripping
- Ideal for treatment room, dental lab, radiography room, etc
- New easy-to-clean finish
- Available in various sizes
- Bright orange and highly recognizable

Ordering Information:

ErgoPlus™ Anti-Fatigue Mat 13"x17"	2317-P
ErgoPlus™ Anti-Fatigue Mat 20"x32"	2332-P
ErgoPlus™ Anti-Fatigue Mat 20"x48"	2348-P

* Data on file



[SCAN OR CLICK HERE TO LEARN MORE](#)



1170-B



1170

Ergo-Step™ Stool

Light-weight, stackable step stools

- Rubber top for added comfort
- Anti-skid feet keep stools firmly in place
- Can be stacked
- Interlocking feature creates a safe working platform
- Bright orange and highly visible

Ordering Information:

Ergo-Step™ Stool - 2 Top Stools, 2 Bottom Stools (Each 5.25" x 18" x 14")	1170
Ergo-Step™ Stool - 4 Bottom Stools (Each 5.25" x 18" x 14")	1170-B
Ergo-Step™ Top Stools (Each 5.25" x 18" x 14")	1170-T
Ergo-Step™ Stool Cart - Holds a Maximum of 6 Stools	1170-C



[SCAN OR CLICK HERE TO LEARN MORE](#)

NOTE: Top step stools (1170-T) should only be stacked on bottom step stools (1170-B). 1170-B can be stacked interchangeably as either a top or bottom stool

STAFF ERGONOMICS (CONTINUED)



Sit-Stand™ Stool

A stool designed to relieve pressure and fatigue during long procedures

- Ideal for the O.R. and vet tech station
- Height 20" (min.) to 30" (max.)
- Weight capacity – 330 lbs

Ordering Information:

Sit-Stand Stool (Non-Sterile)

2300



SCAN OR [CLICK HERE](#)
TO LEARN MORE

TABLE STRAPS

Even the most minor of equipment needs to be paid attention to as it can affect the safety outcome of patients.



TSM-SM / TS-10



TS-SS



TS-PD

STAT-STRAP™ Table Straps

Disposable table straps for surgical positioning

- Disposable to minimize the risk of cross-contamination
- Adjustable sizing to fit most animals

Ordering Information:

STAT-STRAP™ Table Strap Standard, Velcro Closure, 2 piece	TS-SM
STAT-STRAP™ Table Strap Long, Velcro Closure, 2 piece	TS-10
STAT-STRAP™ Table Strap Soft, Velcro Closure, 2 piece	TS-SS
STAT-STRAP™ Table Strap Padded, Velcro Closure, 2 piece	TS-PD



SCAN OR [CLICK HERE](#)
TO LEARN MORE

DISPOSABLE SHEETS

When lying down on beds or other surfaces, animals might be at risk of infection from contaminated surfaces. Our antimicrobial disposable sheets helps maintain a clean surgical surface.



Draw Sheets

Surgical table covers and draw sheets for moving sedated animals

- Soft, non-woven, polyester sheets used to move, position or cover animals before, during or after procedures
- Disposable to help prevent cross-contamination
- Standard draw sheet is capable of supporting up to 1,000 pounds of weight in both wet and dry conditions

Ordering Information:

Disposable, Standard Draw Sheet, White (34" x 66")	LFS-66
Disposable, XL Draw Sheet, White (40" x 80")	LFS-4080
Disposable, Dual-Layer Draw Sheet, White (34" x 66")	LFS-66-BF

TO BE UPDATED

FLUID AND TRIP MANAGEMENT

Slips, trips, and falls can occur within surgical environments, especially when cables are strewn around or the surroundings are wet.



Trip-No-More™ Cord Cover

Bright orange cord cover designed to reduce the risk of trips and falls caused by cords and tubes on the floor

- Serves as a “beacon” so healthcare personnel know where cords and tubes are located
- Adhesive strips keep it firmly in place
- Easy to apply, easy to remove
- Disposable, leaves no residue
- Bright orange and highly recognizable

Ordering Information:

Trip-No-More™ Cord Cover (Non-sterile) 8" x 24"	2302
Trip-No-More™ Cord Cover (Non-sterile) 8" x 36"	2303
Trip-No-More™ Cord Cover (Non-sterile) 8" x 125"	2309
Trip-No-More™ Cord Cover Glow in the Dark (Non-sterile) 8" x 24"	2312
Trip-No-More™ Cord Cover Glow in the Dark (Non-sterile) 8" x 36"	2313



[SCAN OR CLICK HERE TO LEARN MORE](#)



DriFloor™ Absorbent Floor Pad

Absorbent floor pads designed to reduce the risk of slips and falls in the operating room or wherever fluids collect on the floor

- Impervious, low-slip backing
- Super absorbent
- Eliminates the use of blankets and towels on the floor
- Helps reduce the risk of cross-contamination
- Available in a variety of sizes

Ordering Information:

Absorbent Floor Pad with Backing (Non-sterile) 16.5"x20"	2419
Absorbent Floor Pad with Backing (Non-sterile) 33"x40"	2440
Absorbent Floor Pad with Backing (Non-sterile) 35"x46"	2446
Absorbent Floor Pad with Backing (Non-sterile) 28"x125"	2425
Absorbent Floor Pad with Backing - Roll (Non-sterile) 33"x100'	2449



[SCAN OR CLICK HERE TO LEARN MORE](#)

INFECTION CONTROL

During surgical procedures, infectious waste fluid needs to be contained quickly. Our infection control solutions makes disposal of blood and other fluids safe.



STAT-SORB™ Solidifier

Fast-acting solidifier for safe disposal of blood and other fluids

- Screw-off cap for easy dispensing
- Can be used for suction canisters, red bags, spills, accidents and more

Ordering Information:

Solidifier for 500cc (16 oz.) of fluid	SS-500
Solidifier for 1200cc (40 oz.) of fluid	SS-1200
Solidifier for 1500cc (50 oz.) of fluid	SS-1500
Solidifier for 2000cc (68 oz.) of fluid	SS-2000
Solidifier for 3000cc (100 oz.) of fluid	SS-3000



[SCAN OR CLICK HERE TO LEARN MORE](#)

SPECIMEN MANAGEMENT

Proper specimen handling is crucial to help reduce the risk of lost and mislabeled specimens.



Specimen Handling System

With Specimen Handling System products, specimens are easily recognized to help prevent accidental disposal or other mismanagement.

- Leak-proof, puncture resistant specimen cups
- Moisture-proof orange labels
- Bright orange, biohazard labeled bags to protect workers from infectious specimens
- Easy-to-handle heart shaped specimen bowl
- Specimen Zone aids identification of tissue samples
- All items sold separately

Ordering Information:

Orange Specimen Zone (4 Count) (Sterile)	3450
Orange Specimen Labels on roll (Non-Sterile) 2"x4"	3455
Orange Specimen Wrap (Sterile) 23 1/2" x 16 1/2"	3460
Orange Specimen Wrap (Non-Sterile) 23 1/2" x 16 1/2"	3460-NS
120cc Cup with Orange Lid (Sterile)	3470
Heart Shaped Orange Specimen Bowl (Sterile)	3475
Heart Shaped Specimen Bowl (Non-Sterile)	3475-NS
Orange Specimen Bag Large (Non-Sterile) 33 3/8" x 39 1/4"	3480

TO BE UPDATED

SURGICAL SITE MARKING

Preventing errors in the O.R. is crucial. Marking the correct surgical site will reduce the risk of wrong-site surgery.



PREP-RESIST™ Skin Markers

Offers strong visibility following skin prep but are easily removable in the time following surgery

- Resistant to chlorhexidine gluconate and isopropyl alcohol skin prep solutions
- Gentian violet-free
- Available in standard length, petite, and 4-in-1 dual-tip/ dual-ink marker with chiseled tip
- Triangular “no-roll” barrel
- 4-in-1 marker comes with TIME OUT® sleeve

Ordering Information:

Petite Prep-Resistant “No Roll” Skin Marker (Non-Sterile)	1015-NNS
Prep-Resistant “No Roll” Skin Marker (Sterile)	1029
Prep-Resistant “No Roll” Skin Marker (Non-Sterile)	1029-NNS
Dual Tip/ Dual Ink Marker (Permanent and Prep-Resistant Skin Ink) with Chiseled Tip, TIME OUT® Sleeve (Sterile)	1061
Dual Tip/ Dual Ink Marker (Permanent and Prep-Resistant Skin Ink) with Chiseled Tip, TIME OUT® Sleeve (Non-Sterile)	1061-NNS



[SCAN OR CLICK HERE TO LEARN MORE](#)



TIME OUT® Skin Markers & Permanent Markers

For marking the correct surgical site to help reduce the risk of wrong-site surgery

- Triangular barrel prevents rolling
- Available in single and dual-tip
- Gentian violet-free ink
- SANDEL® 4-in-1 Marker features a chiseled tip to allow the user to make lines with multiple thickness options
- Permanent marker for medication labeling and general use

Ordering Information:

“No-Roll” Permanent Marker, Black Ink (Non-Sterile)	1000-NNS
Petite “No-Roll” Permanent Marker, Black Ink (4” in Length) (Non-Sterile)	1002-NNS
Petite “No Roll” Skin Marker (4” in Length) (Non-Sterile)	1011-NNS
Petite “No Roll” Skin Marker (4” in Length) (Non-Sterile)	1011-NS1
“No Roll” Skin Marker (Sterile)	1019
“No Roll” Skin Marker with Flexible Ruler (Sterile)	1019-R
“No Roll” Skin Marker (Non-Sterile)	1019-NNS
“No Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink), TIME OUT® Sleeve (Sterile)	1031
“No Roll” Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels, TIME OUT® Sleeve & Flexible Ruler (Sterile)	1031-LR
“No Roll” Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels, TIME OUT® Sleeve & Flexible Ruler (Non-Sterile)	1031-LRNS
“No Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink), TIME OUT® Sleeve (Non-Sterile)	1031-NNS
“No Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink), Skin Marker with Chiseled Tip, TIME OUT® Sleeve (Sterile)	1041
“No Roll” Dual Tip/Dual Ink (Permanent and Skin Ink) Skin Marker with Chiseled Tip, TIME OUT® Sleeve and Flexible Ruler (Sterile)	1041-R
“No Roll” Dual Tip/Dual Tip (Permanent & Skin Ink) Marker with 8 labels & TIME OUT® Sleeve (Sterile)	1041-L
“No Roll” Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels, TIME OUT® Sleeve & Flexible Ruler (Sterile)	1041-LR
“No Roll” Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels, TIME OUT® Sleeve & Flexible Ruler (Non-Sterile)	1041-LRNS
“No Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink), Skin Marker with Chiseled Tip, TIME OUT® Sleeve (Non-Sterile)	1041-NNS
“No-Roll” Dual Tip/Dual Ink Marker (Permanent & Skin Ink) Skin Marker with Chiseled Tip (Sterile)	1041-NS1



[SCAN OR CLICK HERE TO LEARN MORE](#)

MOPS

Keep the O.R. clean after surgical procedures. Prevent cross contamination with our selection of disposable mops



M16-R



M16-RF



STM-MD



STM-MD-RF

Rayon Mop

A variety of mops for cleaning any hard floor surfaces

- Traditional, string mop heads
- Disposable to help prevent cross-contamination

Ordering Information:

Rayan Mop Head 12oz. (#16)	M16-R
Rayan Mop Head 12oz. (#16) with Fan Tail	M16-RF
Rayan Mop Head 14oz.	M14OZ-R
Rayan Mop Head 14oz. with Fan Tail	M14OZ-RF
Rayan Mop Head 16oz	M16OZ-R
Rayan Mop Head 16oz. with Fan Tail	M16OZ-RF

TO BE UPDATED

STAT-MOP™

Specially designed, disposable mop head to help prevent cross-contamination

- Disposable to help prevent cross-contamination
- More absorbent than traditional cotton, rayon and microfiber mops*
- Low-linting
- Made from 100% post-industrial material*

Ordering Information:

Standard mop head	STM -MD
Standard mop head with fan tail	STM-MD-RF

* Data on file



[SCAN OR CLICK HERE TO LEARN MORE](#)

Microfiber Mops

- Disposable microfiber mops use less chemicals and water than traditional mops
- Designed to attract and trap dust, dirt and bacteria
- Highly absorbent and non-linting

Ordering Information:

Microfiber Flat Mop, Blue	MMBT479RM-18
Microfiber & Cotton Flat Mop, White	MMW-616-P

TO BE UPDATED



MMBT479RM-18



MMW-616-P

WASTE BAGS

Dispose waste materials and biohazard waste in containers that are easily identifiable and are tough to prevent rips and tears.



STAT-BAG™ Heavy-duty Waste Bags

Disposable bags for trash, biohazard waste and laundry

- Specially formulated resin to help prevent tears
- Accommodate irregular shaped items
- Available with drawcord or drawtape for easy, secure closure
- Additional sizes and colors available

Ordering Information:

Clear Bag with Drawtape; Large	BCD-LG
Clear Bag with Drawcord; Large	BCC-LG
Clear Kick Bucket Bag	KBC
Red Bag with Drawtape; Large with Biohazard symbol	BRD-LG
Red Bag with Drawcord; Large with Biohazard symbol	BRC-LG
Red Kick Bucket Bag; with Biohazard symbol	KBR
Blue Bag with Drawtape; Large	BBD-LG
Blue Bag with Drawcord; Large	BBC-LG
Black Bag with Drawcord; Large	BBRC-LG



[SCAN OR CLICK HERE TO LEARN MORE](#)



Z-Suture Saddle Bag™ Waste Bags

Clear bag for disposal of sutures or other items

- Clear for easy identification of contents
- Can be draped over table
- Perforated so bags can be placed in separate locations

Ordering Information:

20" for use on the Mayo Stand or Back Table (Sterile)	1130
---	------



[SCAN OR CLICK HERE TO LEARN MORE](#)



Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: info@ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2023 Ansell Limited. All Rights Reserved.