



Isolator Sleeve/Glove Assembly Unit

BioClean Isolator Sleeve/Glove Assembly Unit is used for fixing glove to isolator sleeve.

- PVC C Ring & Silicone Ring material
- Double bagged in EasyTear Packaging
- Enables fixing of glove to isolator sleeve





BBCO 100

Isolator Glove Assembly Unit

TECHNICAL DATA SHEET

Product Reference	BBCO 100
Material	PVC, Silicone
Color	Grey, Red
Packaging Overview	One assembly unit consisting of 1 x C ring and 1 x O ring per inner PE bag; one inner bag per outer bag.
Storage Instructions	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Country Of Origin	Vietnam
Available Sizes	11
Cleanroom Class	Class 10/ISO 4
Shelf Life	Five (5) years from date of manufacture.



BBCO 100

Isolator Glove Assembly Unit

ORDERING INFORMATION

	PORT SIZE	
	SIZE	11
BBCO 100	REORDER NO.	BBCO-100

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2022 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

