

BIOCLEAN NOTEBOOKS

Cleanroom Notebook
BNBC

About

BioClean Spiral Bound Cleanroom Notebooks are the ideal choice for use in cleanrooms. The notebooks are constructed from latex-free, low-linting, heat resistant, recyclable cleanroom paper, and the spiral bound translucent covers are constructed from low-linting polypropylene. BioClean Notebooks are available in A4, A5, and A6 sizes (50 pages per book).

Specifications

WEIGHT:

80gsm \pm 7g

COMPATIBILITY:

ISO Class 4

CONSTRUCTION:

Spiral bound

COLOUR:

Blue Paper & Blue Cover

Features

- Polypropylene spiral & covers
- Latex-free, low linting & heat resistant
- Packed in an ISO Class 4 cleanroom
- High tear & tensile strength
- All sizes triple bagged
- Sterilizable & autoclavable

LATEX-FREE

LOW LINTING & HEAT RESISTANT

HIGH TEAR & TENSILE STRENGTH

STERILIZABLE & AUTOCLAVABLE

Quality Standards

- Manufactured in a facility holding ISO 9001:2015 certification

Shelf Life & Storage

Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

TO ORDER

Re-order code

BNBC-A4-L
BNBC-A5-L
BNBC-A6-L

Size

A4 Lined
A5 Lined*
A6 Lined*

Packing

A4 - Two books triple bagged, with each book individually double bagged. 25 packs (50 notebooks) per carton
A5 - Two books triple bagged, with each book individually double bagged. 50 packs (100 notebooks) per carton
A6 - Two books triple bagged, with each book individually double bagged. 100 packs (200 notebooks) per carton
Minimum Order Quantity - 1 carton

*Size not available in USA/Canada

For additional information visit us at www.bioclean.com, or contact us at;

Nitritex Ltd

Ground Floor, 15 Kings Court
Willie Snaith Road
Newmarket, Suffolk
CB8 7SG
United Kingdom
Tel: +44 1638 663338

Ansell US

111 Wood Avenue South
Suite 210
Iselin, NJ 08830
USA
Tel: +1 800 800 0444

Nitritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20
Seksyen U1, Hicom Glenmarie Industrial Park
40150 Shah Alam, Selangor
Malaysia
Tel: +603 5569 3857/3859

BNBC/2107/PDS23

Ansell

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2021 Ansell Limited. All Rights Reserved.