

# BIOCLEAN SOFTFLOW

Face Veil with Headloop

## BFV03

### About

Made from soft non-woven material, the BioClean Softflow Cleanroom Face Veil offers comfort for the wearer and features ultrasonically sealed seams and a fully enclosed malleable nose-band.

### Materials

**VEIL:**

Hydroentangled polyester

**FASTENINGS:**

Tubular polyester headloops

**NOSE-BAND:**

Plastic coated aluminium

**BINDING:**

White hydroentangled polyester

### Features

- Ultrasonically sealed seams
- Fully enclosed malleable nose-band
- ISO Class 4 compatible

ULTRASONICALLY SEALED SEAMS

LATEX-FREE

SOFT & COMFORTABLE



Mr. & Mrs Phil Trayshon© represent  
BioClean's range of GMP compatible  
cleanroom facemasks

## Quality Standards

- Manufactured in a facility holding ISO 9001:2015
- Processed in an NEBB certified ISO Class 4 cleanroom environment
- Sampling and inspection in accordance with ISO 2859-3:1994

## Shelf Life & Storage

Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

## Physical Properties

COLOUR	White
STYLE	Headloops
LENGTH (± 5mm)	380mm
WIDTH (± 5mm)	280mm

## Performance Information

COMPATIBILITY	ISO Class 4
PARTICLE FILTRATION EFFICIENCY (PFE)@ 0.5 MICRON	26%
BACTERIAL FILTRATION EFFICIENCY (BFE) @ 3.0 MICRON	>2.6%
DIFFERENTIAL PRESSURE ( $\Delta P$ ), mm H <sub>2</sub> O/cm <sup>2</sup>	<2.0

## TO ORDER

Re-order code

Size

BFV03

### Packing

100 pieces per sealed inner PE bag; one inner bag per sealed outer PE bag; 20 outer bags per lined carton (2000 pieces)

For additional information visit us at [www.bioclean.com](http://www.bioclean.com), or call us at +44 1638 663338

### Nitritex Ltd

Minton Enterprise Park  
Oaks Drive  
Newmarket, Suffolk  
CB8 7YY  
United Kingdom

BFV03/1804/PDS22

**Ansell**

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of facemasks/veils for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2018 Ansell Limited. All Rights Reserved.