

BIOCLEAN DB

Sterile Pouch-style Facemask with Neck Guard

BDBS-G

About

The BioClean DB Pouch-style Sterile Facemask features a large breathing chamber for increased wearer comfort. Made from cleanroom compatible materials the BioClean DB facemask features an integrated neck guard to provide additional coverage reducing the risk of cross-contamination.

Materials

INNER FACING LAYER:

Non-woven spunbonded polyester
(hygroscopic)

FILTER LAYER:

Meltblown polyester

OUTER FACING LAYER:

Non-woven spunbonded polyester
(hydrophobic)

FASTENINGS:

Tubular knitted polyurethane Spandex
yarn headloops

NECK GUARD:

Non-woven spunbonded polyester

NOSE-BAND:

Plastic coated steel

ISO CLASS 4 COMPATIBLE

LARGE BREATHING CHAMBER

INTEGRAL NECK GUARD

Features

- High bacterial & particle filtration efficiency
- Extra long neck guard
- Large breathing chamber
- Ultrasonically sealed edges
- Fully enclosed malleable nose-band

Quality Standards

- Manufactured in a facility holding ISO 9001:2015
- Processed in an NEBB certified ISO Class 4 cleanroom environment
- Sampling and inspection in accordance with ISO 2859-3:1992

Sterilization

Method: Gamma irradiation
Minimum Dose: 25kGy
Sterility Assurance Level: 10^{-6}

Shelf Life & Storage

Five (5) years from date of manufacture.
Store in a dry, cool place ($<40^{\circ}\text{C}$) away from direct sunlight and fluorescent light.

Physical Properties

COLOUR	White
STYLE	Headloops
LENGTH ($\pm 5\text{mm}$)	Facemask: 301mm, Neck Guard: 245mm
WIDTH ($\pm 5\text{mm}$)	Facemask: 100mm, Neck Guard: 210mm

Performance Information

COMPATIBILITY	ISO Class 4 & EU GMP Grade A
PARTICLE FILTRATION EFFICIENCY (PFE)@ 0.3 MICRON	98.2%
BACTERIAL FILTRATION EFFICIENCY (BFE) @ 3.0 MICRON	97.9%
VIRAL FILTRATION EFFICIENCY (VFE)	98.5%
DIFFERENTIAL PRESSURE (ΔP), mm H ₂ O/cm ²	3.57

TO ORDER

Re-order code

BDBS-G

Packing

One piece per sealed inner PE bag; 20 inner bags per sealed outer PE bag; 10 outer bags per lined carton (200 pieces)

For additional information visit us at www.bioclean.com, or contact us at;

Nitritex Ltd
Ground Floor, 15 Kings
Court
Willie Snaith Road,
Newmarket,
Suffolk, CB8 7SG,
United Kingdom
Tel: +44 1638 663338
BDBS-G/1903/PDS34

Ansell US
111 Wood Avenue South
Suite 210
Iselin, NJ 08830
USA
Tel: +1 800 800 0444

Nitritex (M) Sdn.Bhd
2 Jalan Jurunilai U1/20
Seksyen U1, Hicom Glenmarie Industrial Park
40150 Shah Alam, Selangor
Malaysia
Tel: +603 5569 3857/3859

Ansell

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of facemasks/veils for a specific application.