

BIOCLEAN MEA

Sterile Looped Facemask
MEA210-1

About

Providing high particle and bacterial filtration efficiency, the BioClean™ MEA Cleanroom Sterile Looped Facemask is manufactured from cleanroom compatible materials to reduce contamination into the controlled environment and features loops and clip connector to allow for quick and secure fastening at back of head.

Materials

INNER FACING LAYER:

Non-woven polypropylene/polyethylene (hydroscopic)

FILTER LAYER:

Meltblown polyester

OUTER FACING LAYER:

Non-woven spunbonded polyester (hydrophobic)

FASTENINGS:

Non-latex polyurethane loops with blue plastic clip fastener

NOSE-BAND:

Plastic coated steel

Features

- High bacterial, viral & particle efficiency filtration
- Fully enclosed malleable nose-band
- Ultrasonically sealed edges
- Looped with connector for secure fastening

HIGH BFE, VFE & PFE

ISO CLASS 4 COMPATIBLE

LATEX-FREE

Ansell

Quality Standards

- Manufactured in a facility holding ISO 9001:2015
- Processed in an NEBB certified ISO Class 4 cleanroom environment
- Sampling and inspection in accordance with ISO 2859-3:1994
- Tested against standard ASTM F2101 for Bacterial Filtration Efficiency (BFE)

Physical Properties

COLOUR	White
STYLE	Headloops
LENGTH (\pm 5mm)	210mm
WIDTH (\pm 5mm)	95mm

Performance Information

COMPATIBILITY	ISO Class 4 & EU GMP Grade A
PARTICLE FILTRATION EFFICIENCY (PFE)@ 0.3 MICRON	>97%
BACTERIAL FILTRATION EFFICIENCY (BFE) @ 3.0 MICRON	>96.5%
VIRAL FILTRATION EFFICIENCY (VFE)	97.9%
DIFFERENTIAL PRESSURE (ΔP), mm H ₂ O/cm ²	3.2

TO ORDER

Re-order code

MEA210-1

Packing

One piece (with blue clip) per sealed inner PE bag; 50 inner bags per sealed outer PE bag; 12 outer bags per lined carton (600 pieces)

For additional information visit us at www.bioclean.com, or contact us at;

Nitritex Ltd
Ground Floor, 15 Kings
Court
Willie Snaith Road,
Newmarket,
Suffolk, CB8 7SG
UK
Tel: +44 1638 663338
MEA210-1/1903/PDS31

Ansell US
111 Wood Avenue South
Suite 210
Iselin, NJ 08830
USA
Tel: +1 800 800 0444

Nitritex (M) Sdn.Bhd
2 Jalan Jurunilai U1/20
Seksyen U1, Hicom Glenmarie Industrial Park
40150 Shah Alam, Selangor
Malaysia
Tel: +603 5569 3857/3859

Ansell

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of facemasks/veils for a specific application.