

BIOCLEAN CLEARVIEW

Sterile Visor Facemask
VFM

About

Providing good particle and bacterial filtration efficiency, the latex-free BioClean VFM Sterile Visor Facemask combines comfort and breathability. The optically clear, distortion-free anti-fog visor provides additional protection reducing the risk of contamination entering the controlled environment.

Materials

INNER FACING LAYER:

Non-woven polypropylene/polyethylene (hydrophilic)

FILTER LAYER:

Meltblown polyester

OUTER FACING LAYER:

Non-woven spunbonded polyester (hydrophobic)

FASTENINGS:

Ties: Non-woven Polypropylene

Loops: Polyurethane (with blue plastic clip fastener)

NOSE-BAND:

Plastic coated steel

Features

- Anti-fog coating
- Distortion free visor
- Latex-free
- High bacterial & particle filtration efficiency
- Looped or tie-on versions
- ISO Class 4 & EU GMP Grade A compatible

OPTICALLY CLEAR ANTI-FOG VISOR

GOOD BFE & PFE

LATEX-FREE



Ansell

Quality Standards

- Manufactured in a facility holding ISO 9001:2015 and PPER 2016/425 Module D
- Processed in an NEBB certified ISO Class 4 cleanroom environment
- Sampling and inspection in accordance with ISO 2859-3:1993

Sterilization

Method: Gamma irradiation
Minimum Dose: 25kGy
Sterility Assurance Level: 10^{-6}

Shelf Life & Storage

Five (5) years from date of manufacture.
Store in a dry, cool place ($<40^{\circ}\text{C}$) away from direct sunlight and fluorescent light.

Physical Properties

COLOUR	White
STYLE	Tie-on or Headloops
LENGTH ($\pm 5\text{mm}$)	210mm
WIDTH ($\pm 5\text{mm}$)	95mm

Performance Information

COMPATIBILITY	ISO Class 4 & EU GMP Grade A
PARTICLE FILTRATION EFFICIENCY (PFE)@ 0.3 MICRON	$>96\%$
BACTERIAL FILTRATION EFFICIENCY (BFE) @ 3.0 MICRON	$>94\%$
DIFFERENTIAL PRESSURE (ΔP), mm H ₂ O/cm ²	3.2

TO ORDER

Re-order code

VFM210-L
VFM210-T-SLOT

Packing

VFM210-L (Looped with blue clip fastener);
One piece per sealed inner PE bag; 25 inner bags per sealed outer PE bag; four outer bags per lined carton (100 pieces)
VFM210-T-SLOT (Tie-on);
One piece per sealed inner PE bag; 25 inner bags per sealed outer PE bag; four outer bags per lined carton (100 pieces)

For additional information visit us at www.bioclean.com, or contact us at;

Nitritex Ltd
Ground Floor, 15 Kings
Court
Willie Snaith Road,
Newmarket,
Suffolk, CB8 7SG,
United Kingdom
Tel: +44 1638 663338
VFM/1811/PDS30

Ansell US
111 Wood Avenue South
Suite 210
Iselin, NJ 08830
USA
Tel: +1 800 800 0444

Nitritex (M) Sdn.Bhd
2 Jalan Jurunilai U1/20
Seksyen U1, Hicom Glenmarie Industrial Park
40150 Shah Alam, Selangor
Malaysia
Tel: +603 5569 3857/3859

Ansell

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of facemasks/veils for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2019 Ansell Limited. All Rights Reserved.