

About

Constructed from super-soft, lightweight silicone rubber to provide user comfort and enable prolonged use, BioClean Clearview Autoclavable Cleanroom Goggles feature an indirect ventilation system to reduce the risk of contamination entering the controlled environment, and have an anti-fog polycarbonate lens.

Specifications

MATERIAL:

Silicone rubber frame, polycarbonate lens, latex-free silicone head band with polypropylene hooks

COLOUR:

Translucent white frame and strap

CLASSIFICATION:

No Classification

COMPATIBILITY:

ISO Class 4 & EU GMP Grade A

Features

- Anti-fog lens
- Super-soft frame
- Indirect ventilation system
- Non-linting latex-free head band
- Tested to withstand 50 autoclave cycles of 30 minutes duration at 121°C/250°F (under laboratory conditions)
- Anti-fog performance remaining for up to 25 cycles with no signs of degradation (under laboratory conditions)

BIOCLEAN CLEARVIEW Autoclavable Goggles BCAH

ANTI-FOG LENS

SUPER-SOFT FRAME

INDIRECT VENTILATION SYSTEM



Ansell





Quality Standards

- Manufactured in a facility holding ISO 9001:2015
- Processed in an NEBB certified ISO Class 4 cleanroom environment

STERILIZATION:

Supplied non-sterile. To be autoclaved at 121°C/250°F for 30 minutes.

Shelf Life & Storage

Five (5) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

Physical Properties

Total Weight (± 2g)	106g
Lens Thickness (± 0.3mm)	1 mm
Bridge Width (± 0.5mm)	50mm
Bridge Height (±0.5mm)	20mm
Overall Width (hinge to hinge)(±3mm)	180mm
Overall Height (± 3mm)	92mm

Explanation of Markings

REQUIRED MARKINGS	MARKINGS ON FRAME	MEANING
Manufacturer's Identification Mark	BioClean	BioClean™





TO ORDER

Re-order code Size

BCAH Goggle

BCAH-SL Replacement Lens
BCAH-SS Replacement Strap
BCAH-ACC1 Accessory Set

Packing

Goggle (BCAH):

One piece per sealed PE bag (to be removed prior to autoclaving); 12 pieces per box

Replacement Lens (BCAH-SL):

One piece per sealed PE bag; 10 bags per box (10 pieces)

Replacement Strap (BCAH-SS):

25 pieces per sealed PE bag; two bags per box (50 pieces)

Accessory Set (BCAH-ACC1):

One set of accessories per bag; 100 sets per outer bag; 100-300 sets per carton.

For additional information visit us at www.bioclean.com, or contact us at;

Nitritex Ltd

Ground Floor, 15 Kings Court Willie Snaith Road Newmarket, Suffolk CB8 7SG United Kingdom

Tel: +44 1638 663338

BCAH/1907/PDS43

Ansell US

111 Wood Avenue South Suite 210 Iselin, NJ 08830 USA

Tel: +1 800 800 0444

Nitritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20 Seksyen U1, Hicom Glenmarie Industrial Park 40150 Shah Alam, Selangor Malaysia

Tel: +603 5569 3857/3859



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of goggles for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2020 Ansell Limited. All Rights Reserved.