

# BIOCLEAN DUAL

Disposable Sterile Overshoes  
**S-BDOS**

## About

BioClean Dual™ Disposable Sterile Overshoes are practical and durable, featuring a heavy-duty textured cast polyethylene slip-resistant sole for a secure footing.

## Specifications

### MATERIAL:

Top: Spunbonded non-woven polypropylene fabric, Sole: Embossed cast polyethylene film, Elastic: Latex-free

### HEIGHT:

152mm

### LENGTH:

400mm

### COLOUR:

White

### STERILIZATION:

Gamma Irradiation

## Features

- Non-woven spunbonded polypropylene
- Heavy duty slip-resistant sole
- Practical & durable
- Latex-free elastic

**PRACTICAL & DURABLE**

**ISO CLASS 4 COMPATIBLE**

**SLIP-RESISTANT SOLE**



Billy Bioclean©



## Quality Standards

---

- Manufactured in a facility holding ISO 9001:2008 certification
- Processed in an NEBB certified ISO Class 4 environment



## Sterilization

---

Method: Gamma irradiation

Minimum Dose: 25kGy

Sterility Assurance Level:  $10^{-6}$

## Shelf Life & Storage

---

Five (5) years from date of manufacture.  
Store in a dry, cool place ( $<40^{\circ}\text{C}$ ) away  
from direct sunlight and fluorescent  
light.

## TO ORDER

### Re-order code

S-BDOS

### Size

16"

---

### Packing

One pair per sealed inner PE bag; 15 inner bags per sealed outer PE bag; 10 outer bags per lined carton (150 pairs)

For additional information visit us at [www.bioclean.com](http://www.bioclean.com), or call us at +44 1638 663338

#### Nitritex Ltd

Ground Floor, 15 Kings Court  
Willie Snaith Road,  
Newmarket,  
Suffolk, CB8 7SG,  
United Kingdom

S-BDOS/1807/PDS11



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2018 Ansell Limited. All Rights Reserved.