

# BIOCLEAN SUPATAC Tacky Mats STB

## About

BioClean Supatac™ Tacky Mats are an effective particle control solution, manufactured from low density polyethylene coated with an anti-bacterial water based acrylic adhesive.

## Specifications

### MATERIAL:

Low density polyethylene with anti-bacterial finish

### COMPATIBILITY:

ISO Class 5

### CONSTRUCTION:

Upper layers: PE film coated with acrylic adhesive

Lower layers: LPDE backing

Tensile strength: 170kg/cm<sup>2</sup>(min)

Upper layer adhesive force: 180g/25mm  
±40g

Lower layer adhesive force: 600g/25mm  
±100g

Thickness: 0.045mm ± 0.005mm

Elongation: 250% Min

### COLOUR:

Blue

30 NUMBERED LAYERS

ANTI-BACTERIAL

VARIETY OF SIZES

## Features

- 30 numbered layers
- Anti-bacterial water based acrylic adhesive
- Available in variety of sizes
- Also available in white - see product code STW datasheet for details





## Quality Standards

---

- Manufactured in a facility holding ISO 9001:2008 certifications
- Processed in an NEBB certified ISO Class 5 environment



## Shelf Life & Storage

---

Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

## TO ORDER

Re-order code	Size
STB1845	18"x45" Blue
STB1836*	18"x36" Blue
STB2436*	24"x36" Blue
STB2445*	24"x45" Blue
STB2645	26"x45" Blue
STB3636*	36"x36" Blue
STB3645	36"x45" Blue
STB3660*	36"x60" Blue
STB3672*	36"x72" Blue
STB4048*	40"x48" Blue

---

### Packing

Four mats (each with 30 numbered layers) per lined carton

\* Size only available in USA/Canada

For additional information visit us at [www.bioclean.com](http://www.bioclean.com), or call us at +44 1638 663338

**Nitritex Ltd**  
Ground Floor, 15 Kings Court  
Willie Snaith Road,  
Newmarket,  
Suffolk, CB8 7SG,  
United Kingdom  
STB/1804/PDS12



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2018 Ansell Limited. All Rights Reserved.