

### About

BioClean Pharma Covers provide product protection from environmental impurities during storage, processing and transportation with the latex-free elastic providing a secure seal. They minimize the need for autoclave tape during wrapping and save time when prepping glassware, tubing etc for autoclaving. Easy removal minimizes particulate generation and facilitates aseptic technique. Available in a range of diameters and bespoke sizing available on request.

### Specifications

**MATERIAL:**

Polyethylene, IPA resistant coloured thread and high quality elastic

**COMPATIBILITY:**

ISO Class 4

**COLOUR:**

White

### Features

- Breathable microbial barrier
- Lint-free, moisture and puncture resistant
- IPA resistant coloured thread
- Form fitted
- Fast drying time
- Easy to apply and remove
- Autoclavable at 121°C to 127°C

LINT-FREE

MOISTURE AND PUNCTURE RESISTANT

LATEX-FREE ELASTIC

AUTOCLAVABLE



## Quality Standards

- Processed in an NEBB certified ISO Class 4 environment
- Manufactured in a facility holding ISO 9001:2015 and PPER 2016/425 Module D

## Shelf Life & Storage

Three (3) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

## Test Results

Typical Helmke Drum Test Results			
size(μm)	Average counts per minute		
	BPC-008	BPC-041	BPC-061
≥0.3	104	184	223
≥0.5	93	155	208
≥5.0	36	49	66

Endotoxin Test Results	
Method, USP <85> (Kinetic Turbidimetric Method)	<20 EU/Sample
Bacterial Filtration Efficiency (BFE)	
Method, ASTM F2101	98.8%

## TO ORDER

Re-order code	Size
BPC-005	5cm (2.0") Ø
BPC-008	8cm (3.1") Ø
BPC-010	10cm (3.9") Ø
BPC-013	13cm (5.1") Ø
BPC-018	18cm (7.0") Ø
BPC-023	23cm (9.0") Ø
BPC-030*	30cm (11.8") Ø
BPC-041	41cm (16.1") Ø
BPC-051	51cm (20.0") Ø
BPC-061	61cm (24.0") Ø
BPC-076	76cm (29.9") Ø
BPC-091	91cm (35.8") Ø
BPC-122*	122cm (48.0") Ø
BPC-152*	152cm (59.8") Ø
BPC-183*	183cm (72.0") Ø
BPC-229*	229cm (90.0") Ø

### Packing

BPC-005/BPC-008/BPC-010/BPC-013/BPC-018/BPC-023/BPC-030\*/BPC-041/BPC-051:  
50 pieces per sealed inner PE bag; 10 inner bags per sealed outer PE bag; One outer bag per lined carton (500 pieces)

BPC-061/BPC-076/BPC-091/BPC-122\*:  
10 pieces per sealed inner PE bag; 10 inner bags per sealed outer PE bag; One outer bag per lined carton (100 pieces)

BPC-152\*/BPC-183\*/BPC-229\*:  
Two pieces per sealed inner PE bag; 10 inner bags per sealed outer PE bag; One outer bag per lined carton (20 pieces)

\* Sizes not available in USA/Canada

For additional information visit us at [www.bioclean.com](http://www.bioclean.com), or call us at +60 3 5569 3857/3859

#### Niritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20  
Seksyen U1, Hicom Glenmarie Industrial Park  
40150 Shah Alam, Selangor  
Malaysia

BPC/1906/PDS41/APAC

The Ansell logo consists of the word "Ansell" in a bold, blue, sans-serif font. A green swoosh underline is positioned beneath the letters "n", "s", and "e".

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2019 Ansell Limited. All Rights Reserved.