

# SANDEL® Safety Scalpels

S H A R P S   S A F E T Y

## Disposable weighted and unweighted safety scalpels to help reduce the risk of sharps injuries

Studies show the use of safety engineered devices such as safety scalpels can significantly reduce the rate of sharps injuries.<sup>1</sup> SANDEL Safety Scalpels have a retractable safety sheath that covers the blade while the scalpel is being passed. Available weighted and unweighted, SANDEL Safety Scalpels feature a stainless steel Feather blade along with a variety of blade and handle sizes.

### KEY FEATURES & BENEFITS

- Helps meet OSHA Bloodborne Pathogens Standard and AORN sharps safety guidelines
- Weighted safety scalpels are weighted to the feel of a traditional scalpel handle
- Safety shield retracts completely to provide full visibility to the surgical site
- Disposable
- Easy one-hand activation and can be used by a right or left-handed user
- Easy-grip textured surface
- Weighted safety scalpel available with #3 handle with #10, #11, or #15 blade and #4 handle with #20 blade
- Unweighted safety scalpel available with #3 handle with #10, #11 or #15 blade



SANDEL Weighted Safety Scalpels™



SANDEL Non-Weighted Safety Scalpels

### SANDEL Anticipate the Unexpected

Developed in collaboration with clinical professionals, SANDEL products are designed to address key high-risk areas in the healthcare environment. Our signature bright orange products help reduce the risk of preventable injuries, errors and infections. SANDEL products equip caregivers with solutions to do their job safely while improving patient outcomes.

<sup>1</sup>Tosini W, Ciotti C, Goyer F et al., Needlestick injury rates according to different types of safety – engineered devices: results of a French multicenter study. *Infect Control Hosp Epidemiol* 2010; 31(4): 402-407.

# SANDEL® Safety Scalpels

S H A R P S   S A F E T Y

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2210-N	Weighted Safety Scalpel	#3 Handle with #10 Blade, Orange Shield and TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2211-N	Weighted Safety Scalpel	#3 Handle with #11 Blade, Blue Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2215-N	Weighted Safety Scalpel	#3 Handle with #15 Blade, Yellow Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2220-N	Weighted Safety Scalpel	#4 Handle with #20 Blade, Orange Shield and TIME OUT Sleeve (Sterile)	12/Bx 4 Bx/Cs	48
2230-N	Non-Weighted Safety Scalpel	#3 Handle with #10 Blade, Orange Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2231-N	Non-Weighted Safety Scalpel	#3 Handle with #11 Blade, Blue Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96
2235-N	Non-Weighted Safety Scalpel	#3 Handle with #15 Blade, Yellow Shield with TIME OUT Sleeve (Sterile)	12/Bx 8 Bx/Cs	96

✦ **Contact your Ansell representative for ordering or more information.**

**North America**  
US Tel: 866-764-3327  
818-534-2500  
Email: sandel-cs@ansell.com

**Latin America & Caribbean**  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
Tel: +32 (0)2 528 74 00  
Email: info@ansell.eu

**Asia Pacific**  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

**Australia & New Zealand**  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com