

# SANDEL<sup>®</sup> DISARM-IT-ALL<sup>™</sup>



SHARPS SAFETY

## The only scalpel blade disarmer that safely removes traditional scalpel blades, mini scalpel blades and syringe needles

The SANDEL Disarm-It-All is a convenient and versatile scalpel blade disarmer that also acts as a sharps counting and disposal box.

### KEY FEATURES & BENEFITS

- Available in magnet/magnet or magnet/foam
- Provides a safe resting place for sharps in between uses
- Two separate, removable halves with adhesive tape to provide stability on any surface
- Closure tab slides easily into place to lock
- Bright orange and highly recognizable



### SANDEL<sup>®</sup> Anticipate the Unexpected

Developed in collaboration with clinical professionals, SANDEL products are designed to address key high-risk areas in the healthcare environment. Our signature bright orange products help reduce the risk of preventable injuries, errors and infections. SANDEL products equip caregivers with solutions to do their job safely while improving patient outcomes.



# SANDEL<sup>®</sup> DISARM-IT-ALL<sup>™</sup>

S H A R P S   S A F E T Y

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2243	Disarm-It-All (Magnet & Foam)	Foam Needle Counter (40 Count), Scalpel & Specialty Blade Remover, Syringe Needle Remover, Disposal Box	12/Bx 4 Bx/Cs	48
2244	Disarm-It-All (Magnet & Magnet)	Magnetic Needle Counter (30 Count), Scalpel & Specialty Blade Remover, Syringe Needle Remover, Disposal Box	12/Bx 4 Bx/Cs	48

➤ **Contact your Ansell representative for ordering or more information.**

**North America**  
Tel: (866) 764-3327  
(818) 534 2500  
Email: sandel-cs@ansell.com

**Latin America & Caribbean**  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
Tel: +32 (0)2 528 74 00  
Email: info@ansell.eu

**Asia Pacific**  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

**Australia & New Zealand**  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com