



SANDEL® PREP-RESIST™ MARKERS

NEVER EVENT PROTECTION

Introducing PREP-RESIST™ surgical site markers with prep-resistant ink that doesn't overstay its welcome

SANDEL PREP-RESIST surgical site markers are resistant to chlorhexidine gluconate and isopropyl alcohol skin prep solutions. Specially formulated to show up on varying colors of skin, PREP-RESIST markers offer good visibility after prep but largely fade after a day so patients aren't burdened with lingering surgical site markings.

KEY FEATURES & BENEFITS

- Resistant to chlorhexidine gluconate and isopropyl alcohol skin prep solutions
- Gentian violet-free
- Available in standard length, petite, and 4-in-1™ dual-tip/dual-ink marker with chiseled tip
- Triangular “no-roll” barrel
- Bright orange and highly recognizable
- 4-in-1 marker comes with TIME OUT® sleeve



Post Prep: PREP-RESIST Marker vs. Traditional Skin Marker

SANDEL®

Anticipate the Unexpected

Developed in collaboration with clinical professionals, SANDEL products are designed to address key issues in high-risk areas in the healthcare environment. Our signature bright orange products help reduce the risk of preventable injuries, errors and infections. SANDEL products equip caregivers with solutions to do their job safely while improving patient outcomes.

SANDEL® PREP-RESIST™ MARKERS

NEVER EVENT PROTECTION

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
1015-NNS	PREP-RESIST Petite Skin Marker	Petite Prep-Resistant "No Roll" Skin Marker (Non-Sterile)	50/Bag 10 Bags/Cs	500
1029	PREP-RESIST Skin Marker	Prep-Resistant "No Roll" Skin Marker (Sterile)	25/Bag 4 Bx/Cs	100
1029-NNS	PREP-RESIST Skin Marker	Prep-Resistant, "No Roll" Skin Marker (Non-Sterile)	200/Bag 2 Bags/Cs	400
1061	PREP-RESIST 4-in-1 Marker	Dual Tip/Dual Ink Marker (Permanent and Prep-Resistant Skin Ink) with Chiseled Tip, TIME OUT Sleeve (Sterile)	25/Bx 4 Bx/Cs	100
1061-NNS	PREP-RESIST 4-in-1 Marker	Dual Tip/Dual Ink Marker (Permanent and Prep-Resistant Skin Ink) with Chiseled Tip, TIME OUT Sleeve (Non-Sterile)	100/Bag 4 Bags/Cs	400

➤ **Contact your Ansell representative for ordering or more information.**

North America
Tel: (866) 764-3327
(818) 534-2500
Email: sandel-cs@ansell.com