

SANDEL® STAT-BLOC™ ANTIMICROBIAL LINENS

O . R . T U R N O V E R S O L U T I O N S

Disposable, antimicrobial table covers, arm board covers and headrest covers to help prevent cross-contamination

STAT-BLOC™ table sheets provide three layers of protection:

1. Quilted, breathable top layer allows for fluid absorption and minimizes wrinkling to help prevent skin breakdown
2. Strong concentration of Super Absorbent Polymers (SAP) for maximum fluid control
3. Impervious backing with proprietary antimicrobial technology offers complete protection against bacterial migration¹

STAT-BLOC™ arm board covers feature dual, secured straps for added safety and to prevent straps from falling to the floor and becoming contaminated. STAT-BLOC™ headrest covers feature an integrated flap for full coverage.

KEY FEATURES & BENEFITS


- AMT™ Antimicrobial Technology is proven 99.9% effective against MRSA, E. coli and CRE and 99% effective against COVID-19¹
- Shown in testing to be more absorbent than leading table sheets and absorb 3x faster than the leading table sheet¹
- Table sheets are made of 49% renewable bio-based material which is carbon neutral when incinerated¹
- Also available in orange as a visible indicator of the antimicrobial backing
- Less lint than standard reusable linens to help prevent cross-contamination
- Draw sheet and arm board covers come folded inside the table sheet
- Table sheet and combo-fold packaging bags made of 20% recycling content

FEATURED TECHNOLOGIES



ansell.com/STAT-BLOC

¹ Data on file.



SANDEL® Anticipate the Unexpected

Through collaboration with clinical professionals, SANDEL® products are designed to address key high-risk issues affecting healthcare workers. Its signature high-visibility orange products and premium antimicrobial technology equip caregivers in the healthcare environment to do their job safely while improving outcomes.

SANDEL® STAT-BLOC™ ANTIMICROBIAL LINENS

O . R . T U R N O V E R S O L U T I O N S

Stocked and Ready To Order

ITEM #	PRODUCT NAME	DESCRIPTION	CASE QTY.
AS-92-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent, Quilted, Standard Table Sheet	20
AS-92-QZF	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent, Quilted, Fitted Table Sheet / Stretcher Sheet	20
AS-92-QXLF	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent, Quilted, Fitted Table Sheet with XL Absorbent Area	20
AS-100-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent, Quilted, XL Table Sheet	20
AS-120-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent, Quilted, XXL Table Sheet	20
AS-6030-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent, Quilted, Top and Bottom Section Table Sheet - Cystoscopy Table	20
CF-AS6030-LFS-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent Quilted Cysto Table Sheet, 2 pieces, with Disposable Standard White Draw Sheet	20
CF-AS9-LS-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent Quilted Table Sheet with Disposable Standard White Draw Sheet	20
CF-AS9-LS-QZP-C	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent, Quilted Table Sheet, Orange with Sunshine Pattern (Includes Standard Draw Sheet)	20
ABC-10	STAT-BLOC™ Arm Board Cover	Disposable, Antimicrobial Arm Board Cover with 2 Straps	200
ABC-XL	STAT-BLOC™ Arm Board Cover	Disposable, Antimicrobial, XL Arm Board Cover with 2 Straps	200
DPC-10	STAT-BLOC™ Pillow Cover	Disposable, Antimicrobial Pillow Cover	100
HR-10	STAT-BLOC™ Headrest Cover	Disposable, Antimicrobial Headrest Cover	250
LFS-4080	Draw Sheet	Disposable, XL Draw Sheet	50
LFS-66	Draw Sheet	Disposable, Standard Draw Sheet	50
LFS-66-BF	Draw Sheet	Disposable, Butterfly Draw Sheet	50
WGP-01	STAT-BLOC™ Gurney / Wagon Cover	Disposable, Antimicrobial Cover for Gurney or Wagon, Orange with Sunshine Pattern	50

✦ **Contact your Ansell representative for ordering or more information.**

North America
 US Tel: 866-764-3327
 818-534-2500
 Email: sandel-cs@ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking
 © 2024 Ansell Limited. All rights reserved.