

# SANDEL® STAT-BLOC™ ANTIMICROBIAL PEDIATRIC LINENS

SAFETY SOLUTIONS

## Reduce the risk of cross-contamination and healthcare-associated infections

### Disposable Antimicrobial Table Sheet

#### STAT-BLOC™ SHEETS

Printed with a sunshine pattern, specially designed for use with children

- 1 Quilted, breathable, hydrophilic layer**  
Allows for fluid absorption; eliminates wrinkling which can lead to skin breakdown
- 2 Super absorbent layer**  
Strong concentration of Super Absorbent Polymers (SAP)
- 3 Impervious, antimicrobial layer**  
Prevents bacteria from colonizing on the linen;
  - 99.9% effective against MRSA, CRE and E.coli\*
  - Kills 99% of COVID-19 which comes into contact with the antimicrobial backing\*



### Disposable Antimicrobial Wagon & Gurney Cover

#### STAT-BLOC™ WAGON & GURNEY COVER

Helps reduce the risk of cross-contamination between patients



- Impervious, antimicrobial layer
- High resistance to slipping
- Accommodates various wagon configurations
- Fits a standard gurney mattress
- Fun, child-friendly sunshine print
- Wagons provided with purchase of wagon covers



# SANDEL® STAT-BLOC™ ANTIMICROBIAL LINENS

SAFETY SOLUTIONS



Standard or custom packs are available with a variety of components to meet the needs of any department or specialty including the operating room, labor and delivery, cystoscopy, catheterization lab, radiology, endoscopy, and more.

- Improve consistency, efficiency and turnover time
- Enhance infection control
- Decrease lint to help reduce the risk of cross-contamination
- Help maintain patient skin integrity
- Eliminate costs related to reusable textiles
- Comply with guidelines from AORN
- Orange color is a visible indicator of the antimicrobial and provides separation from blue sterile drapes
- Extensive offering of turnover pack components to choose from including:
  - » STAT-STRAP™ Patient Positioning Straps
  - » STAT-BAG™ Heavy Duty Waste Bags
  - » STAT-SORB™ Solidifier
  - » Rayon, Microfiber and Strip Mops
  - » Non-Woven Wipes
  - » Cord Covers
  - » Patient Transfer Sheets
  - » Absorbent Floor Pads
  - » Disposable Patient Positioners



ITEM #	PRODUCT NAME	DESCRIPTION
CF-AS9-LS-QZP-C	STAT-BLOC™ Premium Table Sheet	Antimicrobial, Absorbent, Premium Table Sheet with Sunshine Pattern, Standard Lift Sheet
WGP-01	STAT-BLOC™ Wagon and Gurney Cover	Antimicrobial Wagon and Gurney Cover with Sunshine Pattern

**For more information, visit [www.ansell.com/sunny](http://www.ansell.com/sunny)**

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking)  
© 2022 Ansell Limited. All Rights Reserved.

\*Independently tested under controlled laboratory conditions. AATCC Test Method 100 and ISO ISO18184:201(E). Data on File. Lit. #130432 Rev. B.