

SANDEL® Z-SLIDER™ PATIENT TRANSFER SHEET

STAFF & PATIENT ERGONOMICS

A disposable patient transfer sheet designed to prevent back injuries to healthcare workers

The Z-Slider disposable patient transfer sheet helps prevent disabling back injuries that can result from lifting, transferring and repositioning patients. Z-Slider is placed under the draw sheet and acts to reduce friction while moving the patient.

Use of Z-Slider complies with recommendations by OSHA and the AORN for safe patient handling, as well as “No-Lift” laws in effect in many states.^{1,2}

To enhance staff compliance, Ansell offers several point-of-use dispensers for use in patient rooms or the O.R..

KEY FEATURES & BENEFITS

- Reduces the risk of physical strain to the back, shoulders, neck and arms caused by moving patients
- Single-patient use reduces the risk of cross-contamination
- Cost-effective
- Easy to use - product instructions are printed on the product
- Recyclable (LDPE 4)
- Bright orange and highly recognizable
- Dimensions: 33" x 39" (84 cm x 99 cm)



SANDEL Anticipate the Unexpected

Through collaboration with clinical professionals, SANDEL® products are designed to address key high-risk issues affecting healthcare workers. Its signature high-visibility orange products and premium anti-microbial technology equip caregivers in the healthcare environment to do their job safely while improving outcomes.

¹ Hospital eTool: HealthCare Wide Hazards Module: Ergonomics, U.S. Department of Labor, Occupational Safety and Health Administration (OSHA), www.osha.gov

² Pocket Reference Guide: Safe Patient Handling and Movement in the Perioperative Setting. AORN Toolkit. Association of periOperative Nurses. 2014

SANDEL® Z-SLIDER™ PATIENT TRANSFER SHEET

STAFF & PATIENT ERGONOMICS

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2101	Z-Slider - Single Pack	Patient Transfer Sheet. Individually Wrapped (Non-Sterile)	24/Bx 2 Bx/Cs	48
2101-H	Point-of-Use Clear Wall Holder	Point-of-Use Clear Wall Dispenser for Use with Z-Slider Single Pack	4/Cs	4
2102	Z-Slider - Roll	Patient Transfer Sheet (roll dispenser box) (Non-Sterile)	30/Bx 2 Bx/Cs	60
2103	Z-Slider Wall Dispenser (Large)	Patient Transfer Sheet packaged for Mounting on Any Surface (Non-Sterile)	15/Bx 4 Bx/Cs	60
2103-H	Point-of-Use Clear Wall Holder	Point-of-Use Clear Wall Mounting Holder for Use With 2103 (12" x 20.5")	5/Cs	5
2104	Z-Slider Wall Dispenser (Small)	Patient Transfer Sheet Refill for Point-of-Use Wall Mount Holder (Non-Sterile)	10/Bx 6 Bx/Cs	60
2104-H	Point-of-Use Clear Wall Holder	Point-of-Use Clear Wall Mounting Holder for the 2104 Dispenser Box (11.5" x 11")	5/Cs	5



2104
Z-Slider Refill for Compact
Point-of-Use Wall Mount Holder



2104-H
Z-Slider Compact Point-of-Use
Wall Mounting Holder



2102
Z-Slider Roll



2103
Z-Slider Refill for Full-Size
Point-of-Use Wall Mount Holder



2101
Z-Slider Single-Pack
Individually wrapped



2101-H
Z-Slider Single-Pack Point-
of-Use Wall Mounting Holder



2103-H
Z-Slider Full-Size Point-of-Use
Wall Mounting Holder

➤ **Contact your Ansell representative for ordering or more information.**

North America
US Tel: 866-764-3327
818-534-2500
Email: sandel-cs@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0)2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com