Ansell

ENCORE[®] Latex Moisturizing

SURGICAL GLOVES

Latex, powder-free surgical glove with a built-in moisturizer

KEY FEATURES AND BENEFITS:

- Contains glycerol and dimethicone as moisturizing agents to help retain moisture and rehydrate skin
- Ideal thickness to optimize tactile sensitivity and protection
- Proprietary glove formulation delivers a soft yet a robust surgical glove
- The DERMASHIELD[™] Technology internal coating together with the built-in moisturizer, facilitates dry and damp donning as well as intraoperative changes

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- Damp donning
- Extra comfort during long, sustained procedures
- Use as an underglove when double gloving

FEATURED TECHNOLOGIES

A.R.T Pinhole Barrier Technology

nology HYDRASOFT





ENCORE

and some

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

Ansell

ENCORE[®] Latex Moisturizing

Material	Natural Rubber Latex							
Powder Content	Powder-Free							
Color	Natural							
Thickness	Standard							
Cuff Style	Beaded with SUREFIT™ Technology							
External Glove Surface	Smooth							
Internal Glove Surface	Coated with HYDRASOFT™ Technology							
Allergy Prevention	None							
Grip Level	Moderate							
Double Gloving Recommendation	Underglove							
Tested For Use With Chemotherapy Drugs	No, not tested in accordance with ASTM D6978							
Freedom from Holes (Inspection level I)	0.65 AQL							
Sterile	Yes							
Product Standards	ASTM D3577, EN 420:2003 + A1:2009, ISO 10282, JIS T9107							
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001							
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case							
Shelf Life	2 years							
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.							
Vulcanization Chemical Accelerators	Zinc Dibut	yldithiocarb	amate (ZDB	C), Zinc Diet	hyldithiocar	bamate (ZDI	EC)	
ORDERING INFORMATIO	N							
Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	2018655	2018660	2018665	2018670	2018675	2018680	2018685	2018690

Contact your Ansell representative for ordering or more information.

North America US Tel: 800 952 9916 CA Tel: 844 494 7854 Email: insidesalesus@ansell.com Sandel US Tel: 866 764 3327 Sandel email: sandel-cs@ansell.com Central & South America Tel: +52(442) 296 20 50 Email: cslac@ansell.com Brazil Tel: +55-11-3356-3100 Email: <u>latam.medical@ansell.com</u>

Europe, Middle East & Africa Tel: +32 (0) 2 528 74 00 Email: info@ansell.eu Asia Pacific Tel: +603 8310 6688 Email: apac.medical@ansell.com

India Tel: +91 22 38124500 Email: info.india@ansell.com Australia & New Zealand Tel: +61 3 9270 7270 Email: protection@ap.ansell.com China Tel: +86(21) 3827 5005 Email: infochina@ansell.com



Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking. © 2022 Ansell Limited. All Rights Reserved.