

AlphaTec® 55-303 55-306 55-308



Chemical & Liquid Protection



Neoprene



IDEAL APPLICATIONS

- Biotechnologies
- Blending of liquids
- Chemical & Medical research
- Chemical handling
- Chemical industry
- Chemicals mixing/handling/transferring/pouring/blending
- Lab work and analysis, especially when resistance to acids is required
- Pharmaceuticals
- Product handling - chemical/biological
- Product protection applications

PRIMARY INDUSTRIES



STANDARDS AND CERTIFICATIONS

EN388:200 3	EN374:200 3	EN374:200 3	EN388:201 6	EN ISO 374- 1:2016Type A	EN ISO 374-5:2016
3101		ACJKL	3101A	JKLOPT	VIRUS

DESCRIPTION

- Effective protection for heavy duty chemical environments, delivery exceptional resistance to most acids, alcohols, oils, lubricants and hydrocarbons
- Excellent physical protection against abrasion and good protection against cuts and punctures
- The heavyweight thickness increases durability, augmenting chemical and physical resistance
- Resistant to Type I allergies

PRODUCT DETAILS

Coating material	Neoprene
Grip design	Smooth Finish
Cuff style	Gauntlet
Colour	Black
AQL	0.65
EN Sizes	9,10
Length (mm)	55-303: 813 (Port Size 200mm) 55-306: 813 (Port size 250mm) 55-308: 813 (Port size 300mm)
Packaging	10 pairs in black sealed bags per shipper carton

CATEGORY III

EXTRA FEATURES

Print date 05-12-2019

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J

Boulevard International 55, 1070 Brussels, Belgium

Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03

<http://www.ansell.eu> • E-mail info@ansell.eu