

MICRO-TOUCH™ Coated

M u l t i - P u r p o s e


Powder-free, latex medical exam gloves, coated for quick and frequent donning

- **Easy donning:** Polyacrylate inner coating ensures easy and fast donning, ideal for users who must switch gloves frequently
- **Elasticity:** Ansell's natural rubber latex formulation and optimized film thickness offer high levels of elasticity, flexibility, sensitivity and comfort
- **Medical exam grade:** CE Marked as a Class I Medical Device, it is additionally certified and CE Marked as Category III Personal Protective Equipment, ensuring more consistent quality and providing protection for a wide range of applications
- **Low-allergenicity latex:** The low-allergenicity latex reduces risk of contact dermatitis due to absence of thiurams

CAUTION: This product contains natural rubber latex, which may cause allergic reactions

Recommended For

- Brief, low-risk examination procedures
- Damp donning



MICRO-TOUCH™ Everyday, Everywhere Protection

MICRO-TOUCH™ premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH™ Coated

M u l t i - P u r p o s e

PRODUCT INFORMATION					
Material	Natural Rubber Latex		Not Made From Natural Rubber Latex	No	
Color	White		Cuff Style	Beaded	
Cuff Length	Standard		External Glove Surface	Textured Fingers	
Freedom from Holes (Inspection level I)	1.5 AQL		Finger Thickness (mm/mil)	0.14 / 5.5	
Palm Thickness (mm/mil)	0.11 / 4.3		Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	
Tested For Use With Chemotherapy Drugs	No		Silicone Free	Yes	
Sterile	No		Antistatic	Not Tested	
Glove Length (mm/inches)	240 / 9.45		Product Segmentation	Multi-Purpose	
ORDERING INFORMATION					
Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	553301	553302	553303	553304	553305
PACKAGING AND STORAGE					
Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *80 gloves per dispenser (Size XL) *800 gloves per case (Size XL)				
Shelf Life	3 Years				
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.				
STANDARDS AND CERTIFICATIONS					
CE 0493, Category III, Medical Device Regulation (EU) 2017/745 Class I, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 11193-1:2008, ISO 13485, ISO 9001, ISO 10993-10:2010, UKCA					

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Region
BR +55 (11) 3356-3100
latam.medical@ansell.com
MX: +52 442 296 2050
latam.medical@ansell.com
CO: +57 1 288 3247
latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com