

MICRO-TOUCH™ Nitrile (White)

Standard Thin

White nitrile medical exam glove with textured fingertips

- **Heightened tactility:** Due to its lightweight design, this glove provides enhanced tactility and dexterity
- **Reduced allergy risk:** Nitrile material prevents Type I (latex) allergies
- **Smart packaging:** 150-count packaging reduces waste
- **Secure grip:** Textured fingertips ensures secure handling of medical instruments
- **Easy donning and doffing:** The polymer inner coating reduces friction for easy donning and doffing
- **Specialized protection:** Tested and approved for use with chemotherapy drugs*



*see product packaging for details

Recommended For

- Tested for use with chemotherapy drugs
- Standard, moderate risk examination procedures
- Protection from Type I latex allergy in HCW's or patients
- Environmentally friendly facilities

FEATURED TECHNOLOGIES



MICRO-TOUCH™ Everyday, Everywhere Protection

MICRO-TOUCH™ premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH™ Nitrile (White)

Standard Thin

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	White	Cuff Style	Beaded
Cuff Length	Standard	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.11/4.3
Palm Thickness (mm/mil)	0.07/2.8	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Silicone Free	No	Sterile	No
Antistatic	Not Tested	Glove Length (mm/inches)	240 / 9.45
Product Segmentation	Standard Thin		

PACKAGING AND STORAGE

Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, CE 0493, Category III, Medical Device Regulation (EU) 2017/745 Class I, EN 455 1-4, EN 421:2010, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600-US Food Contact Approved, ISO 9001, ISO 10993-10:2002, ISO 13485, Food Type: Rubber, UKCA

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean Region
BR +55 (11) 3356-3100
latam.medical@ansell.com
MX: +52 442 296 2050
latam.medical@ansell.com
CO: +57 1 288 3247
latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com