

## Sterile, three-layered preparation mat, with leak-proof properties and superb absorption

- **Advanced defenses:** The BioClean™ Chem Prep Mat S-BCPM is composed of three layers – a mesh upper layer allows even liquid distribution, while the middle layer provides excellent absorption, and the impermeable bottom layer protects against leakages
- **Reduced allergy and contamination risks:** This latex-free, low-linting BioClean™ Chem Prep Mat S-BCPM is sterilized by gamma irradiation and individually wrapped to minimize contamination risks
- **Assured quality:** This preparation mat is manufactured in a facility meeting ISO 9001 standards, for product quality you can depend on



### KEY FEATURES & BENEFITS

- ◊ **Three-layered construction:** For highly leak-proof absorbency
- ◊ **Gamma irradiation sterilization:** Suitable for use in an ISO 5 cleanroom
- ◊ **Manufactured in an ISO 9001 facility:** For product quality assurance

### Industries

- Controlled and Critical Environments
- Laboratory and Research
- Production and Manufacturing



## S-BCPM

Sterile Prep-Mat

### TECHNICAL DATA SHEET

|                                   | Product Information                                                                                                                                                                                                                                                                                                                              |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material                          | Top layer - Polyolefin apertured film + Middle layer - Air-laid paper + Bottom layer - Polyethylene                                                                                                                                                                                                                                              |
| Color                             | Blue                                                                                                                                                                                                                                                                                                                                             |
| Manufacturing/QMS Audit Standards | Manufacturing QMS Audit Standards ISO 9001                                                                                                                                                                                                                                                                                                       |
| Regulatory/Standard Compliance    | ISO 9001:2015                                                                                                                                                                                                                                                                                                                                    |
| Packaging                         | S-BCPM4025B: One piece per sealed inner PE bag; one inner bag per second PE bag; 50 double bags per sealed outer PE bag; four outer bags per lined carton (200 pieces)<br>S-BCPM5640B: One piece per sealed inner PE bag; one inner bag per second PE bag; 25 double bags per sealed outer PE bag; four outer bags per lined carton (100 pieces) |
| Storage                           | Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.                                                                                                                                                                                                                                                              |
| Available sizes                   | 40cm x 25cm, 56cm x 40cm                                                                                                                                                                                                                                                                                                                         |
| Sterilization Method              | GAMMA irradiation (25 kGy)                                                                                                                                                                                                                                                                                                                       |
| Sterilization Minimum Dose        | 25kGy                                                                                                                                                                                                                                                                                                                                            |
| Sterility Assurance Level         | 10 <sup>-6</sup>                                                                                                                                                                                                                                                                                                                                 |
| Cleanroom Class                   | Class 100/ISO Class 5 & EU GMP Grade A/B and other sterile cleanrooms                                                                                                                                                                                                                                                                            |
| Shelf Life                        | 3 years                                                                                                                                                                                                                                                                                                                                          |



## S-BCPM

Sterile Prep-Mat

### ORDERING INFORMATION

| SIZE        | 40cm x 25cm  | 56cm x 40cm  |
|-------------|--------------|--------------|
| REORDER NO. | S-BCPM-4025B | S-BCPM-5640B |

### PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

#### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

#### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

#### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

#### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

#### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

#### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

### Technology



Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:  
[www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

