

GAMMEX™ PI Hybrid

SURGICAL GLOVES

Hybrid synthetic surgical gloves offering the comfort of polyisoprene and strength of neoprene

- **Enhanced comfort and durability:** Made with Ansell's proprietary synthetic blend technology, these non-latex surgical gloves combine the latex-like comfort of polyisoprene (PI) with the durability of neoprene so wearers experience the best of both worlds
- **Reduced allergy risks:** Being a non-latex surgical glove, GAMMEX™ PI Hybrid eliminates the risk of a latex allergy while a range of allergens* known to cause allergic contact dermatitis are eliminated
- **Certified ergonomic performance:** This PI surgical glove is certified by United States Ergonomics with a reported 32.3% reduction of hand muscle effort**
- **Double gloving advantage:** This semi-transparent glove is an excellent outer glove for double gloving with a colored under glove to offer easy breach detection
- **Chemo-tested:** ASTM D6978 tested and may be used with a range of chemo drugs***

*Free from natural rubber latex proteins and do not contain Zinc 2-mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG) and cetylpyridinium chloride (CPC)


** When compared to bare hand performance

***Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- Protection from Type I latex allergy in HCW's or patients
- Use as an outer glove when double gloving
- Easy identification of a glove breach

FEATURED TECHNOLOGIES



GAMMEX™ Innovation that Protects

Driven by innovation, GAMMEX™ surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX™ PI Hybrid

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Polyisoprene, Neoprene
Color	Natural
Thickness	Standard
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.4 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, ASTM F739, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC), Sodium di-n-butyldithiocarbamate (SDBC), Xanthogen Polysulphides

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340063055	340063060	340063065	340063070	340063075	340063080	340063085	340063090

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico
 Tel: +52 442 296 2050
 Email: latam.medical@ansell.com
Ansell Brasil Ltda
 +55 (11) 3356-3100
 Email: latam.medical@ansell.com
Ansell Colombia SAS
 Tel: +57 1 288 3247
 Email: latam.medical@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu
Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com