

GAMMEX™ PI Plus Glove-in-Glove™ System

S U R G I C A L G L O V E S

Polyisoprene inner and outer surgical gloves, pre-donned and ready for double gloving

- **Effortless double gloving:** The GAMMEX® PI Plus Glove-in-Glove™ System is a unique two-in-one system with both inner and outer polyisoprene surgical gloves already pre-donned and aligned at the fingertips, enabling wearers to double glove in one don and save half the donning time
- **Sustainable packaging:** Being pre-donned, the inner and outer sterile gloves are packed in one inner wrap for a 50% reduction of paper
- **Advanced allergy prevention:** These non-latex gloves present no risk of latex allergies and do not contain chemical accelerators known to cause allergic contact dermatitis
- **Ergonomic comfort:** Certified by United States Ergonomics for superb wearer comfort and performance
- **Effortless donning and doffing:** Ansell's donning aid technology facilitates easy donning and doffing of the outer glove during intraoperative changes



Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Easy identification of a glove breach

FEATURED TECHNOLOGIES



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



GAMMEX™ PI Plus Glove-in-Glove™ System

SURGICAL GLOVES

PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	White, Green
Thickness	Standard
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	High
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D7160, ASTM F739, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	2 pairs per pouch; 50 pairs per box (25 pouches); 4 boxes per carton/case; 200 pairs per carton/case (100 pouches)
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthates, Zinc Dialkyldithiocarbamate

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340082055	340082060	340082065	340082070	340082075	340082080	340082085	340082090

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico
Tel: +52 442 296 2050
Email: latam.medical@ansell.com
Ansell Brasil Ltda
+55 (11) 3356-3100
Email: latam.medical@ansell.com
Ansell Colombia SAS
Tel: +57 1 288 3247
Email: latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2025 Ansell Limited. All Rights Reserved.