

# GAMMEX™ PI Radiation Attenuation

SPECIALTY GLOVES

## Polyisoprene sterile gloves, protecting hands from risks of scattered radiation exposure


- **Enhanced scattered radiation protection:** Suitable for use with C-Arms, these lead-free gloves reduce skin-dose exposure to scattered radiation by 63.0%, compared to just 1%\*\* with standard surgical gloves
- **Assured dexterity and tactile sensitivity:** Ansell's proprietary, soft, malleable polyisoprene formulation ensures the dexterity, tactile sensitivity and comfort needed for image-guided interventional procedures (fluoroscopy)
- **Reduced allergy risks:** Being a non-latex surgical glove, GAMMEX® PI Radiation Attenuation eliminates the risk of a latex allergy
- **Strengthened grip:** These black specialty gloves also feature a textured finish for good grip control

\*\*Wagner LK, Mulhern OR. Radiation-attenuating surgical gloves: effects of scatter and secondary electron production. Radiology.1996;200(1):45-48

**CAUTION:** These radiation attenuation gloves are designed to protect hands from scattered radiation exposure originating from the X-ray beam during fluoroscopic procedures. They are not intended for use in direct or primary X-ray beams.

### Recommended For

- Protecting hands from image guided interventional procedures (fluoroscopy)
- Protection from Type I latex allergy in HCW's or patients



### GAMMEX™ Innovation that Protects

Driven by innovation, GAMMEX™ surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



# GAMMEX™ PI Radiation Attenuation

SPECIALTY GLOVES

## PRODUCT INFORMATION

|                                         |                                                                                                          |
|-----------------------------------------|----------------------------------------------------------------------------------------------------------|
| Material                                | Polyisoprene with Tungsten                                                                               |
| Color                                   | Black                                                                                                    |
| Thickness                               | Standard                                                                                                 |
| Cuff Style                              | Straight                                                                                                 |
| External Glove Surface                  | Textured Fingertips                                                                                      |
| Internal Glove Surface                  | Chlorinated                                                                                              |
| Allergy Prevention                      | Latex (Type I)                                                                                           |
| Grip Level                              | Moderate                                                                                                 |
| Double Gloving Recommendation           | Underglove                                                                                               |
| Tested For Use With Chemotherapy Drugs  | No, not tested in accordance with ASTM D6978                                                             |
| Freedom from Holes (Inspection level I) | 0.65 AQL                                                                                                 |
| Sterile                                 | Yes                                                                                                      |
| Product Standards                       | ASTM F2547, EN 455 Part 2, EN 61331-1                                                                    |
| Quality / Environmental Standards       | ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001                                                         |
| Packaging                               | 5 pairs/box; 1 box/case                                                                                  |
| Shelf Life                              | 3 years                                                                                                  |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |

## ORDERING INFORMATION

| Size         | 6.5      | 7        | 7.5      | 8        | 8.5      | 9        |
|--------------|----------|----------|----------|----------|----------|----------|
| Product Code | 20873565 | 20873570 | 20873575 | 20873580 | 20873585 | 20873590 |

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 844 494 7854  
Email: insidesalesus@ansell.com  
Sandel US Tel: 866 764 3327  
Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico  
Tel: +52 442 296 2050  
Email: latam.medical@ansell.com  
Ansell Brasil Ltda  
+55 (11) 3356-3100  
Email: latam.medical@ansell.com  
Ansell Colombia SAS  
Tel: +57 1 288 3247  
Email: latam.medical@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu  
Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

India  
Tel: +91 22 38124500  
Email: info.india@ansell.com  
Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com  
China  
Tel: +86(21) 3827 5005  
Email: infochina@ansell.com

**ansell.com**

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking).  
© 2025 Ansell Limited. All Rights Reserved.