

Offers reliable multi-purpose protection for everyday activities with extra durability

- **Reliable protection:** Meets ASTM F1671 requirements for protection against viruses, bacteria, and other micro-organisms
- **Medical exam grade:** Ideal for safely cleaning, disinfecting, and treating wounds — helps protect against unseen hazards during cleaning and sanitization of commonly used spaces
- **Secure grip:** Textured fingertips provide a great grip
- **Genuine versatility:** Offers reliable, multi-purpose protection for everyday tasks, including virus/bacteria protection, first aid, cleaning and food preparation



Industries

- Food Processing
- Janitorial / Sanitation
- Healthcare
- Automotive
- Life Sciences

Recommended For

- Food processing and handling
- Standard patient examinations
- Laboratory analysis
- Catering/Meal preparation
- Light assembly tasks
- Light duty maintenance
- Multipurpose use: Dish Washing, Cleaning, Pet Care & Grooming, Food Preparation, Gardening & many more!

TECHNICAL DATA SHEET

| | Product Information |
|------------------------|--|
| Material | Nitrile |
| External Glove Surface | Textured Fingers |
| Audit Standards | ISO 13485, ISO 9001 |
| Product Certification | Medical Device Regulation (EU) 2017/745 Class I, Personal Protective Equipment Regulation (EU) 2016/425 Category III risks, ASTM D6319, FDA21 CFR 177-2600-US Food Contact Approved, ISO 13485, ISO 9001 |
| Packaging Overview | 100 gloves per dispenser/box 10 dispensers/box per carton/case |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |
| Country Of Origin | China, Malaysia |
| Product Segmentation | Standard Thin |
| Antistatic | No |
| Silicone Free | Yes |


PHYSICAL PROPERTIES

| | Typical Values | | Testing Method |
|---|----------------|-------------|------------------|
| Length (mm/inches) | 240 / 9.5 | | ASTM D6319 |
| Freedom from Holes (Inspection level I) | 1.5 AQL | | ASTM D5151 |
| Palm Thickness (mm/mils) | 0.07 / 2.8 | | ASTM D6319 |
| Finger Thickness (mm/mils) | 0.11 / 4.3 | | ASTM D6319 |
| | BEFORE AGING | AFTER AGING | |
| Ultimate Tensile Strength (MPa) | 14 | 14 | ASTM D412 & D573 |
| Elongation at Break (%) | 500 | 400 | ASTM D412 & D573 |
| Force at break (N) | 6 | 6 | EN 455-2 |

ORDERING INFORMATION

| Size | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|-------------|-------------|-------------|---------------|
| Product Code | 10733070 | 10733080 | 10733090 | 10733100 |

Performance Standards and Regulatory Compliance



For additional information visit us at www.ansell.com, or call us at

| | | |
|--|---|--|
| Europe, Middle East & Africa Region Ansell Healthcare Europe NV T: +32 (0) 2 528 74 00 F: +32 (0) 2 528 74 01 Asia Pacific Region MY +603 8310 6688 apac.medical@ansell.com CN: +86 21 38275000 infochina@ansell.com JP: +813 5549 8151 info.medical.jp@ansell.com | North America Region Ansell Healthcare Products LLC T: +1 800 800 0444 F: +1 800 800 0445 Latin America & Caribbean Region BR +55 (11) 3356-3100 latam.medical@ansell.com MX: +52 442 296 2050 latam.medical@ansell.com CO: +57 1 288 3247 latam.medical@ansell.com | Australia Ansell Limited T: +61 1800 337 041 F: +61 1800 803 578 Russia Ansell PYC Ten. +7 495 258 13 16 |
|--|---|--|

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

