

MICRO-TOUCH™ Coated

M u l t i - P u r p o s e


Powder-free, latex medical exam gloves, coated for quick and frequent donning

- **Easy donning:** Polyacrylate inner coating ensures easy and fast donning, ideal for users who must switch gloves frequently
- **Elasticity:** Ansell's natural rubber latex formulation and optimized film thickness offer high levels of elasticity, flexibility, sensitivity and comfort
- **Medical exam grade:** CE Marked as a Class I Medical Device, it is additionally certified and CE Marked as Category III Personal Protective Equipment, ensuring more consistent quality and providing protection for a wide range of applications
- **Low-allergenicity latex:** The low-allergenicity latex reduces risk of contact dermatitis due to absence of thiurams

CAUTION: This product contains natural rubber latex, which may cause allergic reactions

Recommended For

- Brief, low-risk examination procedures
- Damp donning



MICRO-TOUCH™ Everyday, Everywhere Protection

MICRO-TOUCH™ premium medical examination gloves are scientifically designed to complement every need of medical professionals.

Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH™ Coated

M u l t i - P u r p o s e

PRODUCT INFORMATION

Material	Natural Rubber Latex	Not Made From Natural Rubber Latex	No
Color	White	Cuff Style	Beaded
Cuff Length	Standard	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.14 / 5.5
Palm Thickness (mm/mil)	0.11 / 4.3	Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Tested For Use With Chemotherapy Drugs	No	Silicone Free	Yes
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.45	Product Segmentation	Multi-Purpose

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6-6.5)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	303001XS	303001SM	303001MD	303001LG	303001XL

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser; 10 dispensers/1000 gloves per case; *XL 80 gloves per dispenser; 800 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

CE 0493, Category III, Medical Device Regulation (EU) 2017/745 Class I, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 11193-1:2008, ISO 13485, ISO 9001, ISO 10993-10:2010, UKCA

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean Region
BR +55 (11) 3356-3100
latam.medical@ansell.com
MX: +52 442 296 2050
latam.medical@ansell.com
CO: +57 1 288 3247
latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com