

Previously known as: MICROGARD® 1500 PLUS Lab Coat - Model 204

### Breathable, silicone-free lab coat made from anti-static, SMS fabric.

- **Optimized breathability:** AlphaTec® 1500 PLUS Lab Coat Bound – Model 204 delivers comfortable protection in general industrial environments, thanks to its air and water vapor permeable as well as anti-static, SMS material
- **Reduced contamination risks:** As this protective outer garment is silicone-free, it minimizes contamination risks, making it ideal for the life sciences and metal fabrication industries
- **Improved ergonomics:** In addition, this white protective over garment offers an optimized fit, which improves wearer comfort and workplace safety



## Industries

- Life Sciences
- Metal fabrication


## Applications

- Unpacking of bulk goods
- Loading and unloading of process equipment
- Production line support & maintenance
- Assembly & inspection of components


## KEY FEATURES & BENEFITS

- **Air and water vapor permeable:** For comfortable personal protection
- **Silicone-free design:** Barrier lab coat ensuring clean surfaces
- **Optimized body fit:** Protective garment combining comfort and safety


## Material Diagram




## Seam Diagram





Performance Standards

  
0598

EN ISO 13688:2013

  
PARTIAL BODY  
PROTECTION ONLY

  
TYPE PB [6]

  
TYPE PB [6]

Product Information

Product Information	
Available Sizes	M, L, XL, 2XL
Color	White
Country Of Origin	China
Packaging Overview	Individually packaged
Product Material	SMS fabric with anti-static treatment
Product Reference	1500 PLUS - Lab Coat Bound Model 204
Standards Overview	CE CATEGORY III (Type PB [6])
Seam Type	Bound
Shelf Life	5 years
Style Number	W15B204

Suit components & accessories

Suit Components and Accessories	
Design Features	Medium length labcoat with collar, studded front fastening, left breast and right thigh pockets.

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

**Europe, Middle East & Africa Region**  
Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00

**North America Region**  
Ansell Healthcare Products LLC  
T: +1 800 800 0444

**Australia**  
Ansell Limited  
T: +61 1800 337 041

**Asia Pacific Region**  
Ansell Global Trading Center  
T: +603 8310 6688

**Latin America & Caribbean Region**  
Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

