



Extra BLAS

Sterile Disposable Latex Cleanroom Glove

Sterile elbow-length cleanroom gloves, offering comfortable, flexible personal protection

- **Heightened comfort and flexibility:** BioClean™ Extra BLAS elbow-length gloves' soft, stretchy film composition enhances flexibility and comfort
- **Added protection:** Their elbow-length design ensures forearm, wrist and hand protection, limiting cross contamination risks
- **Reduced contamination risks:** BioClean™ Extra BLAS sterile gloves are also aseptically packed, gamma irradiation sterilized (minimum SAL of 10⁻⁶) and manufactured in accordance with ISO standards, further minimizing contamination risks, and they feature non-particulating EasyTear™ packaging
- **Enhanced strength:** These cleanroom gloves' textured surface and beaded cuff enhance glove strength and stability, reducing the risk of tears during donning
- **Increased practicality:** In addition, BioClean™ Extra BLAS controlled environment gloves are designed for easy double-donning, for everyday convenience



CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

KEY FEATURES & BENEFITS

- ◊ **Soft, stretchy film:** For PPE gloves that offer flexible comfort
- ◊ **Elbow-length design:** Extra coverage, fewer cross contamination risks
- ◊ **Sterility assurance level*:** Aseptically packed, gamma irradiation

www.ansell.com

*minimum SAL of 10⁻⁶

Ansell



Extra BLAS

Sterile Disposable Latex Cleanroom Glove

TECHNICAL DATA SHEET

	Product Information
Material	Natural Rubber Latex
Color	Natural
Shape	Ambidextrous
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory/Standard Compliance	CE 0493, EN 455 Part 2, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 421:2010, Food Contact, Category III, UKCA
Packaging	1 pair per inner PE wallet; one wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs) More sustainable packaging: These sterile gloves are packed in recyclable plastic packaging and delivered in recycled cardboard shipper cases. Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Country of Origin	India
Available sizes	S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Chlorinated
Sterilization Method	GAMMA irradiation (25 kGy)
Sterilization Minimum Dose	25kGy
Sterility Assurance Level	10 ⁻⁶
Cleanroom Class	Class 10/ISO Class 4, EU GMP Grade A/B and other sterile cleanrooms
Shelf Life	3 Years
Tested for use with Chemotherapy Drugs	No
Protein Level	50 µg/g or less of total extractable protein
Anti-static	No



Extra BLAS

Sterile Disposable Latex Cleanroom Glove

Physical Properties					Testing Method
Sizes	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	EN ISO 21420
Typical Length (mm/in)	410 / 16				
Palm Width (mm/in)	86/3.38	96/3.78	106/4.17	115/4.53	
Freedom from Holes	1.5 AQL Performance Level 2				EN 374-2
Typical Particle Count ≥0.5µm (counts / cm²)	<1500				IEST-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.18 / 7.09				EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.21 / 8.27				EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.13 / 5.12				EN 455-2
Force at Break (N) During Aging	≥ 9 N				EN 455-2

IONIC CONTENT

Concentration in µg/cm ²	Typical	Concentration in µg/cm ²	Typical
Ammonium	0.333	Nitrate	0.528
Bromide	0.053	Nitrite	Not Detected
Calcium	1.155	Phosphate	0.014
Chloride	4.412	Potassium	0.075
Fluoride	Not Detected	Sodium	0.025
Lithium	0.230	Sulphate	0.267
Magnesium	Not Detected	Zinc	Not Detected

ORDERING INFORMATION

SIZE	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
REORDER NO.	BLAS-S	BLAS-M	BLAS-L	BLAS-XL

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

