

ENCORE™ Non-Latex PI Hybrid

S U R G I C A L G L O V E S

A hybrid surgical glove ideally balancing polyisoprene (PI) comfort and durability of neoprene

- **Optimal comfort and durability:** These surgical gloves' feature Ansell's patented hybrid blend technology enabling wearers to enjoy the best of two synthetic glove materials: the latex-like comfort of PI and the superior durability of neoprene
- **Heightened allergy prevention:** Being a synthetic glove, ENCORE® Non-Latex PI Hybrid eliminates the risk of a latex allergy while a range of allergens* known to cause allergic contact dermatitis are eliminated
- **Enhanced safety:** When used as an outer glove for double gloving, this semi-transparent glove allows for fast and easy breach detection
- **Strengthened grip:** These hybrid surgical gloves feature a smooth surface with a micro-textured finish for secure instrument handling
- **Chemo-tested:** These sterile gloves comply with ASTM D6978 standards, making them suitable for use with various chemotherapy drugs**

* Free from natural rubber latex proteins and do not contain Zinc 2-mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG) and cetylpyridinium chloride (CPC)

** Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- Easy identification of a glove breach
- Protection from Type I latex allergy in HCW's or patients
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES

 **A.R.T.**
Pinhole Barrier Technology

 **DERMASHIELD**
Donning Aid Technology



ENCORE® The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

ENCORE™ Non-Latex PI Hybrid

SURGICAL GLOVES

PRODUCT INFORMATION

Material	PI (Polyisoprene), Neoprene (Polychloroprene)
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM D6978, ASTM F739, CE 0493 (PPE) + CE 2797 (MD) , EN 16523-1, EN 374-2 and -4, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC), Sodium di-n-butylthiocarbamate (SDBC), Xanthogen Polysulphides

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340105055	340105060	340105065	340105070	340105075	340105080	340105085	340105090

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico
 Tel: +52 442 296 2050
 Email: latam.medical@ansell.com
Ansell Brasil Ltda
 +55 (11) 3356-3100
 Email: latam.medical@ansell.com
Ansell Colombia SAS
 Tel: +57 1 288 3247
 Email: latam.medical@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu
Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com