

# GAMMEX™ Non-Latex PI Micro

S U R G I C A L   G L O V E S

## Polyisoprene, skin-friendly surgical gloves, for exceptional dexterity and tactile sensitivity

- **Elevated tactile sensitivity and dexterity:** The combination of Ansell's unique soft, yet strong, polyisoprene formulation and process helps make GAMMEX® Non-Latex PI Micro surgical gloves 15% thinner\*, delivering exceptional surgical finesse and hand protection
- **Skin-friendly:** These non-latex surgical gloves eliminate any risk of developing a Type I latex allergy, while their ground-breaking non-sensitizing polyisoprene (PI) formulation eliminates the use of chemical accelerators known to cause allergic contact dermatitis
- **Excellent grip efficiency:** The unique micro-textured external finish offers safe grip when handling surgical instruments under both dry and damp conditions
- **Chemo-tested:** Having been tested per ASTM D6978, these sterile surgical gloves may be safely used with a range of chemotherapy drugs\*\*\*



\*\*\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

### Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

### FEATURED TECHNOLOGIES



### GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



# GAMMEX™ Non-Latex PI Micro

## SURGICAL GLOVES

### PRODUCT INFORMATION

|                                         |                                                                                                                                                                |
|-----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material                                | PI (Polyisoprene)                                                                                                                                              |
| Powder Content                          | Powder-Free                                                                                                                                                    |
| Color                                   | White                                                                                                                                                          |
| Thickness                               | Micro                                                                                                                                                          |
| Cuff Style                              | Beaded with SUREFIT™ Technology                                                                                                                                |
| External Glove Surface                  | Smooth with Micro-Textured Finish                                                                                                                              |
| Internal Glove Surface                  | Polymer Coated with DERMASHIELD™ Technology                                                                                                                    |
| Allergy Prevention                      | Latex (Type I)                                                                                                                                                 |
| Grip Level                              | Moderate                                                                                                                                                       |
| Double Gloving Recommendation           | Outer glove                                                                                                                                                    |
| Tested For Use With Chemotherapy Drugs  | Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)                                                                                             |
| Freedom from Holes (Inspection level I) | 0.65 AQL                                                                                                                                                       |
| Sterile                                 | Yes                                                                                                                                                            |
| Product Standards                       | AS/NZS 4179, ASTM D3577, ASTM D7160, CE 0493 (PPE) + CE 2797 (MD) , EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, JIS T9107 |
| Quality / Environmental Standards       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP                                                                                           |
| Packaging                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case                                                                                           |
| Shelf Life                              | 3 years                                                                                                                                                        |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                                                       |
| Vulcanization Chemical Accelerators     | Xanthates, Zinc Dialkyldithiocarbamate                                                                                                                         |

### ORDERING INFORMATION

| Size         | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 340057055 | 340057060 | 340057065 | 340057070 | 340057075 | 340057080 | 340057085 | 340057090 |

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 844 494 7854  
Email: insidesalesus@ansell.com  
Sandel US Tel: 866 764 3327  
Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico  
Tel: +52 442 296 2050  
Email: latam.medical@ansell.com  
Ansell Brasil Ltda  
+55 (11) 3356-3100  
Email: latam.medical@ansell.com  
Ansell Colombia SAS  
Tel: +57 1 288 3247  
Email: latam.medical@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu  
Asia Pacific  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com

India  
Tel: +91 22 38124500  
Email: info.india@ansell.com  
Australia & New Zealand  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com  
China  
Tel: +86(21) 3827 5005  
Email: infochina@ansell.com

**ansell.com**

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking).  
© 2025 Ansell Limited. All Rights Reserved.