



P-Zero BPZS

Sterile Disposable Neoprene
(Polychloroprene) Cleanroom Glove

Sterile polychloroprene cleanroom gloves, offering chemical and chemotherapy drug resistance

- **Specialized protection:** BioClean™ P-Zero BPZS are tested against ASTM D6978 and EN standards for chemotherapy drug handling and chemical resistance
- **Reduced contamination risks:** These latex- and powder-free free cleanroom gloves come in EasyTear packaging, limiting potential contamination risks
- **Heightened flexibility and comfort:** Their stretchy polychloroprene formulation ensures excellent flexibility, dexterity and comfort
- **More sustainable packaging:** These polychloroprene gloves' plastic packaging is recyclable* and delivered in cardboard shipper cases made from recycled content

*Always check your local recyclable status as these materials may not be considered suitable for recycling in your location



KEY FEATURES & BENEFITS

- **ASTM and EN standards:** Protection against cytotoxic drugs
- **ESD properties:** For safeguarded component integrity
- **Non-particulating EasyTear packaging:** Reduced contamination risks

Industries

- Biotechnology Manufacturing
- Controlled and Critical Environments
- Laboratory and Research
- Medical Device Manufacturing
- Pharmaceutical Manufacturing
- Production and Manufacturing

WWW.ANSELL.COM

Ansell



P-Zero BPZS

Sterile Disposable Neoprene
(Polychloroprene) Cleanroom Glove

TECHNICAL DATA SHEET

	Product Information
Material	Neoprene (Polychloroprene)
Color	Natural
Shape	Hand Specific
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory	ASTM D6978, CE 0493, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 455 Part 2, Category III, UKCA
Packaging	One pair per inner PE wallet; one wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs) More sustainable packaging: These sterile neoprene gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases. *Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Country of Origin	Indonesia
Available sizes	6, 6.5, 7, 7.5, 8, 8.5, 9
Powder Content	Powder-Free
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated
Sterilization Method	GAMMA irradiation (25 kGy)
Sterilization Minimum Dose	25kGy
Sterility Assurance Level	10 ⁻⁶
Cleanroom Class	Class 10/ISO Class 4 & EU GMP Grade A/B and other sterile cleanrooms
Shelf Life	Five (5) years from date of manufacture
Tested for use with Chemotherapy Drugs	Yes
Protein Level	N/A: contains no natural rubber latex
Anti-static	No
Vulcanization Chemical Accelerators	None



P-Zero BPZS

Sterile Disposable Neoprene
(Polychloroprene) Cleanroom Glove

Physical Properties								Testing Method
Sizes	6	6.5	7	7.5	8	8.5	9	EN ISO 21420
Typical Length (mm/in)	300 / 12							
Palm Width (mm/in)	77/3	83/3.2	89/3.5	95/3.7	102/4	108/4.3	114/4.4	
Freedom from Holes	0.65 AQLPerformance Level 3							EN 374-2
Typical Particle Count ≥0.5µm (counts / cm²)	<1300							IENT-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.15 / 5.91							EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.18 / 7.09							EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.11 / 4.33							EN 455-2
Force at Break (N) During Aging	≥ 9 N							EN 455-2

IONIC CONTENT

Concentration in µg/cm ²	Typical	Concentration in µg/cm ²	Typical
Ammonium	0.005	Nitrate	1.63
Bromide	Available on File	Nitrite	Available on File
Calcium	0.929	Phosphate	Available on File
Chloride	0.25	Potassium	0.42
Fluoride	Available on File	Sodium	0.15
Lithium	Available on File	Sulphate	0.070
Magnesium	Available on File	Zinc	Available on File

ORDERING INFORMATION

SIZE	6	6.5	7	7.5	8	8.5	9
REORDER NO.	BPZS60	BPZS65	BPZS70	BPZS75	BPZS80	BPZS85	BPZS90

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

