

# GAMMEX® Non-Latex PI Underglove

## SURGICAL GLOVES

### Soft, skin-friendly polyisoprene surgical gloves, for heightened tactile sensitivity

- **Superior comfort and protection:** Ansell's proprietary soft yet strong formulation delivers heightened comfort and hand protection when double gloving
- **Easy breach detection:** Green in color, these sterile gloves enable fast and easy breach detection when worn with a cream or white outer glove
- **Optimal allergy protection:** Being non-latex, these polyisoprene surgical gloves eliminate the risk of developing a latex allergy while their skin-friendly formulation helps minimize sensitization to developing a Type IV chemical allergy as chemical accelerators known to cause skin allergies have been eliminated
- **Effortless donning and doffing:** A unique internal coating enables easy donning and doffing under dry or damp hand conditions



### Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Easy identification of a glove breach
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

### FEATURED TECHNOLOGIES



### GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX® Non-Latex PI Underglove

## SURGICAL GLOVES

### PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	Green
Thickness	Micro
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	Low
Double Gloving Recommendation	Underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthates, Zinc Dialkyldithiocarbamate

### ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5
Product Code	340058055	340058060	340058065	340058070	340058075	340058080	340058085

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: [insidesalesus@ansell.com](mailto:insidesalesus@ansell.com)  
 Sandel US Tel: 866 764 3327  
 Sandel email: [sandel-cs@ansell.com](mailto:sandel-cs@ansell.com)

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: [csiac@ansell.com](mailto:csiac@ansell.com)  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: [info@ansell.eu](mailto:info@ansell.eu)  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: [apac.medical@ansell.com](mailto:apac.medical@ansell.com)

**India**  
 Tel: +91 22 38124500  
 Email: [info.india@ansell.com](mailto:info.india@ansell.com)  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: [protection@ap.ansell.com](mailto:protection@ap.ansell.com)  
**China**  
 Tel: +86(21) 3827 5005  
 Email: [infochina@ansell.com](mailto:infochina@ansell.com)