

Previously known as: MICROGARD® 2000 Lab Coat - Model 209

## Superior protection from low hazard liquid spray and fine particulates..

### DESIGN

MICROGARD® 2000 range is designed to allow water vapour (perspiration) to escape from the suit yet will withstand saturation of liquid chemicals and filter 100% of particulates down to 0.01 microns in size.

### BOUND SEAMS

Superior strength, liquid and particle barrier

### LAB COAT FEATURES

- Front zip fastening
- Left breast pocket
- Lower right pocket
- Bound seams

### COLOURS

White

### SIZES

S-5XL



## Industries

- Life Sciences
- Metal fabrication
- Food Processing


## Applications

- Clean room cleaning and preparing
- Production line support & maintenance
- Maintenance of plant and machines
- Veterinary Services
- Body painting including surface inspection
- Cleaning of plant and machines


## KEY FEATURES & BENEFITS

- Proven protective barrier to low concentration liquid chemicals, liquid & particulate biological hazards
- Comfortable to wear featuring moisture vapor permeable fabric ("breathable") to help reduce the risk of heat stress
- Optimized fit, improving wearer comfort and safety
- Low Linting

## Material Diagram



## Seam Diagram



## Performance Standards



PARTIAL BODY  
PROTECTION ONLY



TYPE PB [6]



0598



EN 14126

ASTM F1671 EN ISO 13688:2013

## Product Information

| Product Information |                                             |
|---------------------|---------------------------------------------|
| Available Sizes     | S, M, L, XL, 2XL, 3XL                       |
| Color               | White                                       |
| Country Of Origin   | Sri Lanka                                   |
| Packaging Overview  | Individually Packed                         |
| Product Material    | Microporous polyethylene Laminate non-woven |
| Product Reference   | 2000 Lab Coat - Model 209                   |
| Standards Overview  | CE CATEGORY III (Type PB (6)-B, EN 14126)   |
| Seam Type           | Bound                                       |
| Shelf Life          | 5 years                                     |

Suit components & accessories

| Suit Components and Accessories |                                                                                       |
|---------------------------------|---------------------------------------------------------------------------------------|
| Design Features                 | Labcoat, front zip fastening, left breast pocket, lower right pocket and bound seams. |

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

**Europe, Middle East & Africa Region**  
Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00

**North America Region**  
Ansell Healthcare Products LLC  
T: +1 800 800 0444

**Australia**  
Ansell Limited  
T: +61 1800 337 041

**Asia Pacific Region**  
Ansell Global Trading Center  
T: +603 8310 6688

**Latin America & Caribbean Region**  
Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

