

# Non-sterile latex cleanroom gloves, ensuring elbow length-protection and easy double-donning

- **Heightened flexibility and comfort:** BioClean™ Legion BLA3 cleanroom gloves' soft, stretchy formulation offer a perfect blend of flexibility and comfort
- **Added protection:** These ambidextrous cleanroom gloves are designed to provide coverage up to the elbow, granting wearers hand, forearm and wrist protection and greatly reducing cross contamination risks
- **Reduced contamination risks:** BioClean™ Extra BLA3 non-sterile gloves come in non-particulating EasyTear packaging, for further reduced contamination risks
- **Enhanced strength:** These powder-free gloves' textured surface enhances grip performance, while a beaded cuff bolsters strength
- **Increased practicality:** These BioClean™ gloves' robust formulation also makes them ideal for quick, trouble-free double-donning, while their ambidextrous nature ensures added convenience



## KEY FEATURES & BENEFITS

- ◊ **Soft, stretchy film:** For flexible, comfortable cleanroom gloves
- ◊ **Elbow-length design:** Greater coverage, less cross contamination
- ◊ **Non-particulating EasyTear packaging:** Lesser contamination risks

## Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Lab and R&D



# Legion BLA3

Non-sterile Disposable Latex  
Cleanroom Glove

## TECHNICAL DATA SHEET

|                                        | Product Information                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material                               | Natural Rubber Latex                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Color                                  | Natural                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Shape                                  | Ambidextrous                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Cuff                                   | Beaded                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Manufacturing/QMS Audit Standards      | ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms                                                                                                                                                                                                                                                                                                                                                                                        |
| Regulatory/Standard Compliance         | CE 0493, EN 421:2010, EN 455 Part 2, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, Food Contact, Category III, UKCA                                                                                                                                                                                                                                                                                                                                                                            |
| Packaging                              | <p>100 pieces per sealed inner PE bag; one inner PE bag per sealed outer PE bag; 10 outer bags per lined carton (1000 pieces).</p> <p><b>More sustainable packaging:</b> These gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases.</p> <p>*Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.</p> |
| Storage                                | Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.                                                                                                                                            |
| Country of Origin                      | India                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Available sizes                        | S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| External Glove Surface                 | Textured Fingers and Palm                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Internal Glove Surface                 | Chlorinated                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Cleanroom Class                        | Class 100/ISO 5                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Shelf Life                             | 3 Years                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Tested for use with Chemotherapy Drugs | No                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Protein Level                          | 50 µg/g or less of total extractable protein                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Anti-static                            | No                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |



# Legion BLA3

Non-sterile Disposable Latex  
Cleanroom Glove

| Physical Properties                                |             |             |             |               | Testing Method  |
|----------------------------------------------------|-------------|-------------|-------------|---------------|-----------------|
| Sizes                                              | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) | EN ISO 21420    |
| Typical Length (mm/in)                             | 410 / 16    |             |             |               |                 |
| Palm Width (mm/in)                                 | 86/3.38     | 96/3.78     | 106/4.17    | 115/4.53      |                 |
| Freedom from Holes                                 | 1.5 AQL     |             |             |               | EN 374-2        |
| Typical Particle Count<br>≥0.5µm (counts / cm²)    | <1500       |             |             |               | IEST-RP-CC005.4 |
| Target Single Wall Palm<br>Thickness (mm/mil)      | 0.18 / 7.09 |             |             |               | EN 455-2        |
| Target Single Wall<br>Finger Thickness<br>(mm/mil) | 0.21 / 8.27 |             |             |               | EN 455-2        |
| Target Single Wall Cuff<br>Thickness (mm/mil)      | 0.13 / 5.12 |             |             |               | EN 455-2        |
| Force at Break (N)<br>Before Aging                 | 1           |             |             |               | EN 455-2        |
| Force at Break (N)<br>During Aging                 | ≥9          |             |             |               | EN 455-2        |

## IONIC CONTENT

| Concentration in µg/cm <sup>2</sup> | Typical      | Concentration in µg/cm <sup>2</sup> | Typical      |
|-------------------------------------|--------------|-------------------------------------|--------------|
| Ammonium                            | 0.055        | Nitrate                             | 0.277        |
| Bromide                             | Not Detected | Nitrite                             | Not Detected |
| Calcium                             | 0.813        | Phosphate                           | Not Detected |
| Chloride                            | 1.968        | Potassium                           | 0.095        |
| Fluoride                            | 0.005        | Sodium                              | 0.022        |
| Lithium                             | Not Detected | Sulphate                            | 0.538        |
| Magnesium                           | 0.020        | Zinc                                | Not Detected |

## ORDERING INFORMATION

| SIZE        | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|-------------|-------------|-------------|-------------|---------------|
| REORDER NO. | BLA301-S    | BLA302-M    | BLA303-L    | BLA304-XL     |

## PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:  
[www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

