



BioClean™ 2000 Overboots (Model 406)

Ideal choice for aseptic controlled environments & cleanroom protection

BioClean™ 2000 Overboots feature excellent fabric tensile strength that provides good barrier protection with a comfortable design. Ideal for users seeking China Medical certified body protection coverage and comfortable fit for medical use, CDC, EMS as well as vaccine manufacturing.

- **China Medical Class I Certified (YY/T 1633-2019):** Quality assurance in protection
- **Reduced contamination risks:** Features Microporous polyethylene film laminate with SBPP fabric which is a lightweight material, low-linting that offer good barrier protection and are less susceptible to contamination
- **Secured seam:** 3 thread overlocking stitch with blue taped seams



Key Features and Benefits

- **Non-particulate film packaging**
- **Produced and packed in NEBB ISO 7 certified cleanroom facility**
- **Recyclable/ recycled content packaging materials:** BioClean 2000 (Model 406) are packed in recyclable* plastic packaging and delivered in recycled cardboard shipper cases

Note: *Always check the recyclable status of HDPE locally, as this plastic material may not be considered suitable for recycling in your country

Industries

- Controlled and Critical Environments
- Pharmaceutical Manufacturing
- Biotechnology Manufacturing
- Medical Device Manufacturing
- Laboratory Research and Development



BioClean™ 2000 Overboots (Model 406)

TECHNICAL DATA SHEET

PRODUCT INFORMATION

Material	Microporous polyethylene laminate non-woven
Color	White
Audit Standards	ISO 14001, ISO 9001
Standards	GB24541-2009
Packaging Overview	One pair per sealed inner PE bag; 100 inner bags per lined carton (100 pairs).
Storage	Store in a dry cool place 40°C away from direct sunlight and fluorescent light.
Country Of Origin	China
Shelf Life	5年
Construction	三线锁边缝线, 带蓝色压胶接缝
Characteristics	*NOTE: The material is static dissipative with charge decay time of T 50 < 4 seconds. (Tested in accordance with EN 1149-3:2006)

SIZE CHART

Universal



BioClean™ 2000 Overboots (Model 406)

MATERIAL PERFORMANCE TEST RESULTS

TEST	RESULT
Tensile Strength Cross Direction (CD)	58N
Tensile Strength Machine Direction (MD)	140N

ORDERING INFORMATION

	SIZE	2000-406
Overboots (Model 406)	REORDER NO.	WH20-T-BC-406-00-NS1

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

Please see product validation pack or contact Ansell customer service for specific data on use of garments with cytotoxic drugs. Garments used for protection against such drugs must be selected specifically for the type of chemicals used.

