

DISPOS-A-GLOVE™ Sterile Single

Essential Series

These single-packed synthetic medical exam gloves are ideal for brief, sterile applications

- **Sterile:** Packaged as a single sterile glove per pouch, these powder-free gloves are suitable for sterile medical exam (non-surgical) procedures
- **Reduced allergy risk:** This synthetic glove is made from copolymer materials (ethylene vinyl acetate film), making it suitable for use by individuals with allergies to natural rubber latex (Type I allergies)
- **Heightened sensitivity:** This glove also offers good tactile sensitivity to the wearer
- **Double-gloving friendly:** DISPOS-A-GLOVE™ can suitably be used as a glove liner for increased protection
- **Cost-effective:** This synthetic glove for sterile use is offered at a competitive price point
- **Easy donning:** This glove is packaged on backing paper, making it easy to don

Recommended For

- Sterile procedures outside of the O.R.
- Protection from Type I latex allergy in HCW's or patients
- Brief, low-risk examination procedures



DISPOS-A-GLOVE

DISPOS-A-GLOVE™ Sterile Single

Essential Series

PRODUCT INFORMATION

Material	Ethylene Vinyl Acetate	Not Made From Natural Rubber Latex	Yes
Color	Clear	Cuff Style	Straight
Cuff Length	Standard	External Glove Surface	Smooth
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.03 / 1.2
Palm Thickness (mm/mil)	0.03 / 1.2	Allergy Prevention	Latex (Type I)
Available Sizes	S (6.5 - 7), M (7.5 - 8), L (8.5 - 9)	Tested For Use With Chemotherapy Drugs	No
Sterile	Yes	Antistatic	Not Tested
Glove Length (mm/inches)	278 / 10.9	Product Segmentation	Essential Series

PRODUCT INFORMATION

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)
Product Code	MDG601	MDG701	MDG801

PACKAGING AND STORAGE

Packaging	1 glove/paper pouch 100 gloves per box 14 boxes per carton/case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

CE Number, EN 455 Part 1, EN 455 Part 3, ISO 13485, ISO 10993-10:2002, Medical Device Regulation (EU) 2017/745 Class I



Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: csllac@ansell.com
Latin America & Caribbean Region
BR +55 (11) 3356-3100
latam.medical@ansell.com
MX: +52 442 296 2050
latam.medical@ansell.com
CO: +57 1 288 3247
latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
MY +603 8310 6688
apac.medical@ansell.com
CN: +86 21 38275000
infochina@ansell.com
JP: +813 5549 8151
info.medical.jp@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com