

# MICRO-TOUCH™ Nitrile E.P.

P e r f o r m a n c e   S e r i e s

## Durable design with an extended cuff for high-risk protection within healthcare environments

- **Enhanced grip:** Textured for better grip in both wet and dry clinical settings
- **Added protection:** The extended 295mm/11.6" cuff provides added protection over the wrist and forearm
- **Barrier and chemical protection:** Offering effective barrier protection against punctures and a wide range of chemicals
- **Reduced allergy risk:** Prevents latex Type I allergies
- **Specialized protection:** Tested for use with chemotherapy drugs\*

\*see product packaging for details

### Recommended For

- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients
- Sustained, high-risk examination procedures
- Extra protection over the wrist and arm

### FEATURED TECHNOLOGIES



## MICRO-TOUCH™ Everyday, Everywhere Protection

MICRO-TOUCH™ premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH™ Nitrile E.P.

Performance Series

## PRODUCT INFORMATION

|                                         |                                                                    |                                        |                                                       |
|-----------------------------------------|--------------------------------------------------------------------|----------------------------------------|-------------------------------------------------------|
| Material                                | Nitrile                                                            | Not Made From Natural Rubber Latex     | Yes                                                   |
| Color                                   | Blue                                                               | Cuff Style                             | Beaded                                                |
| Cuff Length                             | Extended                                                           | External Glove Surface                 | Textured Fingers                                      |
| Freedom from Holes (Inspection level I) | 1.5 AQL                                                            | Finger Thickness (mm/mil)              | 0.16/6.3                                              |
| Palm Thickness (mm/mil)                 | 0.11/4.3                                                           | Allergy Prevention                     | Latex (Type I)                                        |
| Available Sizes                         | XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10) | Tested For Use With Chemotherapy Drugs | Yes, in accordance with ASTM D6978 and US FDA cleared |
| Silicone Free                           | No                                                                 | Sterile                                | No                                                    |
| Antistatic                              | Not Tested                                                         | Glove Length (mm/inches)               | 295 / 11.6                                            |
| Product Segmentation                    | Performance Series                                                 |                                        |                                                       |

## ORDERING INFORMATION

| Size         | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|--------------|-------------|-------------|-------------|---------------|
| Product Code | 6034050      | 6034051     | 6034052     | 6034053     | 6034054       |

## PACKAGING AND STORAGE

|                      |                                                                                                          |
|----------------------|----------------------------------------------------------------------------------------------------------|
| Packaging            | 100 gloves per dispenser 10 dispensers per case 1000 gloves per case                                     |
| Shelf Life           | 3 Years                                                                                                  |
| Storage Instructions | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition. |

## STANDARDS AND CERTIFICATIONS

ASTM D6319, ASTM D6978, EN 374:2003, EN ISO 21420:2020, FDA21 CFR 177-2600-US Food Contact Approved, ISO 13485, ISO 9001, ISO 10993-10:2002

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email:  
insidesalesus@ansell.com

Latin America & Caribbean  
Region  
BR +55 (11) 3356-3100  
latam.medical@ansell.com  
MX: +52 442 296 2050  
latam.medical@ansell.com  
CO: +57 1 288 3247  
latam.medical@ansell.com

Europe, Middle East &  
Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email:  
apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email:  
protection@ap.ansell.com