


## Provides 30% longer, up-to-the-elbow personal protection in high-risk situations

- **Extended protection:** The MICROFLEX™ 93-244 nitrile medical exam gloves provide 30% longer, up-to-the-elbow (390mm long) personal protection for high-risk situations
- **Reduced allergy risk:** Manufactured without accelerators, helping to protect users from Type IV skin allergies
- **Chemical & Chemo protection:** Applicable to a wide range of healthcare situations, these disposable gloves have been tested against commonly used chemicals in the Sterilization Department and chemotherapy drugs



### Industries

- Healthcare

### Recommended For

- Sterile Processing Department
- Preventing the risk of cross-contamination
- Laboratory analysis
- Central Supply
- Biohazard response

## TECHNICAL DATA SHEET

|                        | Product Information                                                                                                                                                                                                                                                                                                                                                          |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material               | Nitrile                                                                                                                                                                                                                                                                                                                                                                      |
| External Glove Surface | Fully Textured                                                                                                                                                                                                                                                                                                                                                               |
| Audit Standards        | ISO 13485, ISO 9001                                                                                                                                                                                                                                                                                                                                                          |
| Product Certification  | Medical Device Regulation (EU) 2017/745 Class I, Personal Protective Equipment Regulation (EU) 2016/425 Category III risks, CE 0493, EU Food Contact Approved, EC 2023/2006, EN 16523-1, EN ISO 374-1:2016 Type B KPT, EN ISO 374-5:2016 Virus, EN 455 1-4, EN ISO 21420:2020, FDA21 CFR 177-2600-US Food Contact Approved, ISO 13485, ISO 9001, UKCA, ASTM D6319, EN 1149-5 |
| Packaging Overview     | 100 gloves per box<br>10 box per case<br>1000 gloves per case                                                                                                                                                                                                                                                                                                                |
| Storage Instructions   | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                                                                                                                                                                                                                                                                     |
| Country Of Origin      | Malaysia                                                                                                                                                                                                                                                                                                                                                                     |
| Product Segmentation   | High Risk                                                                                                                                                                                                                                                                                                                                                                    |
| Antistatic             | Yes                                                                                                                                                                                                                                                                                                                                                                          |
| Silicone Free          | Yes                                                                                                                                                                                                                                                                                                                                                                          |

## PHYSICAL PROPERTIES

|                                         | Typical Values    |             | Testing Method      |
|-----------------------------------------|-------------------|-------------|---------------------|
| Length (mm/inches)                      | 390 / 15.4        |             | EN 455-2,ASTM D6319 |
| Freedom from Holes (Inspection level I) | 1.5 AQL           |             | ASTM D5151,EN 455-1 |
| Palm Thickness (mm/mils)                | 0.13 mm / 5.1 mil |             | EN 455-2,ASTM D6319 |
| Finger Thickness (mm/mils)              | 0.17 mm / 6.7 mil |             | EN 455-2,ASTM D6319 |
|                                         | BEFORE AGING      | AFTER AGING |                     |
| Elongation at Break (%)                 | ≥ 500             | ≥ 400       | ASTM D412 & D573    |
| Force at break (N)                      | ≥ 6               | ≥ 6         | EN 455-2            |

ORDERING INFORMATION

| Size         | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) | XXL (10.5 - 11) |
|--------------|--------------|-------------|-------------|-------------|---------------|-----------------|
| Product Code | 93244060     | 93244070    | 93244080    | 93244090    | 93244100      | 93244110        |

Performance Standards and Regulatory Compliance

CE 0493

MD



EN ISO 374-1:2016  
Type B

  
K P T

EN ISO 374-5

  
VIRUS

ISO

ISO 13485

ISO

ISO 9001

UK

CA

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

**Europe, Middle East & Africa Region**  
Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

**Asia Pacific Region**  
MY: +603 8310 6688  
apac.medical@ansell.com  
CN: +86 21 38275000  
infochina@ansell.com  
JP: +813 5549 8151  
info.medical.jp@ansell.com

**North America Region**  
Ansell Healthcare Products LLC  
T: +1 800 800 0444  
F: +1 800 800 0445

**Latin America & Caribbean Region**  
BR: +55 (11) 3356-3100  
latam.medical@ansell.com  
MX: +52 442 296 2050  
latam.medical@ansell.com  
CO: +57 1 288 3247  
latam.medical@ansell.com

**Australia**  
Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

**Russia**  
Ansell PYC  
Ten. +7 495 258 13 16

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.