

## Durable gloves for heavy-duty jobs requiring extra protection

- **Robust design:** MICROFLEX™ E-Grip™ Max L92 is a medical exam glove with a thick design for enhanced durability, protecting against punctures, rips and tears
- **Comfortable fit:** Made of a soft, comfortable latex for comfort during extended wear
- **Easy donning/doffing:** Because the glove is double-chlorinated, it is easy to don and doff
- **Food handling compliance:** Suitable for a wide range of food-handling applications



**CAUTION:** This product contains natural rubber latex, which may cause allergic reactions

### Industries

- Automotive Aftermarket
- Automotive
- Healthcare
- Fire and Rescue

### Recommended For

- General Purpose Auto Aftermarket
- Standard, moderate risk examination procedures
- Inspection (picking, checking and marking)
- Picking parts and accessories

## TECHNICAL DATA SHEET

|                        | Product Information                                                                                                     |
|------------------------|-------------------------------------------------------------------------------------------------------------------------|
| Material               | Natural Rubber Latex                                                                                                    |
| External Glove Surface | Textured Fingers                                                                                                        |
| Audit Standards        | ISO 13485:2003                                                                                                          |
| Product Certification  | ASTM D3578, EN ISO 21420:2020, FDA21 CFR 177-2600-US Food Contact Approved, US QSR/FDA 510(k) Medical Examination Grade |
| Packaging Overview     | 100 gloves per dispenser<br>10 dispensers per case<br>1000 gloves per case                                              |
| Storage Instructions   | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                |
| Country Of Origin      | Malaysia                                                                                                                |
| Product Segmentation   | Multi-Purpose                                                                                                           |
| Antistatic             | Not Tested                                                                                                              |
| Silicone Free          | No                                                                                                                      |

## PHYSICAL PROPERTIES

|                                         | Typical Values |             | Testing Method      |
|-----------------------------------------|----------------|-------------|---------------------|
| Length (mm/inches)                      | 245 / 9.65     |             | ASTM D3767,EN 420   |
| Freedom from Holes (Inspection level I) | 1.5 AQL        |             | ASTM D5151,EN 455-1 |
| Palm Thickness (mm/mils)                | 0.13/5.1       |             | ASTM D3767,EN 420   |
| Finger Thickness (mm/mils)              | 0.16/6.3       |             | ASTM D3767,EN 420   |
|                                         | BEFORE AGING   | AFTER AGING |                     |
| Elongation at Break (%)                 | ≥ 750          | ≥ 600       | ASTM D412           |

ORDERING INFORMATION

| Size         | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|--------------|-------------|-------------|-------------|---------------|
| Product Code | L920         | L921        | L922        | L923        | L924          |

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

Asia Pacific Region

MY: +603 8310 6688  
apac.medical@ansell.com  
CN: +86 21 38275000  
infochina@ansell.com  
JP: +813 5549 8151  
info.medical.jp@ansell.com

North America Region

Ansell Healthcare Products LLC  
T: +1 800 800 0444  
F: +1 800 800 0445

Latin America & Caribbean Region

BR: +55 (11) 3356-3100  
latam.medical@ansell.com  
MX: +52 442 296 2050  
latam.medical@ansell.com  
CO: +57 1 288 3247  
latam.medical@ansell.com

Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

Russia

Ansell PYC  
Ten. +7 495 258 13 16

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

