



Advance BASL

Sterile Disposable Latex Cleanroom Glove

Sterile latex cleanroom gloves, delivering comfort and flexibility

- **Added protection:** BioClean™ Advance BASL sterile cleanroom gloves' flexible, yet sturdy, design facilitates double-donning, for an additional protective layer
- **Enhanced features:** These cleanroom gloves feature a textured finger surface, improving grip levels, while a beaded cuff delivers added stability and prevents tearing when donning
- **Reduced contamination risks:** BioClean™ Advance BASL powder-free sterile latex gloves are gamma irradiation sterilized and come in non-particulating EasyTear packaging, lessening particulate contamination risks in cleanroom environments



CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions

KEY FEATURES & BENEFITS

- ◊ **Double-donning-friendly design:** For convenient daily hand protection
- ◊ **Textured finger surface:** For sterile gloves with strengthened grip
- ◊ **Non-particulating EasyTear packaging:** Lesser contamination risks

Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Sterile Processing Department

Recommended For

- Clean room cleaning and preparing
- Preventing the risk of cross-contamination
- Transferring liquids and solids
- Chemical Splash
- Pharmaceutical manufacturing
- Medical Device Manufacturing
- Extra comfort during long, sustained procedures
- Biotechnologies Blending operations



Advance BASL

Sterile Disposable Latex Cleanroom Glove

TECHNICAL DATA SHEET

	Product Information
Material	Natural Rubber Latex
Color	Natural
Shape	Hand Specific
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory/Standard Compliance	CE 0493, UKCA, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 455 Part 2, Food Contact, Category III
Packaging	1 pair per inner PE wallet; 1 wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs). More sustainable packaging: These sterile latex gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases. *Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Available sizes	6, 6.5, 7, 7.5, 8, 8.5, 9
External Glove Surface	Textured Fingers
Internal Glove Surface	Chlorinated
Sterilization Method	GAMMA irradiation (25 kGy)
Sterility Assurance Level	10 ⁻⁶
Cleanroom Class	Class 10/ISO Class 4, EU GMP Grade A/B and other sterile cleanrooms
Shelf Life	5 Years
Tested for use with Chemotherapy Drugs	No
Protein Level	50 µg/g or less of total extractable protein
Anti-static	No



Advance BASL

Sterile Disposable Latex Cleanroom Glove

Physical Properties								Testing Method
Sizes	6	6.5	7	7.5	8	8.5	9	EN ISO 21420
Typical Length (mm/in)	295 / 11.61							
Palm Width (mm/in)	77/3	83/3.3	89/3.5	95/3.7	102/4	109/4.3	114/4.5	
Freedom from Holes	1.5 AQL Performance Level 2							EN 374-2
Typical Particle Count ≥0.5µm (counts / cm²)	<2000							IEST-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.21							EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.24							EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.20							EN 455-2
Force at Break (N) During Aging	≥ 9N							EN 455-2

IONIC CONTENT

Concentration in µg/cm ²	Typical	Concentration in µg/cm ²	Typical
Ammonium	< 0.05	Nitrate	< 0.50
Bromide	< 0.01	Nitrite	< 0.01
Calcium	< 0.80	Phosphate	< 0.01
Chloride	< 1.80	Potassium	< 0.08
Fluoride	< 0.01	Sodium	< 0.08
Lithium	< 0.01	Sulphate	< 0.22
Magnesium	< 0.05	Zinc	< 0.60

ORDERING INFORMATION

SIZE	6	6.5	7	7.5	8	8.5	9
REORDER NO.	BASL6.0	BASL6.5	BASL7.0	BASL7.5	BASL8.0	BASL8.5	BASL9.0

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

