

GAMMEX™ PI Hybrid Micro

S U R G I C A L G L O V E S

Polyisoprene and neoprene blend surgical gloves, made with HYBRID™ Synthetic Blend Technology, delivers comfort, durability and enhanced sensitivity for precision performance

- **Enhanced comfort and durability:** Combining polyisoprene and neoprene to deliver a latex-like comfort and heightened durability
- **Heightened sensitivity:** Being up to 24% thinner than Ansell's standard nonlatex surgical gloves, GAMMEX™ PI Hybrid Micro offer exceptional level of tactile sensitivity and dexterity for delicate procedures that requires utmost precision
- **Double-gloving versatility:** This brown micro glove is suitable to be donned either as an outer glove; or inner glove when double gloving for easy breach detection
- **Reduced allergy risks:** This synthetic blend material eliminates Type I latex allergy risks, and is free from Zinc 2-mercaptobenzothiazole-free (ZMBT-free), diphenylguanidine-free (DPG-free) and cetylpyridinium chloride-free (CPC-free) to minimize the risk of chemical Type IV allergies and sensitivities to these chemical accelerators

Recommended For

- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- Use as an underglove when double gloving
- All surgeries, especially delicate, precision performance procedures such as ophthalmology, micro and cardiovascular requiring high tactile sensitivity and dexterity

FEATURED TECHNOLOGIES



GAMMEX™ Innovation that Protects

Driven by innovation, GAMMEX™ surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX™ PI Hybrid Micro

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Polyisoprene, Neoprene
Powder Content	Powder-Free
Color	Brown
Thickness	Micro
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, CE 0493 (PPE) + CE 2797 (MD) , EN 16523-1, EN 420:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC), Sodium di-n-butyldithiocarbamate (SDBC), Xanthogen Polysulphides

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340002055	340002060	340002065	340002070	340002075	340002080	340002085	340002090

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico
 Tel: +52 442 296 2050
 Email: latam.medical@ansell.com
Ansell Brasil Ltda
 +55 (11) 3356-3100
 Email: latam.medical@ansell.com
Ansell Colombia SAS
 Tel: +57 1 288 3247
 Email: latam.medical@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu
Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com