

Non-sterile latex disposable gloves, for hand protection offering comfort and tactility

- **Assured cleanliness:** MICROFLEX™ CE4-200 latex disposable gloves are designed to be compatible with Class 10 (ISO 4) cleanroom environments
- **Reduced contamination risks:** They are ultra-clean processed to produce low surface particulates, ionic extractables and non-volatile residues
- **Elevated comfort and tactility:** These cleanroom gloves' latex formulation ensures a blend of comfort, tactility and dexterity, for easy-to-use everyday hand protection
- **CAUTION:** This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.



*tested in accordance with IEST-RP-CC-005.3 and FTMS 101-C, 4046

KEY FEATURES & BENEFITS

- ◊ **Class 10 (ISO 4) cleanroom compatibility:** Assured cleanliness*
- ◊ **Unique latex formulation:** Dexterity, tactile sensitivity and comfort
- ◊ **Ultra-clean process:** Low surface particulates, ionic extractables and non-volatile residues

*manufactured in an integrated facility certified to ISO 9001 standards

Industries

- Pharmaceutical Manufacturing
- Biotechnology Manufacturing
- Medical Device Manufacturing
- Electronics manufacturing

Recommended For

- Clean room cleaning and preparing
- Spill or leakage cleanup
- Transferring liquids and solids
- Weighing and dispensing of solid and liquid raw materials
- Loading centrifuges and chromatography columns
- Blending, compounding solids and liquids



MICROFLEX™

CE4-200

Non-sterile Disposable Latex
Cleanroom Glove

TECHNICAL DATA SHEET

	Product Information
Material	Natural Rubber Latex
Color	Natural
Shape	Ambidextrous
Cuff	Beaded
Regulatory/Standard Compliance	EN 420:2003 + A1:2009, ISO 13485, ISO 9001
Packaging	100 gloves per inner polybag, 10 inner polybags per lined carton
Available sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)
Powder Content	Powder-Free
Internal Glove Surface	Powder-Free
Cleanroom Class	Class 10/ISO 4
Shelf Life	3 years
Tested for use with Chemotherapy Drugs	No
Anti-static	No

Physical Properties		Testing Method
Typical Length (mm/in)	300 / 12	EN 420/ASTM D3767
Freedom from Holes	1.5 AQL	EN 455-1/ASTM D3577
Typical Particle Count $\geq 0.5\mu\text{m}$ (counts / cm^2)	< 850 particle count/ cm^2	IEST-RP-CC005.4


IONIC CONTENT

Concentration in $\mu\text{g}/\text{cm}^2$	Typical	Concentration in $\mu\text{g}/\text{cm}^2$	Typical
Bromide	<0.20	Nitrate	<0.50
Calcium	<0.50	Phosphate	<0.20
Chloride	<1.00	Potassium	<0.20
Fluoride	<0.20	Sodium	<0.20
Magnesium	<0.03	Sulphate	<0.20

ORDERING INFORMATION

SIZE	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
REORDER NO.	CE4-200-XS	CE4-200-S	CE4-200-M	CE4-200-L	CE4-200-XL

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.