



## Maxima BLLS

Sterile Disposable Latex Cleanroom Glove

# Shoulder-length sterile latex gloves, providing comprehensive personal protection

- **Added protection:** BioClean™ Maxima BLLS are anatomically shaped, shoulder-length sterile latex gloves specially designed to provide extra protection
- **Elevated comfort:** These cleanroom gloves' latex material is inherently comfortable and flexible, for easy-to-wear personal protective equipment
- **Improved grip:** These sterile cleanroom gloves also boast a textured palm and fingers, enhancing grip performance
- **Reduced contamination risks:** These BioClean™ gloves have undergone gamma irradiation sterilization to minimize risk of contamination entering the cleanroom environment
- **Genuine adaptability:** In addition, BioClean™ Maxima BLLS shoulder-length PPE gloves are designed to facilitate double-donning



## KEY FEATURES & BENEFITS

- **600mm/24" length:** For additional arm coverage and protection
- **Latex composition:** Combination of flexibility and comfort
- **Textured surface:** For latex gloves with improved grip

## Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Lab and R&D
- Sterile Processing Department





## Maxima BLLS

Sterile Disposable Latex Cleanroom Glove

### TECHNICAL DATA SHEET

|                                   | Product Information                                                                                                                                   |
|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material                          | Natural Rubber Latex                                                                                                                                  |
| Color                             | Natural                                                                                                                                               |
| Shape                             | Hand Specific                                                                                                                                         |
| Cuff                              | Beaded                                                                                                                                                |
| Manufacturing/QMS Audit Standards | ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms                                    |
| Regulatory                        | CE 0598, EN 420:2003 + A1:2009, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 455 Part 2, Food Contact, Category III                          |
| Packaging                         | One pair per inner PE wallet; one wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 10 outer bags per lined carton (100 pairs) |
| Storage                           | Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.                                                                   |
| Country of Origin                 | Malaysia                                                                                                                                              |
| Available sizes                   | 6, 6.5, 7, 7.5, 8, 8.5, 9                                                                                                                             |
| Powder Content                    | Powder-Free                                                                                                                                           |
| External Glove Surface            | Textured Fingers and Palm                                                                                                                             |
| Internal Glove Surface            | Chlorinated                                                                                                                                           |
| Sterilization Method              | GAMMA irradiation (25 kGy)                                                                                                                            |
| Sterilization Minimum Dose        | 25kGy                                                                                                                                                 |
| Sterility Assurance Level         | 10 <sup>-6</sup>                                                                                                                                      |
| Cleanroom Class                   | Class 10/ISO Class 4 & EU GMP Grade A/B and other sterile cleanrooms                                                                                  |
| Shelf Life                        | Five (5) years from date of manufacture.                                                                                                              |
| Protein Level                     | 50 µg/g or less of total extractable protein                                                                                                          |



# Maxima BLLS

Sterile Disposable Latex Cleanroom Glove

| Physical Properties                                |                             |        |        |        |       |         |         |         | Testing Method  |
|----------------------------------------------------|-----------------------------|--------|--------|--------|-------|---------|---------|---------|-----------------|
| Sizes                                              | 6                           | 6.5    | 7      | 7.5    | 8     | 8.5     | 9       | 10      | EN 420          |
| Typical Length (mm/in)                             | 600 / 24                    |        |        |        |       |         |         |         |                 |
| Palm Width (mm/in)                                 | 77/3                        | 83/3.3 | 89/3.5 | 95/3.7 | 102/4 | 108/4.2 | 115/4.5 | 120/4.7 |                 |
| Freedom from Holes                                 | 1.5 AQL Performance Level 2 |        |        |        |       |         |         |         | EN 374-2        |
| Typical Particle Count<br>≥0.5µm (counts / cm²)    | <1200                       |        |        |        |       |         |         |         | IEST-RP-CC005.4 |
| Target Single Wall Palm<br>Thickness (mm/mil)      | 0.18 / 7.09                 |        |        |        |       |         |         |         | EN 455-2        |
| Target Single Wall<br>Finger Thickness<br>(mm/mil) | 0.20 / 7.87                 |        |        |        |       |         |         |         | EN 455-2        |
| Target Single Wall Cuff<br>Thickness (mm/mil)      | 0.12 / 4.72                 |        |        |        |       |         |         |         | EN 455-2        |
| Force at Break (N)<br>During Aging                 | ≥9N                         |        |        |        |       |         |         |         | EN 455-2        |

## IONIC CONTENT

| Concentration in µg/cm <sup>2</sup> | Typical      | Concentration in µg/cm <sup>2</sup> | Typical      |
|-------------------------------------|--------------|-------------------------------------|--------------|
| Ammonium                            | Not Tested   | Nitrate                             | 0.01         |
| Bromide                             | Not Detected | Nitrite                             | 0.01         |
| Calcium                             | 0.07         | Phosphate                           | Not Detected |
| Chloride                            | 0.12         | Potassium                           | 0.03         |
| Fluoride                            | Not Detected | Sodium                              | 0.01         |
| Lithium                             | Not Detected | Sulphate                            | 0.01         |
| Magnesium                           | 0.02         | Zinc                                | Not Detected |

## ORDERING INFORMATION

| SIZE        | 6      | 6.5    | 7      | 7.5    | 8      | 8.5    | 9      | 10      |
|-------------|--------|--------|--------|--------|--------|--------|--------|---------|
| REORDER NO. | BLLS60 | BLLS65 | BLLS70 | BLLS75 | BLLS80 | BLLS85 | BLLS90 | BLLS100 |

## PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:  
[www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

