

Chemical-resistant sterile nitrile cleanroom gloves, perfect for skilled handling

- **Chemical protection:** BioClean™ Emerald BENS are resistant to a range of chemicals and chemotherapy drugs (PPE Category III, EN ISO Type B, ASTM D6978)
- **Reduced allergy risk:** The gloves are powder-free, latex-free and accelerator-free minimizing allergy risks, notably Type I latex and Type IV chemical allergies and sensitivities
- **Clean processed:** BioClean™ Emerald BENS undergo deionized washing to further reduce surface particulates, ionic extractables and non-volatile residues



KEY FEATURES & BENEFITS

- **Double gloving properties:** Ideal as undergloves for double gloving requirements
- **Enhanced dexterity:** Hand specific design for ergonomic fit and better handling during skilled tasks
- **Sterility Assurance:** Gamma sterilized with sterility assurance level (SAL) of 10⁻⁶

Industries

- Biotechnology Manufacturing
- Controlled and Critical Environments
- Laboratory Research and Development
- Medical Device Manufacturing
- Pharmaceutical Manufacturing
- Production and Manufacturing

Recommended For

- Chemical Splash
- Sterile Processing Department
- Chemical Handling
- Biotechnologies Blending operations
- Transferring liquids and solids
- Sample taking and processing
- Aseptic filling and labelling



Emerald BENS

Sterile Disposable Nitrile Cleanroom Glove

TECHNICAL DATA SHEET

	Product Information
Material	Nitrile
Color	Green
Shape	Hand Specific
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory	ASTM D6978, CE 0493, EN 421:2010, EN 455 Part 2, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, Food Contact, Category III, UKCA
Packaging	1 pair per inner wrap; 1 inner wrap per pouch; 10 pouches per outer bag; 20 outer bags per lined carton (200 pairs) More sustainable packaging: These sterile nitrile gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases. *Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Country of Origin	Malaysia
Available sizes	6, 6.5, 7, 7.5, 8, 8.5, 9
Powder Content	Powder-Free
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Chlorinated
Sterilization Method	Gamma Irradiation
Sterilization Minimum Dose	25kGy
Sterility Assurance Level	10 ⁻⁶
Cleanroom Class	Class 10/ISO Class 4 & EU GMP Grade A/B and other sterile cleanrooms
Shelf Life	Five (5) years from date of manufacture.
Tested for use with Chemotherapy Drugs	Yes
Protein Level	N/A: contains no natural rubber latex
Anti-static	Yes
Vulcanization Chemical Accelerators	<ul style="list-style-type: none">None



Emerald BENS

Sterile Disposable Nitrile Cleanroom Glove

Physical Properties								Testing Method
Sizes	6	6.5	7	7.5	8	8.5	9	EN ISO 21420
Typical Length (mm/in)	300 / 12							
Palm Width (mm/in)	77/3.0	83/3.3	89/3.5	95/3.7	102/4.0	108/4.3	114/4.5	
Freedom from Holes	1.5 AQL Performance							EN 374-2
Typical Particle Count ≥0.5µm (counts / cm²)	<1200							IEST-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.11 / 4.33							EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.15 / 5.90							EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.08 / 3.15							EN 455-2
Ultimate tensile strength (MPa) During Aging	Min. 15							ASTM D412-06a
Force at Break (N) During Aging	≥ 6N							EN 455-2

IONIC CONTENT

Concentration in µg/cm ²	Typical	Concentration in µg/cm ²	Typical
Ammonium	Available On File	Nitrate	<0.25
Bromide	<0.01	Nitrite	Available On File
Calcium	<0.60	Phosphate	<0.01
Chloride	<0.35	Potassium	<0.08
Fluoride	<0.01	Sodium	<0.05
Lithium	<0.01	Sulphate	<0.05
Magnesium	<0.01	Zinc	<0.10

ORDERING INFORMATION

SIZE	6	6.5	7	7.5	8	8.5	9
REORDER NO.	BENS6.0	BENS6.5	BENS7.0	BENS7.5	BENS8.0	BENS8.5	BENS9.0

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

